

OA & INJURY
SOLUTIONS CATALOGUE



ÖSSUR APAC OFFICE

<p>Össur CHINA 2F, W16 B No. 1801 Hongmei Road 200233, Shanghai, China</p> <p>Tel: 86 21 6127 1707 Fax: 82 21 6127 1799 asia@ossur.com WWW.OSSUR.ASIA</p>	<p>Össur ASEAN 60 Paya Lebar Road Paya Lebar Square #07-51 Singapore 409051</p> <p>Tel: 65 6397 7190</p> <p>enquiriessg@ossur.com WWW.OSSUR.SG</p>	<p>Össur INDIA 207B, 2nd Floor, Rectangle one, District Centre, Saket, New Delhi - 110017 India</p> <p>Toll Free: 1800 3000 6220</p> <p>ossurindia@ossur.com WWW.OSSUR.IN</p>	<p>Össur AUSTRALIA Unit 15, 114 Merrindale Drive Croydon Victoria 3136 Australia</p> <p>Tel: 61 3 8761 6408 Fax: 61 3 9761 6067 infomelbourne@ossur.com WWW.OSSUR.COM.AU</p>	<p>Össur KOREA No.207, 2Fl., Samsung IT Valley, 27, Digital-ro 33-gil, Guro-gu, Seoul 08380, Korea</p> <p>Tel: 82 2 868 5482 Fax: 82 2 868 5483 ossurkorea@ossur.com WWW.OSSUR.KR</p>
--	---	--	---	--

FOLLOW ÖSSUR ON



WWW.OSSUR.ASIA

Össur APAC
2F, W16 B
No. 1801 Hongmei Road
200233, Shanghai, China

TEL: 86 21 6127 1707
FAX: 86 21 6127 1799
asia@ossur.com

WELCOME TO THE ÖSSUR ASIA-PACIFIC 2018 OA & INJURY SOLUTIONS CATALOGUE

HIGH QUALITY, HIGH PERFORMANCE PRODUCTS ... FOR A LIFE WITHOUT LIMITATIONS.

Össur's diverse range of OA & Injury Solutions products offers expert treatment for a range of conditions. Combining scientifically advanced products with the highest standards of quality assurance and customer service, our destination is Life Without Limitations.

Your business is important to us. We are committed to offering products and services that will enhance your operation, from highly innovative products backed by exceptional service and educational support, to more common items that represent great value for money. We hope that our OA & Injury Solutions catalogue will guide you through this comprehensive range of products.

Should you need any assistance, please do not hesitate to contact our Customer Service team or visit our website, www.ossur.asia, for additional information.

OA & INJURY SOLUTIONS CATALOGUE

TABLE OF CONTENTS

INTRODUCTION TO OA SOLUTIONS	4
UNLOADER® HIP	6
UNLOADER ONE® LITE	8
UNLOADER ONE®	10
UNLOADER ONE® CUSTOM	13
UNLOADER® CUSTOM	16
CTI® OA.....	18
FORM FIT® OA WRAPAROUND	20
INTRODUCTION TO INJURY SOLUTIONS	22
MIAMI LSO™	24
MIAMI LUMBAR® POSTEO	26
MIAMI LUMBAR TLSO®	27
MIAMI J® ADVANCED	28
MIAMI J®	29
MIAMI JTO® THORACIC EXTENSION	31
OCCIAN® BACK & OCCIAN® BACK JR	32
MIAMI JR®	33
PAPOOSE®	34
PHILADELPHIA® ADJUSTABLE COLLAR	35
PHILADELPHIA® TRACHEOTOMY COLLAR.....	36
UNIVERSAL CERVICAL COLLAR	37
FOAM CERVICAL COLLAR.....	37
NECLOC®	38
NECLOC® KIDS	38
PATRIOT® COLLAR.....	39
FORM FIT® BACK SUPPORT, AIR	40
FORM FIT® BACK SUPPORT	41
BACK SUPPORT SUSPENDERS.....	42
UPPER EXTREMITY	43
INNOVATOR X®	44
FORM FIT® SHOULDER BRACE	45
SHOULDER ABDUCTION SLING	46
PREMIUM SHOULDER IMMOBILIZER.....	46
PADDED STRAP SHOULDER IMMOBILIZER.....	47
PREMIUM CONTACT CLOSURE ARM SLING	47
BUCKLE CLOSURE ARM SLING.....	48
BERGSCHULTZ UNIVERSAL ARM SLING.....	48
PREMIUM PADDED ARM SLING	49
FRONT CLOSURE CLAVICLE SUPPORT.....	49
FIGURE-8 CLAVICLE SPLINT	50

TENNIS ELBOW SUPPORT WITH HOT/COLD GEL THERAPY	50
FORM FIT® TENNIS ELBOW	51
AIRFORM® TENNIS ELBOW SUPPORT.....	51
EXOFORM® WRIST	52
EXOFORM® CARPAL TUNNEL WRIST	52
FORM FIT® WRIST / FORM FIT® WRIST & FOREARM	53
FORM FIT® WRIST & THUMB, UNIVERSAL	54
FORM FIT® THUMB SPICA	55
UNIVERSAL THUMB SPLINT	56
6" THUMB AND WRIST	56
HIP	57
REBOUND® HIP	58
KNEE	59
REBOUND® CARTILAGE	60
REBOUND® PCL.....	63
CTI® CUSTOM	65
CTI®	68
REBOUND® DUAL	71
PARADIGM®	73
FORM FIT® KNEE MCL	74
C180™ ROCKET	75
REBOUND® KNEE BRACE.....	76
INNOVATOR®	77
REBOUND POST OP KNEE	78
EXOFORM® KNEE IMMOBILIZER.....	79
FORM FIT® KNEE IMMOBILIZER	80
FX™ PATELLA STABILIZER	81
FORM FIT® HINGED LATERAL J	82
FORM FIT® KNEE, HINGED	83
FORM FIT® KNEE ROM	84
NEOPRENE HINGED KNEE SUPPORT	85
NEOPRENE WRAPAROUND HINGED KNEE SUPPORT	85
NEOPRENE KNEE SUPPORT WITH STABILIZED PATELLA	86
NEOPRENE KNEE SLEEVE	86
FOOT AND ANKLE	87
REBOUND® DIABETIC WALKER.....	88
REBOUND® AIR WALKER	90
AFO DYNAMIC	92

OA & INJURY SOLUTIONS CATALOGUE

AFO LIGHT.....	94
REBOUND FOOT-UP®	95
FOOT-UP®	96
FORM FIT® WALKER, AIR	97
FORM FIT® WALKER	98
RANGE OF MOTION WALKER.....	99
DH OFFLOADING WALKER®.....	100
DH OFFLOADING POST-OP STYLE SHOE	100
PEDIATRIC WALKER	101
METAL UPRIGHT AIR WALKER	102
REBOUND® ANKLE BRACE.....	103
FORM FIT® ANKLE BRACE	104
FORM FIT® ANKLE WITH SPEEDLACE	105
AIRFORM® UNIVERSAL INFLATABLE ANKLE STIRRUP	106
AIRFORM® INFLATABLE ANKLE STIRRUP	107
AIRFORM® PRE-INFLATED ANKLE STIRRUP	107
FORM FIT® HONEYCOMB ANKLE STIRRUP	108
FORM FIT® FOAM ANKLE BRACE	108
EXOFORM® DORSAL NIGHT SPLINT	109
FORM FIT® NIGHT SPLINT	110
AIRFORM® NIGHT SPLINT	111
AFO LEAF SPRING	111
MESH TOP POST-OP SHOE	112
SOFT TOP POST-OP SHOE.....	112
CANVAS ROCKER BOTTOM CAST SHOE	113
SPECIALITY CARE, CASTING & TRACTION	114
STOCKINETTE	115
TECHFORM® PREMIUM	115
TECHFORM® ORTHO-ROLL®	116
TECHFORM® PRE-CUT SPLINTS.....	116

ÖSSUR OA SOLUTIONS

Enhance quality of life and mobility for people living with osteoarthritis.

These product solutions have been built on decades of experience and are designed to enhance the quality of life of people living with osteoarthritis by reducing pain and improving mobility. Through innovation leadership, OA expertise and partnerships with medical professionals, we offer comprehensive, proven and effective treatment options.

UNLOADER®

The Össur Unloader® brand name has become synonymous with Osteoarthritis bracing, representing OA Solutions' technology as the industry standard, with more than 15 clinical studies demonstrating its efficacy.

Össur strives to maintain its position as the innovation leader for OA Solutions by continuing to develop simple, yet effective clinically proven treatment solutions that increase people's quality of life.



William McGirt
Professional Golfer

HIP OA



UNLOADER® HIP

PAGE 6

KNEE OA



UNLOADER ONE® LITE

PAGE 8



UNLOADER ONE®

PAGE 10



UNLOADER® CUSTOM

PAGE 16

UNLOADER® HIP

Unloader Hip is designed to reduce pain by optimizing load dispersion for patients suffering from mild and moderate osteoarthritis of the hip, thus contributing to hip stability and improved mobility.




INDICATIONS FOR USE

- Mild and moderate osteoarthritis of the hip
- Conditions benefitting from repositioning of the hip and/or compression around the greater trochanters

FEATURES AND BENEFITS

- The Dynamic Rotation Straps with SmartDosing dial allows individual adjustment of the external rotation force, to achieve optimum pain relief for each patient
- Easy donning and doffing thanks to pocket wrap closure system and zipper
- Pulley system for improved stability
- Premium quality lycra for breathability and a snug fit



REQUIRED MEASUREMENTS

- Measure the hip circumference at the anterior superior iliac spine



UNLOADER® HIP

ORDERING INFORMATION

UNLOADER HIP

Size	Hip Circumference	Woman	Men
Small	31.5-35" (80-88cm)	B-465350030	B-465350031
Medium	35-38" (88-96cm)	B-465350032	B-465350033
Large	38-41" (96-104cm)	B-465350034	B-465350035
X-Large	41-45" (104-114cm)	B-465350036	B-465350037
XX-Large	45-49" (114-124cm)	B-465350038	B-465350039

ACCESSORIES/ OPTIONS

REPLACEMENT TROUSERS

Size	Measurements	Woman	Men
Small	31.5-35" (80-88cm)	B-465350077	B-465350078
Medium	35-38" (88-96cm)	B-465350079	B-465350080
Large	38-41" (96-104cm)	B-465350081	B-465350082
X-Large	41-45" (104-114cm)	B-465350083	B-465350084
XX-Large	45-49" (114-124cm)	B-465350085	B-465350086

OTHER ACCESSORIES

Description	Item No
Rotation strap	B-465350064
Washbag	B-465350065
Trochanter pad	B-465350066
Thigh strap S	B-465350067
Thigh strap M	B-465350068
Thigh strap L	B-465350069
Thigh strap XL	B-465350070
Thigh strap XXL	B-465350071
Bilateral kit* S	B-465350072
Bilateral kit* M	B-465350073
Bilateral kit* L	B-465350074
Bilateral kit* XL	B-465350075
Bilateral kit* XXL	B-465350076

*The bilateral kit consists of a rotation control strap, a thigh strap and a trochanter pad

UNLOADER ONE® LITE

Unloader One Lite is a slim and low profile knee brace that relieves pain from knee osteoarthritis and degenerative meniscal tears by means of the clinically proven Dynamic Force System that provides uni-compartmental unloading of the knee joint.

With an integrated Lycra® sleeve and shells with overmold, the brace is an optimal combination of comfort and sleek look.

INDICATIONS FOR USE

- Mild to moderate unicompartmental osteoarthritis
- Degenerative meniscus
- Other knee conditions that may benefit from unicompartmental unloading

FEATURES AND BENEFITS

- 3-Point Leverage System provides clinically relevant pain relief of affected compartment
- Single SmartDosing dial allows fine-tuned pain management
- Helps keep patients mobile during daily activities
- Discreet hinge and integrated sleeve ensure a sleek look on leg and snug fit underneath clothing
- Lightweight and easy to use

REQUIRED MEASUREMENTS

- Circumference measurement: 6"/15cm above and below mid-patella



UNLOADER ONE® LITE

ORDERING INFORMATION

UNLOADER ONE LITE

Size	Circumference 6" (15cm) above mid patella	Circumference 6" (15cm) below mid-patella	Left medial	Right medial	Left lateral	Right lateral
X-Small	12-15" (30-38cm)	10-12" (25-30cm)	UOL23000	UOL13000	UOL24000	UOL14000
Small	15-17.75" (38-45cm)	12-13.75" (30-35cm)	UOL23010	UOL13010	UOL24010	UOL14010
Medium	17.75-22" (45-56cm)	13.75-16.5" (35-42cm)	UOL23020	UOL13020	UOL24020	UOL14020
Large	22-26" (56-66cm)	16.5-20" (42-51cm)	UOL23030	UOL13030	UOL24030	UOL14030
X-Large	26-29" (66-73cm)	20-24" (51-61cm)	UOL23040	UOL13040	UOL24040	UOL14040

Hinge is placed on the side of the affected compartment

ACCESSORIES

Snap-On DFS Strap Kit	Item No
Unloader One Lite DFS Kit X-Small	UNF01302
Unloader One Lite DFS Kit Small	UNF01304
Unloader One Lite DFS Kit Medium	UNF01306
Unloader One Lite DFS Kit Large	UNF01308
Unloader One Lite DFS Kit X-Large	UNF01310

Includes DFS kit, polybag with zip lock 8x10in/20x25cm, Label, QRC

Calf Strap Kit	Item No
Unloader One Lite calf strap kit X-Small	UNF01402
Unloader One Lite calf strap kit Small	UNF01404
Unloader One Lite calf strap kit Medium	UNF01406
Unloader One Lite calf strap kit Large	UNF01408
Unloader One Lite calf strap kit X-Large	UNF01410

Includes calf strap kit, poly bag with zip lock 8x10in/20x25cm, Label GPO Thermal Ossur

Alligator DFS Strap Kit	Item No
Unloader One Lite alligator DFS kit	UNF01502

Includes alligator DFS kit (one size), poly bag with zip lock 8x10in/20x25cm, Label, QRC

UNLOADER ONE®

The Unloader One is a comfortable, lightweight, low-profile brace that provides excellent suspension and pain relief. This lightweight brace preserves the clinically proven Össur Unloader 3-Points of Leverage design by incorporating Dual Dynamic Force Straps that are easily fine-tuned with the SmartDosing® dial. Color-coded Quick Fit Buckles simplify application and removal, while flexible shells and breathable tibial liners, coated with Össur Sensil® Silicone, improve comfort and reduce migration, helping to ensure patient compliance. Recent research demonstrates that the Unloader One reduces pain, improves function and decreases the use of pain medication.



INDICATIONS FOR USE

- Moderate to severe unicompartmental osteoarthritis

FEATURES AND BENEFITS

- Dual Dynamic Force Straps designed to reduce pressure on sensitive bony surfaces of the knee, resulting in reduction of pain
- Tool-free SmartDosing dial enables fine-tuned pain management and can be adjusted even when wearing the brace under clothing
- Color coded quick fit buckles enable easy application and removal
- Proprietary Sensil® silicone liners designed to reduce migration
- Flexible, low-profile upright contours to patient anatomy and allows use underneath clothing
- Quick Connect Strapping for tool-free removal and replacement
- Lightweight frame is perforated for flexibility and ventilation with custom paint options for the Custom version
- Unloader One is available in pre-fabricated and Custom version

REQUIRED MEASUREMENTS

- Circumference measurement: 6" (15cm) below mid-patella

CUSTOM STANDARD COLORS



CUSTOM PAINT



UNLOADER ONE®

ORDERING INFORMATION

UNLOADER ONE WITH SMARTDOSING

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2-14" (31-36cm)	B-240618712	B-240518712	B-240628712	B-240528712
Medium	14-16.75" (36-43cm)	B-240618713	B-240518713	B-240628713	B-240528713
Large	16.75-20" (43-51cm)	B-240618714	B-240518714	B-240628714	B-240528714
X-Large	20-24" (51-61cm)	B-240618715	B-240518715	B-240628715	B-240528715

UNLOADER ONE SHORT WITH SMARTDOSING

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	9.8-12.2" (25-31cm)	B-241618711	B-241518711	B-241628711	B-241528711
Small	12.2-14" (31-36cm)	B-241618712	B-241518712	B-241628712	B-241528712
Medium	14-16.75" (36-43cm)	B-241618713	B-241518713	B-241628713	B-241528713
Large	16.75-20" (43-51cm)	B-241618714	B-241518714	B-241628714	B-241528714
X-Large	20-24" (51-61cm)	B-241618715	B-241518715	B-241628715	B-241528715

UNLOADER ONE WITH RATCHET

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2-14" (31-36cm)	B-240619712	B-240519712	B-240629712	B-240529712
Medium	14-16.75" (36-43cm)	B-240619713	B-240519713	B-240629713	B-240529713
Large	16.75-20" (43-51cm)	B-240619714	B-240519714	B-240629714	B-240529714
X-Large	20-24" (51-61cm)	B-240619715	B-240519715	B-240629715	B-240529715

UNLOADER ONE SHORT WITH RATCHET

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	9.8-12.2" (25-31cm)	B-241619711	B-241519711	B-241629711	B-241529711
Small	12.2-14" (31-36cm)	B-241619712	B-241519712	B-241629712	B-241529712
Medium	14-16.75" (36-43cm)	B-241619713	B-241519713	B-241629713	B-241529713
Large	16.75-20" (43-51cm)	B-241619714	B-241519714	B-241629714	B-241529714
X-Large	20-24" (51-61cm)	B-241619715	B-241519715	B-241629715	B-241529715

Standard version is recommended for patients between 5'5" and 6'4" (1,65 and 1,92 meters tall).
Short version is recommended for patients under 5'5" (1,65 meters tall).
Hinge is placed on the side of the affected compartment

ACCESSORIES

ACCESSORY KITS

Description	SMARTDOSING (STANDARD VERSION)		SMARTDOSING (SHORT VERSION)	
	RM/LL	RL/LM	RM/LL	RL/LM
Strap & Strap Pad & Buckle Clasp Kit	B-705277003	B-705278003	B-705277004	B-705277005
Strap & Strap Pad Kit	B-705270003		B-705277006	
Strap Kit	B-701270013		B-705277007	

UNLOADER ONE®

ACCESSORIES

ACCESSORY KITS

Description	Right medial/Left lateral	Right lateral/Left medial
Buckle Clasp Kit	B-708217003	B-708218003
Alligator Clip Kit (included in all strap kits)	B-704210003	
Extension Stop Kit	B-301040046	
Velcro Twist-Fit Attachment Kit	B-721240003	

ÖSSUR KNEE UNDERSLEEVE

Size	Thigh Circumference	Item No
X-Small	12-15" (30.5-38cm)	B-250150201
Small	15-18" (38-45.5cm)	B-250150202
Medium	18-22" (45.5-56cm)	B-250150203
Large	22-26" (55-66cm)	B-250150204
X-Large	26-29" (66-73.5cm)	B-250150205
XX-Large	29-33" (73.5-84cm)	B-250150206

When using an Össur Knee Undersleeve we advise using the suspension strap kit to minimize migration

SUSPENSION STRAP KIT

Description	Item No
Suspension Strap Kit	UNL00212

SENSIL LINERS

Size	Calf		Thigh	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-419700311	B-419800311	B-409700321	B-409800321
Small	B-419700312	B-419800312	B-409700322	B-409800322
Medium	B-419700313	B-419800313	B-409700323	B-409800323
Large	B-419700314	B-419800314	B-409700324	B-409800324
X-Large	B-419700315	B-419800315	B-409700325	B-409800325

DOESKIN LINERS

Size	Calf		Thigh	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-424703721	B-424803721	B-412701721	B-412801721
Small	B-424703722	B-424803722	B-412701722	B-412801722
Medium	B-424703723	B-424803723	B-412701723	B-412801723
Large	B-424703724	B-424803724	B-412701724	B-412801724
X-Large	B-424703725	B-424803725	B-412701725	B-412801725

When using a doeskin calf liner we advise using a suspension strap kit to minimize migration

SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

Size	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-715217031	B-715218031
Small	B-715217032	B-715218032
Medium	B-715217033	B-715218033
Large	B-715217034	B-715218034
X-Large	B-715217035	B-715218035

UNLOADER ONE® CUSTOM



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the Unloader One. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



DIGITAL MEASUREMENT SYSTEM (DMS)

This system requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for an Unloader One Custom without the need for training.

ORDERING INFORMATION

UNLOADER ONE CUSTOM WITH SMARTDOSING

Color	Thigh Circumference	Calf Circumference	Left	Right
Black	34.3-73.5cm	25-61.5cm	B-115600112	B-115500112
Grey	34.3-73.5cm	25-61.5cm	B-115601712	B-115501712

ACCESSORIES

ÖSSUR KNEE UNDERSLEEVE

Size	Circumference	Grey	Black
X-Small	12-15" (30.5-38cm)	B-250150201	B-250150101
Small	15-18" (38-45.5cm)	B-250150202	B-250150102
Medium	18-22" (45.5-56cm)	B-250150203	B-250150103
Large	22-26" (56-66cm)	B-250150204	B-250150104
X-Large	26-29" (66-73.5cm)	B-250150205	B-250150105
XX-Large	29-33" (73.5-84cm)	B-250150206	B-250150106

When using an Össur Knee Undersleeve we advise using the suspension strap kit to minimize migration

SENSIL LINERS

Size	Calf		Thigh	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-419700311	B-419800311	B-409700321	B-409800321
Small	B-419700312	B-419800312	B-409700322	B-409800322
Medium	B-419700313	B-419800313	B-409700323	B-409800323
Large	B-419700314	B-419800314	B-409700324	B-409800324
X-Large	B-419700315	B-419800315	B-409700325	B-409800325

UNLOADER ONE® CUSTOM

DOESKIN LINERS

Size	Calf		Thigh	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-424703721	B-424803721	B-412701721	B-412801721
Small	B-424703722	B-424803722	B-412701722	B-412801722
Medium	B-424703723	B-424803723	B-412701723	B-412801723
Large	B-424703724	B-424803724	B-412701724	B-412801724
X-Large	B-424703725	B-424803725	B-412701725	B-412801725

When using a doeskin calf liner we advise using a suspension strap kit to minimize migration. For replacement liners, please contact Customer Service.

SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

Size	Grey		Black	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-715217031	B-715218031	B-715217011	B-715218011
Small	B-715217032	B-715218032	B-715217012	B-715218012
Medium	B-715217033	B-715218033	B-715217013	B-715218013
Large	B-715217034	B-715218034	B-715217014	B-715218014
X-Large	B-715217035	B-715218035	B-715217015	B-715218015

SUSPENSION STRAP KIT

Description	Grey	Black
Suspension Strap Kit	UNL00212	UNL00211

SMARTDOSING ACCESSORY (STANDARD VERSION)

Description	Grey RM/LL	Grey RL/LM	Black RM/LL	Black RL/LM
SmartDosing Strap & Strap Pad & Buckle Clasp Kit	B-705277003	B-705278003	B-705277001	B-705278001
Strap & Strap Pad Kit	B-705270003		B-705270001	
Strap Kit	B-701270013		B-701270011	
SmartDosing Complete Strap Set	B-705270003			

SMARTDOSING ACCESSORY KITS (SHORT VERSION)

Description	Grey RM/LL	Grey RL/LM
Short SmartDosing Strap & Strap Pad & Buckle Clasp Kit	B-705277004	B-705277005
Short Strap & Strap Pad Kit	B-705277006	
Short Strap Kit	B-705277007	

ACCESSORY KITS

Description	Grey RM/LL	Grey RL/LM	Black RM/LL	Black RL/LM
Buckle Clasp Kit	B-708217003	B-708218003	B-708217001	B-708218001
Alligator Clip Kit (included in all strap kits)	B-704210003		B-704210001	
Extension Stop Kit	B-301040046		B-301040047	
Velcro Twist-Fit Attachment Kit	B-721240003		B-721240001	

UNLOADER ONE® PLUS

ORDERING INFORMATION

UNLOADER ONE PLUS WITH RATCHET

Size	Thigh Circumference	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Large	21-25" (53-64cm)	14-16.75" (36-43cm)	B-245610004	B-245510004	B-245620004	B-245520004
X-Large	25-29" (64-74cm)	16.75-20" (43-51cm)	B-245610005	B-245510005	B-245620005	B-245520005

UNLOADER ONE PLUS OTS, LONG VERSION

Size	Thigh Circumference	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Large	21-25" (53-64cm)	14-16.75" (36-43cm)	B-244610004	B-244510004	B-244620004	B-244520004
X-Large	25-29" (64-74cm)	16.75-20" (43-51cm)	B-244610005	B-244510005	B-244620005	B-244520005

ACCESSORIES

SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

Size	Right medial/Left lateral	Right lateral/Left medial
Medium	B-715237033	B-715238033
Large	B-715237034	B-715238034
X-large	B-715237035	B-715238035

LINERS

Size	Calf, sensil		Thigh, doeskin	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
Medium	B-413700313	B-413800313	B-414700323	B-414800323
Large	B-413700314	B-413800314	B-414700324	B-414800324
X-Large	B-413700314	B-413800315	B-414700325	B-414800325
XX-Large			B-414700326	B-414800326

ACCESSORY KITS

Description	Right medial/Left lateral	Right lateral/Left medial
Strap & Strap Pad & Buckle Clasp Kit	B-701237008	B-701238008
Strap & Strap Pad Kit	B-701230006	
Strap Kit	B-701230004	
Pad Kit	B-702230003	
Extra Pads	B-702230013	
Buckle Clasp Kit	B-708237003	B-708238003
Alligator Clip Kit (included in all strap kits)	B-704210003	
Gastroc Kit	B-705230001	
Elastic Thigh Strap 8"	B-701237083	B-701238083
Elastic Thigh Strap 10"	B-701237103	B-701238103

UNLOADER® CUSTOM

With a semi-rigid frame and an adjustable hinge, the Unloader Custom accommodates a wide variety of leg shapes and sizes. The Adjustable Dynamic Joint enables changes to the hinge and additional leverage for optimal pain management.



INDICATIONS FOR USE

- Mild to severe unicompartmental osteoarthritis

FEATURES AND BENEFITS

- Custom molded semi-flexible shells are durable, adaptable to leg shape changes, yet rigid enough to contain tissue while force is applied at the knee
- Patented Adjustable Dynamic Joint™ (ADJ) high-carbon steel hinge is strong enough to unload the pressure from the compartment of the knee. Adjustments are made easily
- The Dynamic Force Strap® designed to reduce pressure on sensitive bony surfaces of the knee, resulting in reduction of pain
- Available in standard and in SmartDosing version. Tool-free SmartDosing dial option enables fine-tuned pain management and can be adjusted even when wearing the brace under clothing
- SmartDosing version available with Sensil® silicone calf liner to reduce migration
- Available in 12 standard colors in matte or glossy finish or custom painted


REQUIRED MEASUREMENTS

- Custom-to-measurement

STANDARD COLORS



CUSTOM PAINT



UNLOADER® CUSTOM



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the Unloader One. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

ORDERING INFORMATION

UNLOADER CUSTOM

Description	Item No
Right	B-122500020
Left	B-122600020

ACCESSORIES

Description	Item No
Knee Shield	B-250110000
Hinge cover	B-251030000
Water Ski Strap Lock	B-303010084
Distal Non-Elastic Thigh Strap	B-507010000
Proximal Non-Elastic Thigh Strap	B-507010000
AMS Wrap (custom made)	B-426000150
AMS Sleeve (custom made)	B-426000250

ELASTIC CALF SUSPENSION STRAP

Size	Item No	Size	Item No	Size	Item No
14"	B-504020114	18"	B-504020018	22"	B-504020022
15"	B-504020115	19"	B-504020019	23"	B-504020023
16"	B-504020116	20"	B-504020020	24"	B-504020024
17"	B-504020017	21"	B-504020021		

SLEEVES

Size	Circumference	Sofsleeve	Neoprene undersleeve	Neoprene oversleeve
Small	< 15" (38cm)	20643	20251	20051
Medium	15-18" (38-46cm)	20644	20252	20052
Large	18-19.5" (46-50cm)	20644	20253	20053
X-Large	19.5-21" (50-55cm)	20645	20254	20054
XX-Large	21.5-24.5" (55-62cm)	20645	20255	20055
XXX-Large	24.5-28" (62-71cm)	20646	20256	20056
XXXX-Large	28-34" (71-86cm)	20647	20257	20057

Size	Measuring	Össur knee undersleeve	Size	Measuring	Össur knee undersleeve
X-Small	12-15" (30.5-38cm)	B-250150101	Large	22-26" (56-66cm)	B-250150104
Small	15-18" (38-45.5cm)	B-250150102	X-Large	26-29" (66-73.5cm)	B-250150105
Medium	18-22" (45.5-56cm)	B-250150103	XX-Large	29-33" (73.5-84cm)	B-250150106

CTI® OA

CTi® OA's composite shells can be constructed with up to 7° of varus or valgus correction, offering a high level of protection and stability for clinicians interested in treating osteoarthritis and ligament instability equally. Key features include a Total Support System, Accutrac® hinges that mimic the natural motion of the knee joint and breathable liners, coated with Össur® Sensil® Silicone to improve efficacy and virtually eliminate migration.



INDICATIONS FOR USE

- Mild to severe unicompartmental osteoarthritis
- ACL, MCL, LCL, PCL, rotary and combined instabilities
- Demanding activities of daily living

FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with extension stops
- Proprietary Sensil® silicone liners designed to reduce migration
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Carbon composite, rigid frame construction
- Custom-fabricated to each patient's unique measurements with up to 7° varus/valgus correction built into the frame
- Non-corrosive materials, ideal for water sports
- Standard model for medium to high impact levels. Vapor model for low to medium impact levels. Pro Sport model for highest impact levels

REQUIRED MEASUREMENTS

- Custom-made to measure



CTI® OA

ORDERING INFORMATION

CTI OA

Description	Left	Right
Standard	B-116630010	B-116530010
Pro Sport	B-117630010	B-117530010
Vapor	B-118630010	B-118530010

For accessories see CTi Custom.



CTI® MEASURING DEVICE

Simplified measuring technique saves fitting time. The need for costly casting materials is eliminated. Provides coordinates of patient's leg from which a three-dimensional chart is produced for brace manufacturing. Designed exclusively for use with the CTi OA knee brace. The CTi measuring device is provided free of charge with the initial CTi brace order.

Description	Item No
CTi Custom Measuring Device	42711
Constant Force Caliper	B-312000000



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the CTi OA. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.




DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for an CTi OA without the need for training.

STANDARD COLORS



CUSTOM PAINT



FORM FIT® OA WRAPAROUND

The Form Fit OA Wraparound is a cost effective, low-profile option for patients with mild to moderate unicompartmental knee osteoarthritis. The Dynamic Force Strap tensions as the leg extends, providing pain relief in a single pull motion. The slim polycentric hinge design and breathable material offers a discreet and comfortable fit. Dot-coded color wrap closures and non-slippery rubber tabs allows for a simple donning procedure for the user.



INDICATIONS FOR USE

- Mild to moderate unicompartmental knee osteoarthritis

FEATURES AND BENEFITS

- Easy wraparound design
- Dot coded closure wraps
- Non-slippery rubber tab
- Dynamic Force Strap
- Low-profile design & slim hinge
- Breathable lycra

REQUIRED MEASUREMENTS

- Thigh Circumference: 15cm (6") above mid-patella
- Calf Circumference: 15cm (6") below mid-patella

FORM FIT® OA WRAPAROUND

ORDERING INFORMATION

FORM FIT OA WRAPAROUND

Size	Thigh Circumference	Calf Circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	14-16" (36-41cm)	11-12.5" (28-32cm)	B-254901008	B-254901001	B-254901022	B-254901015
Small	16-18" (41-46cm)	12.5-14" (32-36cm)	B-254901009	B-254901002	B-254901023	B-254901016
Medium	18-21" (46-53cm)	14-16" (36-41cm)	B-254901010	B-254901003	B-254901024	B-254901017
Large	21-23.5" (53-60cm)	16-18" (41-46cm)	B-254901011	B-254901004	B-254901025	B-254901018
X-Large	23.5-26.5" (60-67cm)	18-20" (46-51cm)	B-254901012	B-254901005	B-254901026	B-254901019
2X-Large	26.5-29.5" (67-75cm)	20-22" (51-56cm)	B-254901013	B-254901006	B-254901027	B-254901020
3X-Large	29.5-32" (75-81cm)	22-24" (56-61cm)	B-254901014	B-254901007	B-254901028	B-254901021

ACCESSORIES

SUSPENSION STRAP

Description	Item No
Suspension Strap	OAW00007

ÖSSUR INJURY SOLUTIONS

For people recovering from fractures, ligament injuries or in need of post-operative treatment solutions.


These product solutions have been built on decades of experience and are designed to enhance the healing process of bone and soft tissue injuries. Through innovation leadership, proven biomechanical approaches, and partnerships with leading medical professionals, we strive to offer comprehensive and effective treatment options.

Össur Injury Solutions are globally recognized, primarily under the brand names CTi®, Rebound® and Miami®, all of which are developed to assist healing and increase people's quality of life.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimize healing while maintaining function and mobility. These Functional Healing products incorporate Össur proprietary technology and are designed to enhance clinically relevant outcomes.



FUNCTIONAL HEALING-
WHERE MOBILITY MEANS RECOVERY.



SPINE



MIAMI LSO™

PAGE 24



MIAMI J®

PAGE 29



FUNCTIONAL HEALING



REBOUND® CARTILAGE

PAGE 60



REBOUND® PCL

PAGE 63

KNEES



CTi®

PAGE 65



REBOUND® DUAL

PAGE 71



REBOUND® POST-OP KNEE

PAGE 78

FEET



REBOUND® AIR WALKER

PAGE 90

MIAMI LSO™

Miami LSO is a spinal orthosis designed for protection of injured ligaments or muscles and post-surgical immobilization. It provides low-profile, anatomical immobilization, easy fit and wear, options for modular support and an advanced pulley system with a length-adjustable pulley cord.



INDICATIONS FOR USE

Conditions requiring gross immobilization of the trunk in the lumbar region. Such conditions may include:

- Stable, non-displaced spinal fractures (L1-L5)
- Spinal stenosis
- Herniated discs
- Degenerative spinal pathologies
- Spondylolisthesis
- Spondylolysis

Contra-indication: Unstable, displaced spinal fractures

FEATURES AND BENEFITS

- Slim, breathable material
- Proven Miami double-pulley system with length-adjustable pulley cord
- Intuitive fitting without trimming with fold-over belt arms & Precise-Size™ landmarks
- Optional fit tool for perfect fitting directly on the patient
- Universal belt sizing accommodates waists up to 178cm/70"
- Contoured posterior panels sized S-XL for a precise, anatomical fit
- Quick-Fit pocket enclosure and blue highlights for patient ease of use

REQUIRED MEASUREMENTS


- Waist circumference







MIAMI LSO™

ORDERING INFORMATION




OPTION 1: Order braces for inventory

MIAMI LSO BELT		
		
SIZE	WAIST CIRC.	PART #
6" Universal (Slim)	26-50" (66-127cm)	MLSO-1
8" Universal (Standard)	26-50" (66-127cm)	MLSO-2
6" (Slim) Belt Extension	26-70" (66-177cm)	MLSO-1EXT
8" (Standard) Belt Extension	26-70" (66-177cm)	MLSO-2EXT



CONTOURED POSTERIOR PANEL			CONTOURED POSTERIOR PANEL WITH LATERAL SUPPORTS		
					
SIZE	WAIST CIRC.	PART #	SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86cm)	MLSO-P11	Small (S)	26-34" (66-86cm)	MLSO-P12
Medium (M)	34-38" (86-96cm)	MLSO-P21	Medium (M)	34-38" (86-96cm)	MLSO-P22
Large (L)	38-50" (96-127cm)	MLSO-P31	Large (L)	38-50" (96-127cm)	MLSO-P32
Extra Large (XL)	50-66" (127-167cm)	MLSO-P41	Extra Large (XL)	50-66" (127-167cm)	MLSO-P42

STANDARD POSTERIOR PANEL		ANTERIOR PANEL	
			
SIZE	PART #	SIZE	PART #
Standard	MLSO-P56	6" Panel (Slim)	MLSO-PA6
		8" Panel (Standard)	MLSO-PA8

OPTION 2: Order a brace in a kit WITH ANTERIOR PANEL

KIT CONFIGURATION	POSTERIOR PANEL SIZE	PART # SLIM WIDTH (6")	PART # STANDARD WIDTH (8")
	Standard	MLSO-156-1	MLSO-256-1
	SM	MLSO-111-1	MLSO-211-1
	MD	MLSO-121-1	MLSO-221-1
	LG	MLSO-131-1	MLSO-231-1
	XL	MLSO-141-1	MLSO-241-1
	SM	MLSO-112-1	MLSO-212-1
	MD	MLSO-122-1	MLSO-222-1
	LG	MLSO-132-1	MLSO-232-1
	XL	MLSO-142-1	MLSO-242-1

WITHOUT ANTERIOR PANEL

KIT CONFIGURATION	POSTERIOR PANEL SIZE	PART # SLIM WIDTH (6")	PART # STANDARD WIDTH (8")	ÖSSUR FIT TOOL	PART #
	Standard	MLSO-156-2	MLSO-256-2		B-04010000

Note: The standard posterior panel is always recommended when no other posterior panel is being used.

MIAMI LUMBAR® POSTEO

Designed with patient comfort and compliance in mind, the Miami Lumbar Posteo provides functional relief of vertebral fracture pain, and facilitates return to mobility. The lumbar compression system offers both pain relief via active vertebral offloading and better grip of patient's hips for greater ease of use and compliance. The Miami Lumbar Posteo includes an innovative approach to ease of use and consistent donning.



INDICATIONS FOR USE

- Back pain caused by stable vertebral compression fractures (VCF)

FEATURES AND BENEFITS

- Active offloading via thoracic extension and lumbar compression
- Integrated pulley system for customized fit
- Trimmable axillary straps with adjustable tension for optimal stabilization
- Cold-moldable aluminum PTE to fit different anatomies
- Extra cushioning around the aluminum strut for minimized skin pressure

REQUIRED MEASUREMENTS

- Back length measured from the top of the sacrum up to C7 vertebrae, following the spinal curvature

Please note that the Miami Lumbar Posteo includes a universal lumbar belt which can be trimmed down to the appropriate size. Calculation of the appropriate back length size, measured from the tip of the sacrum to the C7 vertebra, following the curvature of the spine.

ORDERING INFORMATION

MIAMI LUMBAR POSTEO

Size	Back Length	Item No
X-Small	14.5"-17" (37-44cm)	LMB-POSTEO-XS
Small	17"-20" (44-51cm)	LMB-POSTEO-SM
Medium	19"-22" (49-56cm)	LMB-POSTEO-MD
Large	21"-24" (55-61cm)	LMB-POSTEO-LG

ACCESSORIES

Description	Item No
Strap Pad	LMB-POSTEO-STPD



MIAMI LUMBAR TLSO®

Miami Lumbar TLSO is an easy-to-use lumbar orthosis for post-surgical immobilization,. The superior anatomical design follows the contours of the spine while ensuring patient comfort. The modular system enables step-down treatment options for every stage of the recovery process. A TLSO attachment can be used to provide gross immobilization of the thoracic and lumbar spines.



INDICATIONS FOR USE

Miami Lumbar belt with added rigid panels and TLSO attachment:

- Post-surgical immobilization
- Herniated disc
- Spondylolisthesis
- Spondylolysis
- Degenerative spinal fractures

Contraindications: Unstable, displaced fractures

FEATURES AND BENEFITS

- Anatomical design follows the contours of the spine perfectly, ensuring better immobilization and patient comfort
- Designed for maximum breathability and patient compliance by incorporating aerospace linings and ventilated panels
- Modular system enables step-up, step-down treatment options
- All belts can be trimmed down to the appropriate size allowing for reduced inventory levels
- Compression system eases application
- Optional air bladder for increased customization

REQUIRED MEASUREMENTS

- Waist circumference

ORDERING INFORMATION

ACCESSORIES

Item	Item No
TLSO system with ATE & PTE	LMB-T01
PTE and shoulder straps	LMB-T02
Air Bladder	LMB-A01
Pad replacement set	LMB-A02

KITS

Description	Item No
Standard Miami Lumbar TLSO Kit	LMB-XXYYNPT
Slim Miami Lumbar TLSO Kit	LMBSL-0XXYYNPT

XX = 01 (XS); 02 (S); 03 (M); 04 (L); 05 (XL); 06 (2XL); 07 (3XL); 08 (4XL); 09 (5XL); 10 (6XL)
YY = NL (Without Lateral Supports); LS (With Lateral Supports)

MIAMI J® ADVANCED

The Miami J Advanced is based on the same proven principles as the Miami J cervical collar, widely acknowledged to be the best c-spine immobilizer available. Helping to provide immobilization of the injured area, it also minimizes skin breakdown and patient discomfort, thus resulting in optimal patient compliance. This unique adjustable collar fits various phenotype sizes, ensuring an accurate fit while enabling reduced stock.




INDICATIONS FOR USE

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc
- Post-operative
- Motor Neurone Disease

FEATURES AND BENEFITS

- Unique sizing mechanism allows for easy modification to suit individual patients, simultaneously achieving an accurate fit and reduced stock levels
- Bioengineered to minimize pressure points in key known areas: Chin, occiput, trapezius and clavicle
- The shape of the front panel avoids pressing on the cervical veins
- Sternal pad disperses ambulatory pressure and is fully removable, enabling the front of the collar to be taken off when cleaning intubated patients, without the need to remove the tube
- Tamper resistant quick-lock buttons prevent unintended height adjustment
- Patented Flex-Edge® technology delivers soft over-molded edges that gently conform to the contours of the anatomy of the back of the skull. The chin insert on the front of the collar ensures even pressure distribution MRI, CT and X-ray lucent



REQUIRED MEASUREMENTS

- Adjustable to Miami J phenotype sizes, except Super Short and XS

ORDERING INFORMATION

MIAMI J ADVANCED

Description	Collar set*	Collar
Miami J Advanced	MJAR-101	MJA-101

*The collar set comes with one additional set of replacement pads.

ACCESSORIES/OPTIONS

REPLACEMENT PADS

Description	Item No
Replacement Pads	MJAP-101
Occian Back	ACB-100

MIAMI J®

Known as the superior c-spine immobilizer through multiple independent studies, the Miami J's patented design meets the specialized needs of all those afflicted with cervical injuries. The Miami J offers features that enhance compliance while providing immobilization, inhibition of skin breakdown and patient comfort.



INDICATIONS FOR USE

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc
- Post-operative
- Motor Neurone Disease

FEATURES AND BENEFITS

- Sorbatex™ padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- Bioengineered to minimize pressure points in key known areas: Chin, occiput, trapezius and clavicle
- Angulation buttons for occiput and sternal pad is removable when supine, lying and swallowing
- MRI, CT and X-ray lucent, MR safe

REQUIRED MEASUREMENTS

- Sizing based on 6 phenotypes to ensure correct fit and optimum immobilization

Phenotype-driven sizing system makes size selection simple while drastically reducing margin of error in fitting a patient. Each collar has a bioengineered front with a specific back designed to work together for superior immobilization. For paediatric sizing, please see Miami Jr in this Catalog.

Specialty Sizes (40% of population)

Standard Sizes (60%)

200S Kyphotic (10%)	250 XSmall (5%)	200L 200B Stout (15%)	500 Tall (10%)	400 Regular (30%)	300 Short (30%)
<ul style="list-style-type: none"> • kyphotic • chin-on-chest • geriatric • osteoporitic and/or ankylosing spondylitic 	<ul style="list-style-type: none"> • Very short, thin neck • Principally female • More prevalent in the Asian population 	<ul style="list-style-type: none"> • Very large neck circumference • Obese • "No neck" • Massive shoulders 	<ul style="list-style-type: none"> • Long, tall neck • "Swan neck" • Young Women • Adolescents 	<ul style="list-style-type: none"> • Female • Standard size for mature women • Thin, mature male 	<ul style="list-style-type: none"> • Standard size for mature men • Short-necked women

*For individuals with neck circumference > 25" (63.5cm), the Bariatric collar is recommended.

MIAMI J®

ORDERING INFORMATION

MIAMI J

Size	Collar set*	Collar	Collar Back	Collar Front
Super Short	MJR-200S	MJ-200S	MJ-400BK	MJ-200SFRT
Stout	MJR-200L	MJ-200L	MJ-200LBK	MJ-200LFRT
XS	MJR-250	MJ-250	MJ-250BK	MJ-250FRT
Short	MJR-300	MJ-300	MJ-300BK	MJ-300FRT
Regular	MJR-400	MJ-400	MJ-400BK	MJ-400FRT
Tall	MJR-500	MJ-500	MJ-500BK	MJ-500FRT
Bariatric	MJR-200B	MJ-200B		

*Collar sets come with one additional set of replacement pads.
Miami J Front and Back Panels may be interchanged to accommodate patient anatomy, e.g. a Short Front may be used with a Tall Back.
Choose collar size based on patient neck height, not patient body size. Super Short collars are packaged with Regular Back Panels

ACCESSORIES

REPLACEMENT PADS

Size	Item No
Super Short	MJP-100
Stout	MJP-200L
XS	MJP-250
Short	MJP-100
Regular	MJP-100
Tall	MJP-100
Bariatric	MJP-200B

MIAMI J PATIENT PACK

Description	Item No
Miami J Patient Pack	MJ-PTPK2

MIAMI J STERNAL PAD

Description	Item No
Miami J Sternal Pad	MJ-SPK

MIAMI JTO® THORACIC EXTENSION

The Miami JTO Thoracic Extension is designed in conjunction with leading spine surgeons, orthotists and nurses. The orthosis is used exclusively with all Miami J® collars when even greater stability at the cervical junctions is indicated.



INDICATIONS FOR USE

- | | |
|--|--|
| <p>Lower C-Spine disruptions</p> <ul style="list-style-type: none"> • Level one instabilities • Isolated compartment fracture • Lateral mass fracture • Spinous Process fracture | <p>Upper C-Spine disruptions</p> <ul style="list-style-type: none"> • For increased stability over the collar alone • Stable Jefferson or Hangman's Fracture • Type I Odontoid Fracture |
|--|--|

FEATURES AND BENEFITS

- Provides more control than a collar alone can provide
- Innovative device designed to treat cervical and high-thoracic injuries
- Patented design with 'floating back' permits maximum immobilization without the risk of compensatory motion
- Snap on and off without tools
- Sorbatex™ II padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- MRI, CT and X-ray lucent, MR compatible

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

MIAMI JTO THORACIC EXTENSION

Description	Item No
Miami JTO	MT-900
Miami JTO with replacement softgoods	MTR-900

Compatible with all sizes of the Miami J® and the P3 Miami Jr®. Miami J collar is not included with Miami JTO Thoracic Extension. Miami J collar need to be purchased also

ACCESSORIES

Description	Item No
Replacement softgoods	MTP-950

OCCIAN® BACK & OCCIAN® BACK JR

The Occian Back is designed exclusively for the specialised needs of the multi-system trauma patient in the ICU. Universally sized, specifically engineered and designed with advanced, Intuitech™ pressure-relieving memory foam, the Occian Back helps caregivers manage c-spine patients who are at high risk of occipital breakdown.



INDICATIONS FOR USE

- Patients not awake and alert within 24-48 hours after admission
- Patients with pre-existing skin injuries and/or lacerations to the head/neck
- Burn patients who must wear a collar
- Elective surgery patients desiring night time post-operative comfort
- Patients who remain in a supine position for long periods of time related to respiratory, neurological, or haemodynamic instability

FEATURES AND BENEFITS

- Intuitech pressure relieving memory foam virtually eliminates occipital breakdown
- Maintains C spine stability
- Provides unique pressure-relieving comfort by replacing the back section of a Miami J Collar.
- Clinically proven to virtually eliminate occipital breakdown without sacrificing c-spine stability.
- Easy to clean
- MRI, CT & X-ray lucent, MR safe

REQUIRED MEASUREMENTS

- Universal, one size


ORDERING INFORMATION

Description	Item No
Occian Back	ACB-100
Occian Back Junior	ACB-200

Compatible with all sizes of the Miami J®, Miami J® Advanced and the Miami Jr®. Collar front panel not included.

MIAMI JR®

The Miami Jr is the breakthrough collar designed exclusively for children of 12 years and under. The Miami Jr line, which spans the growth cycle, is sized according to age. Össur is the only company that has a paediatric line of collars based on the anatomical differences between adults and children. We realise that children are not just "little adults"-they require innovative design.



INDICATIONS FOR USE

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Post-operative
- Motor Neurone Disease
- Herniated cervical disc

FEATURES AND BENEFITS

- Innovative collar designed exclusively for children 12 years and under anatomy pediatric
- Sorbatex™ II pads make the collar "skin friendly" and protect skin integrity during extended wear
- MRI, CT and X-ray lucent, MR safe

REQUIRED MEASUREMENTS

- Age

ORDERING INFORMATION

MIAMI JR

Size	Collar set*	Collar
0-6 months	MJR-P0	MJ-P0
6-24 months	MJR-P1	MJ-P1
2-6 years	MJR-P2	MJ-P2
6-12 years	MJR-P3	MJ-P3

*Collar sets come with one additional set of replacement pads.

ACCESSORIES

REPLACEMENT PADS

Size	Item No
0-6 months	MJP-P0
6-24 months	MJP-P1
2-6 years	MJP-P2
6-12 years	MJP-P3

PAPOOSE®

The Papoose is an innovative design for the youngest and most challenging patients-newborn infants. The Papoose cradles the baby's body securely while gently immobilizing the head and spine and maintaining airway and spine alignment to help achieve optimal patient outcome.



INDICATIONS FOR USE

- Immobilization of the spine
- Spinal fractures

FEATURES AND BENEFITS

- Securely and gently immobilizes the head and spine in the proper anatomical airway and spinal alignment
- Occipital shape protects against positional plagiocephaly
- True cervical spine immobilizer for infants up to 3 months
- Allows for continued intimacy between caregiver and child
- MRI, CT and X-ray lucent, MR compatible

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

PAPOOSE

Description	Item No
Papoose without front collar	PO-100
Papoose without front with replacement pads	POR-100
Papoose with front collar	PA-100
Papoose with front with replacement pads	PAR-100

ACCESSORIES

REPLACEMENT PADS

Description	Item No
Papoose replacement pads (with front)	PAP-100
Papoose replacement pads (without front)	POP-100

PHILADELPHIA® ADJUSTABLE COLLAR

INDICATIONS FOR USE


- Post-trauma
- Herniated cervical disc
- Cervical spondylitis
- Post-operative
- Cervical spondylosis
- RA and OA cervical spine

FEATURES AND BENEFITS

- Circumference and height adjustable
- Permanent lock on height adjustment
- Anatomically shaped Plastazote® foam (Not made with natural rubber latex, non-toxic and hypoallergenic to reduce skin irritation), leaving clearance for cervical veins
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy.
- Well ventilated
- X-ray and CT lucent. MR conditional.

REQUIRED MEASUREMENTS

- Refer to the sizing chart



ORDERING INFORMATION

PHILADELPHIA ADJUSTABLE TRACHEOTOMY COLLAR

Size	Neck circumference
Small (2 tears)	10-13" (25.4-33cm)
Medium (1 tear)	13-16" (33-40.6cm)
Large (no tear)	16-19" (40.6-48.3cm)

NOTE: The table above represents common sizing only. Sizing must be verified by the healthcare provider for each patient individually.

Item No	Size / Description
PHP-AA	Conforms to sizes 2, 3 and 4 of the original Philadelphia Tracheotomy Collar

PHILADELPHIA® TRACHEOTOMY COLLAR

INDICATIONS FOR USE

- Post-trauma
- Herniated cervical disc
- Cervical spondylitis
- Post-operative
- Cervical spondylosis
- Motor Neurone Disease
- RA and OA cervical spine

FEATURES AND BENEFITS

- Not made with natural rubber latex, non-toxic and hypoallergenic
Plastazote® foam reduces skin irritation
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy.

REQUIRED MEASUREMENTS

- Neck circumference and height from tip of chin to sternal notch

ORDERING INFORMATION

Size	Neck circumference	6cm / 2.25" height	8cm / 3.25" height	11cm / 4.25" height	13cm / 5.25" height
X-Small	up to 10" (< 25cm)	PHP-T2XS			
Small	10-13" (25-33cm)	PHP-T2S	PHP-T3S	PHP-T4S	PHP-T5S
Medium	13-16" (33-41cm)	PHP-T2M	PHP-T3M	PHP-T4M	PHP-T5M
Large	16-19" (41-48cm)	PHP-T2L	PHP-T3L	PHP-T4L	PHP-T5L
X-Large	19" up (48cm)	PHP-T2XL	PHP-T3XL	PHP-T4XL	PHP-T5XL

Size	Neck circumference	Height	Item No
Infant	6-8" (15-20cm)	1.25" (3cm)	PHP-TI
Paediatric	8-11" (20-28cm)	1.75" (4.5cm)	PHP-TP

ACCESSORIES

PHILADELPHIA COOLMAX COTTON LINER

Size	Neck circumference	Item No
Infant	6-8" (15-20cm)	PHP-LCI
X-Small	up to 10" (< 25cm)	PHP-LCXS
Small	10-13" (25-33cm)	PHP-LCS
Medium	13-16" (33-41cm)	PHP-LCM
Large	16-19" (41-48cm)	PHP-LCL
X-Large	19" up (> 48cm)	PHP-LCXL

PHILADELPHIA TRACHEOTOMY ATTACHMENTS

Size	Stabilizer	Halo
Pediatric	PHP-OSP	PHP-OHP
Adult	PHP-OS	PHP-OH



UNIVERSAL CERVICAL COLLAR

INDICATIONS FOR USE

- Neck extension injuries (Whiplash)
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc

FEATURES AND BENEFITS

- Medium density foam
- Contact closure for easy fit

REQUIRED MEASUREMENTS

- Mandible to clavicle height

ORDERING INFORMATION

Height	Item No
2" (5cm)	203020
3" (7.5cm)	203000
4" (10cm)	203010

FOAM CERVICAL COLLAR

INDICATIONS FOR USE

- Neck extension injuries (Whiplash)
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc

FEATURES AND BENEFITS

- Easy to apply

REQUIRED MEASUREMENTS

- Neck circumference and mandible to clavicle height
- Ordering Information

Size	Circumference	Neck size	Item No
Small	3" x 18.5" (7cm x 47cm)	11-15" (28-38cm)	203013
Small/Long	3" x 22" (7cm x 56cm)	15-20" (38-51cm)	203014
Medium	3.5" x 20" (9cm x 51cm)	13-18" (33-46cm)	203015
Medium/Long	3.5" x 23" (9cm x 58cm)	16-21" (41-53cm)	203016
Large	4" x 23" (10cm x 58cm)	15-20" (38-51cm)	203017
Large/Long	4" x 25" (10cm x 64cm)	16-22" (41-56cm)	203018



NECLOC®

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Ergonomic two-piece clamshell design and retaining strap secure flexion during application
- Five unique phenotype-driven sizes make it easy to select the correct size
- Minimal manipulation of patient during application
- X-ray and CT lucent. MR conditional



ORDERING INFORMATION

Size	Item No
Stout	NL-200EL
XS	NL-250E
Small	NL-300E
Medium	NL-400E
Large	NL-500E

NECLOC® KIDS

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Ergonomic two-piece clamshell design and retaining strap secure flexion during application
- Three pediatric sizes ensure optimal airway and cervical spine immobilization
- Minimal manipulation of patient during application to children
- Meets the unique developing anatomical requirements of children
- X-ray and CT lucent. MR conditional



REQUIRED MEASUREMENTS

- Age

ORDERING INFORMATION

Age	Item No
0-2 years	NL-P1
2-6 years	NL-P2
6-12 years	NL-P3

PATRIOT® COLLAR

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Multi-height adjustment-Adult (4 positions) and Pediatric (3 positions)
- Open rear ventilation panel for fluid drainage, airflow and cervical spine palpation
- Universal sizing optimizes inventory requirements



REQUIRED MEASUREMENTS

- Neck circumference

ORDERING INFORMATION

Size	Neck circumference	Blue
Child	8-18" (20-46cm)	PHP-PP
Adult	11-23" (28-58cm)	PHP-PA

FORM FIT® BACK SUPPORT, AIR

INDICATIONS FOR USE

- Muscular low back pain
- Lumbosacral Strains
- Chronic low back Pain
- Sciatica

FEATURES AND BENEFITS

- Intuitive integrated pneumatic system for more customized compression and support
- Anti-Migration system reduces slippage
- Double wrap design for improved compression
- Finger loops for easy application
- Hot/Cold Gel pad available separately

REQUIRED MEASUREMENTS

- Waist circumference



ORDERING INFORMATION

Size	Waist circumference	Standard
Small	25-30" (63-76cm)	209133
Medium	30-34" (76-86cm)	209135
Large	34-38" (86-97cm)	209137
X-Large	38-42" (97-107cm)	209138
XX-Large	42-56" (107-142cm)	209139

ACCESSORIES

Description	Item No
Gel pad	GB28

Description	Item No
Back Support Suspenders Clip On	10912

FORM FIT® BACK SUPPORT

INDICATIONS FOR USE

- Lumbar strains
- General posture support
- Lumbosacral strains
- Sciatica

FEATURES AND BENEFITS

- Muscular low back pain
- Lumbosacral strains
- Chronic low back pain
- Sciatica

REQUIRED MEASUREMENTS

- Waist circumference



ORDERING INFORMATION

Size	Waist Circumference	Standard
Small	25-30" (63-76cm)	209143
Medium	30-34" (76-86cm)	209145
Large	34-38" (86-97cm)	209147
X-Large	38-42" (97-107cm)	209148
XX-Large	42-56" (107-142cm)	209149
XXX-Large	56-60" (142-152cm)	209151

ACCESSORIES

Description	Item No
Back Support Suspenders Clip On	10912
Hot/Cold Therapy Gel Pad	GB28

BACK SUPPORT SUSPENDERS

FEATURES AND BENEFITS

- Suspenders can be added to each support upon request

REQUIRED MEASUREMENTS

- Please specify sewn-in or removable suspender design



ORDERING INFORMATION

Description	Item No
Clip-on	10912



INNOVATOR X®

Benefiting from our proven, patented technologies, Innovator X, with extending struts, can help improve patient outcomes. With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, Innovator X is designed to increase patient satisfaction and compliance.



INDICATIONS FOR USE

- Controlled elbow flexion/extension following surgery or trauma

FEATURES AND BENEFITS

- Easy-to-use Innovator hinge
- Slide-to-Size™ struts lead to a more customized fit
- Sensil® silicone padding on the strapping system reduces slipping and malleable aluminium cuffs allow for customization
- Patented Form Fit® 3-dimensional padding and lycra lining provide superior comfort and compliance
- Comfort pad at the neck and shoulder
- Flexion range can be quickly adjusted to 0, 10, 20, 30, 45, 60, 75, 90, 105 and 120. Extension range can be adjusted to 0, 10, 20, 30, 45, 60, 75 and 90
- Optional anti-supination/pronation arm bar kit is available

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

INNOVATOR X

Side	Item No
Left	W-823101
Right	W-823100

ACCESSORIES

Description	Item No
Innovator X Bar Kit	W-A82310
Replacement Softgoods Left	W-823111
Replacement Softgoods Right	W-823110

FORM FIT® SHOULDER BRACE

The Form Fit Shoulder Brace is designed to provide immobilization following shoulder injury/surgery. With its streamlined fit, sporty look and breathability it offers excellent comfort and ease of use. Same as its predecessor, the SmartSling®, the Form Fit Shoulder Brace also provides abduction of the shoulder (when the abduction pillow is attached), or external rotation of the shoulder (when the external rotation pillow is attached).

Optional axilla pillow helps distribute weight over a larger surface area and decreases pressure on the ulnar nerve.



INDICATIONS FOR USE

- Immobilization for rotator cuff repairs
- Bankart Repairs
- Capsular Shift
- Gleno humeral dislocation/subluxations

FEATURES AND BENEFITS

- Breathable material
- Quick-release buckles
- Offloading strap
- Extra padding
- Universal thumb loop



Axilla Pillow



External rotation pillow



REQUIRED MEASUREMENTS

- Olecranon to end of closed fist

ORDERING INFORMATION

FORM FIT SHOULDER BRACE WITH ABDUCTION

Size	Fits Length	Part
Small	10.5-13" (up to 29cm)	SSS
Medium	13.25-15.75" (29-40cm)	SSR
Large	16-17.5" (40-44.5cm)	SSL

FORM FIT SHOULDER BRACE

Size	Fits Length	Part
Small	10.5-13" (up to 29cm)	SSS-SO
Medium	13.25-15.75" (29-40cm)	SSR-SO
Large	16-17.5" (40-44.5cm)	SSE-SO

FORM FIT SHOULDER BRACE WITH EXTERNAL ROTATION

Size	Fits Length	Part
Small	10.5-13" (up to 29cm)	SSS-ER
Medium	13.25-15.75" (29-40cm)	SSR-ER
Large	16-17.5" (40-44.5cm)	SSE-ER

ACCESSORIES

Description	Part
Axilla Pillow	SSAXP
External Rotation Pillow	SSE

SHOULDER ABDUCTION SLING

INDICATIONS FOR USE

- Shoulder injuries
- Post-operative

FEATURES AND BENEFITS

- Provides shoulder abduction at 10° with support
- Cushion rests on waist to support weight of arm
- Breathable cotton provides a comfortable fit
- Fits left and right

REQUIRED MEASUREMENTS

- Measure elbow crease to fingertip



ORDERING INFORMATION

Size	Elbow crease to fingertip	Item No
Small	up to 12.5" (< 32cm)	204053
Medium	12.5-15.5" (32-39cm)	204055
Large	15.5" and up (39cm +)	204057

PREMIUM SHOULDER IMMOBILIZER

INDICATIONS FOR USE

- Post-trauma
- Post-operative

FEATURES AND BENEFITS

- Comfortable fit
- Easily adjusts
- One size fits all (max circumference 45")



ORDERING INFORMATION

Description	Item No
Universal	204100

PADDED STRAP SHOULDER IMMOBILIZER

FEATURES AND BENEFITS

- Padded pressure sensitive straps improve comfort
- Waist strap is easily adjusted or can be removed
- Deep pocket design
- Two adjustment points for each strap allows for better fit
- Hook and loop closure make it easy to adjust



ORDERING INFORMATION

Size	Length	Height	U/M	Item No
X-Small	11" (28cm)	7.5" (44cm)	1 ea	204120
Small	13.5" (34cm)	7.5" (44cm)	1 ea	204122
Medium	15.5" (39cm)	7.5" (44cm)	1 ea	204123
Large	17.5" (44cm)	7.5" (44cm)	1 ea	204125
X-Large	20" (51cm)	7.5" (44cm)	1 ea	204127

PREMIUM CONTACT CLOSURE ARM SLING

INDICATIONS FOR USE

- Post upper limb trauma
- Post upper limb surgery

FEATURES AND BENEFITS

- Envelope style with thumb support
- Individually boxed
- Includes double "D" ring and slide buckle adjustment
- Canvas Material



ORDERING INFORMATION

Size	Length	Depth	Strap length	Standard	With Comfort Pad
X-Small	12" (30cm)	4" (10cm)	28" (71cm)	204112	P204112
Small	14" (35cm)	5" (13cm)	30" (76cm)	204113	P204113
Medium	17" (43cm)	6.75" (17cm)	35" (89cm)	204115	P204115
Large	19" (48cm)	6.75" (17cm)	38" (97cm)	204117	P204117
X-Large	19.5" (50cm)	8" (20cm)	43" (109cm)	204118	P204118

BUCKLE CLOSURE ARM SLING

INDICATIONS FOR USE

- Post-trauma
- Post-surgery

FEATURES AND BENEFITS

- Contact closure shoulder strap with clip buckle swathe
- Available in 6 pack only



ORDERING INFORMATION

Size	Length	Depth	Strap length	Item No
Child	10.5" (27cm)	5" (13cm)	24" (61cm)	204020
X-Small	12.25" (31cm)	4" (10cm)	32" (81cm)	204022
Small	14.5" (37cm)	5.25" (13cm)	35" (89cm)	204023
Medium	17.5" (45cm)	6.5" (17cm)	39.5" (100cm)	204025
Large	19" (48cm)	7.5" (19cm)	43.5" (110cm)	204027

BERGSCHULTZ UNIVERSAL ARM SLING

INDICATIONS FOR USE

- Post-trauma
- Post-operative

FEATURES AND BENEFITS

- Designed by Jim Bergschultz. R.O.T.
- Blue heavy denim material with contact closure

REQUIRED MEASUREMENTS

- Universal, one size



ORDERING INFORMATION

Description	Item No
Bergschulz Universal Arm Sling	204070

PREMIUM PADDED ARM SLING

INDICATIONS FOR USE

- Post-trauma
- Post-surgery

FEATURES AND BENEFITS

- Blue envelope style
- Canvas Material



ORDERING INFORMATION

Size	Length	Depth	Strap length	Item No
Small	13" (33cm)	7" (18cm)	31" (79cm)	204313
Medium	16" (41cm)	8" (20cm)	35" (89cm)	204315
Large	19" (48cm)	9.5" (24cm)	36" (91cm)	204317
X-Large	22" (56cm)	9.5" (24cm)	40" (102cm)	204318

FRONT CLOSURE CLAVICLE SUPPORT

INDICATIONS FOR USE

- Clavicle fractures and dislocations
- Postural support
- Soft tissue injuries

FEATURES AND BENEFITS

- Features front closure design for convenience in application
- Comfortable padding over clavicle area

REQUIRED MEASUREMENTS

- Chest circumference at the sternum area



ORDERING INFORMATION

Size	Circumference	Item No
X-Small	20-24" (51-61cm)	205032
Small	24-30" (61-76cm)	205033
Medium	30-36" (76-91cm)	205035
Large	36-42" (91-107cm)	205037
X-Large	42-48" (107-122cm)	205038
XX-Large	48-54" (122-137cm)	205039

FIGURE-8 CLAVICLE SPLINT

INDICATIONS FOR USE

- Clavicle fractures and dislocations
- Postural support
- Soft tissue injuries

FEATURES AND BENEFITS

- Prong buckle closure ensures correct fit without slippage

REQUIRED MEASUREMENTS

- Chest circumference



ORDERING INFORMATION

Size	Circumference	Item No
X-Small	20-24" (51-61cm)	205002
Small	24-30" (61-76cm)	205003
Medium	30-36" (76-91cm)	205005
Large	36-42" (91-107cm)	205007
X-Large	42-48" (107-122cm)	205008

TENNIS ELBOW SUPPORT WITH HOT/COLD GEL THERAPY

INDICATIONS FOR USE

- Epicondylitis

FEATURES AND BENEFITS

- Universal 2" / 5cm wide with loop-lock Velcro closure
- Loop-lock contact closure
- Twin blue cushions provide counter force, reducing traction on the tendons and inflammation

REQUIRED MEASUREMENTS

- Circumference of the elbow



ORDERING INFORMATION

Description	Size	Item No
Tennis Elbow Support with hot/cold gel therapy	Universal (fits elbow circumference 8-14"/20-36cm)	E310000

FORM FIT® TENNIS ELBOW

INDICATIONS FOR USE

- Epicondylitis

FEATURES AND BENEFITS

- Patented ergonomically-shaped plastic eliminates tourniquet effect associated with traditional designs
- Twin blue cushions provide counterforce, reducing traction on the tendons and thus inflammation
- Ribbed surfaces reduce vibration to the tendon, increasing comfort
- Durable design and materials enable everyday use for both work and athletic activities

REQUIRED MEASUREMENTS

- Circumference of the elbow



ORDERING INFORMATION

Size	Circumference	Item No
Small	8.7-10.4" (22-26cm)	B-247000002
Medium	10.4-12.8" (27-32cm)	B-247000003
Large	12.8-15" (33-38cm)	B-247000004

AIRFORM® TENNIS ELBOW SUPPORT

INDICATIONS FOR USE

- Epicondylitis

FEATURES AND BENEFITS

- Focused pneumatic compression
- Conforms and comfortably cushions sensitive injured anatomy
- Easy to apply, simple design that fits most, reduces inventory needs

REQUIRED MEASUREMENTS

- Universal: Fits elbow circumference 20-36cm / 8-14"



ORDERING INFORMATION

Description	Size	Item No
AirForm Tennis Elbow Support	Universal (fits elbow circumference 8-14"/20-36cm)	320000

EXOFORM® WRIST

INDICATIONS FOR USE

- Mild to Moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other select soft tissue injuries

FEATURES AND BENEFITS

- Ergonomically contoured padding
- Rigid plastic supports contour to the hand and wrist
- Adjustable aluminum palmar stay allows for fine-tuning when needed

REQUIRED MEASUREMENTS

- Wrist circumference



ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
X-Small	5.25-6.25" (13-16cm)	507082	507072
Small	6.25-7.25" (16-18.5cm)	507083	507073
Medium	7.25-8" (18.5-20.5cm)	507085	507075
Large	8-9" (20.5-23.5cm)	507087	507077
X-Large	9-10" (23.5-25.5cm)	507088	507078

EXOFORM® CARPAL TUNNEL WRIST

INDICATIONS FOR USE

- Carpal tunnel syndrome
- Tendonitis

FEATURES AND BENEFITS

- Allows for full MCP flexion
- Contoured for anatomical fit
- Provides radial and ulnar styloid relief
- Low profile
- Preshaped into a neutral position for comfortable limited motion/ flexibility

REQUIRED MEASUREMENTS

- Wrist and palm circumference



ORDERING INFORMATION

Size	Wrist Circumference	Palm Circumference	Left	Right
Small	5.25-6.5" (13-16.5cm)	6.75-8" (17-20.5cm)	517083	517073
Medium	6.5-7.75" (16.5-19cm)	8-9.25" (20.5-23.5cm)	517085	517075
Large	7.75-10" (19-25.5cm)	9.25-10.5" (23.5-26cm)	517087	517077

FORM FIT® WRIST / FORM FIT® WRIST & FOREARM

Form Fit Wrist, and Form Fit Wrist & Forearm, provides immobilization and protection for patients suffering from tendonitis, select soft tissue injuries-including mild to moderate sprains and strains, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.



INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

FEATURES AND BENEFITS

- Breathable material
- Contour fit
- Contact Closure Straps
- Customizable Stays

REQUIRED MEASUREMENTS

- Wrist circumference

ORDERING INFORMATION

FORM FIT WRIST & FOREARM

Size	Wrist Circumference	8" (20cm) LENGTH FORM FIT® WRIST		10" (25cm) LENGTH FORM FIT® WRIST & FOREARM	
		Left	Right	Left	Right
XSmall	5.25-6.25" (13-16cm)	317082	317072	617082	617072
Small	6.25-7.25" (16-18.5cm)	317083	317073	617083	617073
Medium	7.25-8" (18.5-20.5cm)	317085	317075	617085	617075
Large	8-9" (20.5-23.5cm)	317087	317077	617087	617077
XLarge	9-10 (23.5-25.5cm)	317088	317078	617088	617078

ACCESSORIES

GEL PAD FOR HOT/COLD THERAPY

Description	Left	Right
Form Fit Wrist 8" (20cm) version*	GB20L	GB20R
Form Fit Wrist & Forearm 10" (25cm) version	GB22L	GB22R



* Fits XSmall/Small. For Medium-XLarge, please use Form Forearm 10" part numbers

FORM FIT® WRIST & THUMB, UNIVERSAL

Form Fit Universal Wrist and Thumb products are designed to restrict motion of the wrist and the thumb. Providing immobilization and protection from select soft tissue injuries, mild to moderate sprains and sprains, including gamekeeper's thumb and post cast healing. The lightweight, universal design also reduces inventory needs.



INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

FEATURES AND BENEFITS

- Quick-Pull close strap
- Positional dorsal stay
- Customizable thumb stay
- Malleable palmar stay

REQUIRED MEASUREMENTS

- Wrist circumference

ORDERING INFORMATION

FORM FIT WRIST

Size	Wrist Circumference	Left	Right
Pediatric	3.5-5" (9-13cm)	B-252603301	B-252603300
8" (20cm) Universal	4.5-9.5" (11-24cm)	B-252603109	B-252503109
8" (20cm) X-Large	7.5-11" (19-28cm)	B-252603106	B-252603105

FORM FIT WRIST & FOREARM UNIVERSAL

Size	Wrist Circumference	Left	Right
10" (25cm) Universal	4.5-9.5" (11-24cm)	B-252603210	B-252603209

FORM FIT THUMB SPICA UNIVERSAL

Size	Wrist Circumference	Left	Right
Pediatric	3.5-5" (9-13cm)	B-253603301	B-253603300
8" (20cm) Universal	4.5-9.5" (11-24cm)	B-253603109	B-253503109
8" (20cm) X-Large	7.5-11" (19-28cm)	B-253603106	B-253603105

FORM FIT® THUMB SPICA

The Form Fit® Thumb Spica is designed to restrict wrist and thumb motion. It provides immobilization and protection for patients suffering from tendonitis (e.g. DeQuervain's Tendonitis), select soft tissue injuries-including mild to moderate sprains and strains and Gamekeeper's Thumb -, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.



INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
 - e.g., DeQuervain's Tendonitis
- Post-cast healing
- Other selected soft tissue injuries
 - e.g., Gamekeeper's Thumb

FEATURES AND BENEFITS

- Breathable material
- Contour fit
- Contact Closure Straps
- Customizable Stays

REQUIRED MEASUREMENTS

- Wrist circumference

ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
XSmall	5-6" (13-16cm)	3020	3120
Small	6-7" (16-18.5cm)	3030	3130
Medium	7-8" (18.5-20.5cm)	3050	3150
Large	8-9" (20.5-23.5cm)	3070	3170
XLarge	9-10" (23.5-25.5cm)	3080	3180

ACCESSORIES

GEL PAD FOR HOT/COLD THERAPY

Description	Left	Right
Form Fit Wrist 8" (20cm) version*	GB20L	GB20R
Form Fit Wrist & Forearm 10" (25cm) version	GB22L	GB22R

* Fits XSmall/Small. For Medium-Xlarge, please use Form Forearm 10" part numbers



UNIVERSAL THUMB SPLINT

INDICATIONS FOR USE

- Gamekeeper's thumb
- Torn collateral ligaments
- DeQuervain Tendonitis

FEATURES AND BENEFITS

- Elastic contact closures
- Restricts movement of the MP/CMC
- Light weight and breathable

REQUIRED MEASUREMENTS

- Universal, one size



ORDERING INFORMATION

Description	Item No
Universal Thumb Splint	30000

6" (20CM) THUMB AND WRIST

INDICATIONS FOR USE

- Gamekeeper's thumb
- Torn collateral ligaments
- DeQuervain Tendonitis

FEATURES AND BENEFITS

- Contact closure straps
- Contoured fit
- Malleable stays

REQUIRED MEASUREMENTS


- Wrist circumference



ORDERING INFORMATION

Description	Item No	Description	Item No
Wrist Brace 6" (20cm) Right Sm	417073C	Thumb Spica 6" (20cm) Right Sm	30304C
Wrist Brace 6" (20cm) Right Md	417075C	Thumb Spica 6" (20cm) Right Md	30504C
Wrist Brace 6" (20cm) Right Lg	417077C	Thumb Spica 6" (20cm) Right Lg	30704C
Wrist Brace 6" (20cm) Right Xlg	417078C	Thumb Spica 6" (20cm) Right Xlg	30804C
Wrist Brace 6" (20cm) Left Sm	417083C	Thumb Spica 6" (20cm) Left Sm	31304C
Wrist Brace 6" (20cm) Left Md	417085C	Thumb Spica 6" (20cm) Left Md	31504C
Wrist Brace 6" (20cm) Left Lg	417087C	Thumb Spica 6" (20cm) Left Lg	31704C
Wrist Brace 6" (20cm) Left Xlg	417088C	Thumb Spica 6" (20cm) Left Xlg	31804C

HIP



REBOUND® HIP

Össur's Rebound Hip brace has a low-profile, lightweight design that supports patient comfort and compliance. It delivers gentle controlling forces and effective functional ROM restriction in extension/flexion and abduction/adduction to promote optimal post-surgical outcomes. The universal frame and softgood are easy to stock, assemble and quickly fit on patients.



INDICATIONS FOR USE

Hip conditions that could benefit from motion restriction following surgery:

- Femoral Acetabular Impingement arthroscopic surgery
- Gluteus medius repair surgery
- Hamstring repair surgery

FEATURES AND BENEFITS

- Low-profile, lightweight design that supports patient comfort and compliance
- Quick-release pegs easily adjust the desired range-of-motion needed for the patient's optimal rehab
- Quick-fit with sliding mechanism on frame to adjust for patient's height
- Optional, adjustable ROM hinge from -10°, 0°, 10°, or 20°

REQUIRED MEASUREMENTS

- Circumference at iliac crest height
- Length from greater trochanter to mid patella

ORDERING INFORMATION

REBOUND HIP

Description	Item No	
	Preset to Left	Preset to Right
Rebound Hip	980130RB	980135RB
Rebound Hip with adjustable hinge	980120RB	980125RB

*Frame: Universal; Softgood: Universal fitting to 48" (122cm) hip circumference

ACCESSORIES

Description	Item No
OR Kit	980301
Belt Extender (Fitting 49-67"/124-170cm hip circumference)	980121RB
Fitting Tool	B-040100000



KNEE



REBOUND® CARTILAGE

The success of an articular cartilage repair is dependent on patient compliance to the rehabilitation process. The protocol is based on restricted cartilage load combined with gradual increase in activity. Many rehabilitation protocols suggest crutches for partial weight-bearing and it is the standard of care for successful outcomes. However, patient compliance to partial weight-bearing instructions with crutches has been shown to be inconsistent.

The Rebound Cartilage is a protective functional solution designed to support regeneration of the knee cartilage after repair procedures such as Microfracture, OATS and (M)ACI. Based on the clinically-proven 3 Points of Leverage System and patented Cartilage Protection Straps, joint unloading is maintained in flexion for protection of the cartilage during the healing process.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimise healing while maintaining function and mobility.



FUNCTIONAL HEALING-WHERE MOBILITY MEANS RECOVERY

INDICATIONS FOR USE

- Articular cartilage defect repair, such as (M)ACI, OATS and Microfracture
- Cartilage knee trauma
- Meniscal repair
- Avascular necrosis
- Condylar bone marrow lesions (bone bruises)

FEATURES AND BENEFITS

- Dual Cartilage Protection Straps™ (CPS) maintain joint unloading in flexion for protection of the cartilage during the healing process
- Breathable wrap around sleeve integrated with Sensil® silicone coated calf liners and doeskin thigh liners provide increased comfort while ensuring adequate suspension
- Color coded Quick Fit buckles enable easy application and removal
- Optional modular Flexion Control kit available when needed per prescribed rehab protocol for range of movement restriction

REQUIRED MEASUREMENTS

- Circumference measurement: 6" (15cm) below mid-patella

CUSTOM STANDARD COLORS



CUSTOM PAINT



REBOUND® CARTILAGE

ORDERING INFORMATION



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the Rebound PCL. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for a Rebound PCL without the need for training.

REBOUND CARTILAGE

Size	Circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2-14" (31-35.5cm)	B-254610002	B-254510002	B-254620002	B-254520002
Medium	14-16.75" (35.5-42.5cm)	B-254610003	B-254510003	B-254620003	B-254520003
Large	16.75-20" (42.5-51cm)	B-254610004	B-254510004	B-254620004	B-254520004
X-Large	20-24" (51-61.5cm)	B-254610005	B-254510005	B-254620005	B-254520005

REBOUND CARTILAGE-WITH FLEXION CONTROL

Size	Circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2-14" (31-35.5cm)	B-254810002	B-254710002	B-254820002	B-254720002
Medium	14-16.75" (35.5-42.5cm)	B-254810003	B-254710003	B-254820003	B-254720003
Large	16.75-20" (42.5-51cm)	B-254810004	B-254710004	B-254820004	B-254720004
X-Large	20-24" (51-61.5cm)	B-254810005	B-254710005	B-254820005	B-254720005

REBOUND CARTILAGE CUSTOM

Description	Item No
Rebound Cartilage Custom Right	B-154500001
Rebound Cartilage Custom Left	B-154600001

Note: Circumference limitation: Thigh circumferences smaller than 34.3cm or larger than 73.5cm. Calf circumferences smaller than 25cm or larger than 61.5cm.

REBOUND CARTILAGE CUSTOM -WITH FLEXION CONTROL

Description	Item No
Rebound Cartilage Custom-With flexion control Right	B-154700001
Rebound Cartilage Custom-With flexion control Left	B-154800001

REBOUND® CARTILAGE

ACCESSORIES

WRAPAROUND SLEEVE (WITH INTEGRATED SILICONE CALF LINER AND DOESKIN THIGH LINER)

Size	Right medial/Left lateral	Right lateral/Left medial
Small	B-254900072	B-254900082
Medium	B-254900073	B-254900083
Large	B-254900074	B-254900084
X-Large	B-254900075	B-254900085
Optional sleeveless liner kit	Right medial/Left lateral	Right lateral/Left medial
Silicone calf liner, doeskin thigh liner, strap pads	B-254900070	B-254900080
Doeskin calf and thigh liners, suspension strap kit, strap pads	B-427705000	B-427705001

STRAP KITS

Size	Complete strap kit (incl. 4 straps)	CPS strap kit (incl. 2 straps)
Small	B-701290002	B-700290002
Medium	B-701290003	B-700290003
Large	B-701290004	B-700290004
X-Large	B-701290005	B-700290005

UNIVERSAL FLEXION CONTROL KIT

Description	Item No
Left medial	B-705296100
Right medial	B-705295100
Left lateral	B-705296200
Right lateral	B-705295200

OTHER ACCESSORIES

Description	Item No
Flexion Stop Kit	B-729000000
Velcro Twist-Fit Attachment Kit	B-705290000
Suspension Strap	UNL00213

REBOUND® PCL

The PCL, compared to the ACL, has the intrinsic ability to heal and regain continuity after an injury. The PCL experiences variable tensile forces during knee flexion. Insufficient healing after a PCL rupture leads to significant pathologies, such as, chronic anterior knee pain or knee OA. Healing of the PCL in an elongated position can lead to chronic instability and disability. Therefore, physiologically correct positioning of the PCL during conservative or post-surgical rehabilitation is important.

The custom made Rebound PCL is the world's first Dynamic Force PCL brace. It applies a physiologically correct dynamic force to the PCL during rehabilitation.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimise healing while maintaining function and mobility.

FUNCTIONAL HEALING-WHERE MOBILITY MEANS RECOVERY

INDICATIONS FOR USE

- Non-surgical treatment of PCL ruptures
- Post-surgical rehabilitation for
- Partial PCL tears
- PCL reconstruction / augmentation
- PCL revision

FEATURES AND BENEFITS

- Dynamic Tension System (DTS) is designed to increase the load on the tibia as the knee goes through flexion, providing an anterior drawer force necessary to decrease load on the healing PCL to allow it to heal 3 options to customize the standardized, biochemically correct posterior tibial load for each patient by using shear knobs. And having the possibility to re-adjust the posterior tibial load, e.g. in case of swelling
- Lightweight brace with proprietary Sensil® silicone calf liner and doeskin thigh liner reduce migration while providing comfort to the patient
- Custom molded posterior frame for optimized fit and support, with open anterior for ease of donning/doffing
- Available in 12 colors

REQUIRED MEASUREMENTS

- Rebound PCL Custom: Össur's SmartMeasure or DMS (Digital Measurement System)
- Rebound PCL:
 - * M/L measurement: In a standard, weight bearing position
 - Using a Caliper: M/L at joint width
 - Using a Constant Force Caliper: M/L at knee axis line
 - ** Calf circumference measurement: 15cm below mid-patella
 - *** Thigh circumference measurement: 15cm above mid-patella

STANDARD COLORS



REBOUND® PCL

ORDERING INFORMATION

REBOUND PCL

Description	Right	Left	Caliper M/L*	Calf Circ.**	Thigh CIRC.***
Small	B-124501158	B-124501161	8.9-10.2cm	29-34.5cm	42-45cm
Medium	B-124501159	B-124501162	10.2-11.4cm	34.5-38cm	45-51cm
Large	B-124501160	B-124501163	11.4-12.7cm	38-42cm	51-58.5cm

REBOUND PCL CUSTOM

Description	Left	Right	Measurement
Rebound PCL Brace Custom	B-124501155	B-124501154	SmartMeasure iOS App

ACCESSORIES

STRAPS

Size	Item No
X-Small-Large	PCL00061
X-Large-XX-Large	PCL00062

OTHER ACCESSORIES

Description	Item No
Extension Stop Kit	B-705130042
Flexion Stop Kit	B-705130043
Shear Knob Kit	B-465350087
Condyle Pad Kit	B-705130046

LINER KITS

STANDARD UPPER LINER KIT

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
Left	B-124501170	B-124501171	B-124501172	B-124501173	B-124501174	B-124501175
Right	B-124501164	B-124501165	B-124501166	B-124501167	B-124501168	B-124501169

Includes: doeskin thigh liner, doeskin femoral cuff liner (softgood on lower thigh strap), doeskin pre-femoral liner (softgood on upper thigh strap), Velcro (4x round coins, 8x rec tab)

STANDARD LOWER LINER KIT

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
Left	B-124601176	B-124601177	B-124601178	B-124601179	B-124601180	B-124601181
Right	B-124501176	B-124501177	B-124501178	B-124501179	B-124501180	B-124501181

Includes: silicone calf liner / AMS, doeskin tibia cuff liner (softgood on lower calf strap), Velcro (2x round coin, 4x rec tab)

DOESKIN LOWER LINER KIT

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
Left	B-124601182	B-124601183	B-124601184	B-124601185	B-124601186	B-124601187
Right	B-124501182	B-124501183	B-124501184	B-124501185	B-124501186	B-124501187

Includes: doeskin calf liner / AMS, doeskin tibia cuff liner (softgood on lower calf strap), tibia strap pad, Velcro (2x round coins, 4x rec tab)

CTI® CUSTOM

Introduced in 1983, CTi remains the cornerstone of our ligament bracing product line today. While the original CTi has evolved over time, our basic philosophy remains unchanged—a static support system consisting of a rigid frame will hold the bones in proper alignment and provide high level of control to prevent displacement. The constant development and improvement of this classic design is why the CTi remains the top choice of athletes and leading physicians around the world.

The rigid, carbon composite CTi knee brace frame captures the tibia, providing the brace-to-bone contact needed for maximum stabilization of the knee joint. The non-elastic cruciate strapping combined with a rigid frame that truly captures the tibia, results in a static support system.

The CTi Custom is generally used for high levels of support where comfort and strength are not compromised with weight and bulkiness. It can also offer greater advantages for accurate fitting where leg shape abnormalities present a problem.



INDICATIONS FOR USE

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with extension stops
- Proprietary Sensil® silicone liners designed to reduce migration
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Pre-preg carbon fiber, rigid frame construction
- Non-corrosive materials, ideal for water sports
- Custom paint option available
- Standard model for moderate to high impact levels
- Pro Sport model for the highest impact levels
- Vapor model for low to moderate impact levels
- Available in 12 colors


REQUIRED MEASUREMENTS

- Custom-made to measure

STANDARD COLORS



CUSTOM PAINT



CTI® CUSTOM

ORDERING INFORMATION

CTI CUSTOM

Description	Left	Right
Standard	B-116600010	B-116500010
Pro Sport	B-117600010	B-117500010
Vapor	B-118600010	B-118500010



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the CTi. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for a CTi without the need for training.



CUSTOM MEASURING DEVICE (CMS)

Simplified measuring technique saves fitting time. The need for costly casting materials is eliminated. Provides coordinates of patient's leg from which a three-dimensional chart is produced for brace manufacturing. Designed exclusively for use with the CTi knee brace. The CTi measuring device is provided free of charge with the initial CTi brace order.

Description	Item No
CTi Custom Measuring Device	42711
Constant Force Caliper	B-312000000

Contact your local representative for technical support

ACCESSORIES

CALF LINERS

Description	Item No
CTi Custom liner	B-411500115

STRAPS KITS

Size	Item No
Small, Medium	B-701130104
Large, X-Large, XX-Large	B-701130108

PCL KITS

Size	Item No
Small, Medium	B-705000132
Large, X-Large, XX-Large	B-705000235

CTI® CUSTOM

ACCESSORIES

STRAP PADS

Size	Item No
Small	B-304010057
Medium, Large	B-304010058
X-Large, XX-Large	B-304010059

AMS (ANTI-MIGRATION SYSTEM) WRAPS

Description	Item No
AMS Wrap EVA	22095
AMS Wrap Neo	22094
AMS Wrap Neo XL	22152

MX GEAR

Description	Item No
Patella Cup	B-711130011
Gear Guards	B-705000003
MX Kit S (patella cup, gear guards, undersleeve, stickers)	B-711140002
MX Kit M (patella cup, gear guards, undersleeve, stickers)	B-711140003
MX Kit L (patella cup, gear guards, undersleeve, stickers)	B-711140004
MX Kit XL (patella cup, gear guards, undersleeve, stickers)	B-711140005

SLEEVES

Size	Circumference 15cm above mid patella	Softsleeve	Undersleeve	Oversleeve	Sportsleeve
Small	<15" (38cm)	20643	20251	20051	20682
Medium	15-18" (38-45cm)	20644	20252	20052	20683
Large	18-19.5" (45-49.5cm)	20644	20253	20053	20684
X-Large	19.5-21.5" (49.5-54.5cm)	20645	20254	20054	20685
XX-Large	21.5-24.5" (54.5-62cm)	20645	20255	20055	N/A
XXX-Large	24.5-28" (62-71cm)	20646	20256	20056	N/A
XXXX-Large	28-34" (71-86cm)	20647	20257	20057	N/A

OTHER ACCESSORIES

Description	Item No
CTi D-Ring Kit Standard 2x 1" (2.5cm) and 2x 1.5" (3.8cm) d-rings	B-704132600
CTi D-Ring Kit 2" (5cm) 2x 2" (5cm) d-rings	B-704132602
CTi Condyle Pad Kit	B-704130150
Extension Stop Kit	22130
Flexion Stop Kit	19254

CTI®

The proven technology of the CTi knee brace is available with off-the-shelf convenience, providing adjustability and an excellent out-of-the-box fit while maintaining the level of support and protection that you expect from a CTi knee brace.

Outfitted with anatomically-correct Accutrac® hinges and flexible subshells and d-rings, the rock-solid CTi carbon frame provides a rigid exoskeleton that stabilizes the knee joint. Breathable liners coated with bio-inert Össur Sensil® Silicone ensure CTi OTS remains properly positioned on the leg, for optimal functional benefit.



INDICATIONS FOR USE

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

FEATURES AND BENEFITS

- Anatomically-correct Accutrac hinges with extension stops to glide with the knee
- Flexible subshell and Adjustable Tibial Strap capture the tibia and allow adjustments in hinge depth for a more customized fit
- Breathable liners coated with bio-inert Össur Sensil Silicone to reduce migration and skin irritation
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Non-corrosive materials, ideal for water sports
- Available in 2 models: Standard model for moderate to high impact levels, Pro Sport model for the highest impact levels

REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella



CTI®

ORDERING INFORMATION

CTI

Size	Caliper	Circumference	Standard Left	Standard Right	Pro Sport Left	Pro Sport Right
Small	8.9-10.2cm	28-34cm	B-238600112	B-238500112	B-239600112	B-239500112
Medium	10.2 -11.4cm	34-38cm	B-238600113	B-238500113	B-239600113	B-239500113
Large	11.4-12.7cm	38-41cm	B-238600114	B-238500114	B-239600114	B-239500114
X-Large	12.7-14.0cm	41-44cm	B-238600115	B-238500115	B-239600115	B-239500115
XX-Large	14.0-15.2cm	44-48cm	B-238600116	B-238500116	B-239600116	B-239500116

Note: M/L measurements are the preferred and most accurate form of measurement for the CTi. Circumferential measurements are approximate.

ACCESSORIES

LINERS

Size	Silicone tibia cuff	Silicone calf	Silicone thigh	Lycra thigh
Medium	B-304030030	B-425580002	B-425571002	B-425570002
Large		B-425580004	B-425571004	B-425570004
X-Large, XX-Large		B-425580005	B-425571004	B-425570004

STRAPS KITS

Size	Item No
Small, Medium	B-701140103
Large, X-Large, XX-Large	B-701140104

PCL KITS

Size	Item No
Small, Medium	B-705000132
Large, X-Large, XX-Large	B-705000235

STRAP PADS

Size	Item No
Small	B-304010057
Medium, Large	B-304010058
X-Large, XX-Large	B-304010059

AMS (ANTI-MIGRATION SYSTEM) WRAPS

Description	Item No
Silicone AMS Wrap Small, Medium	B-715260002
Silicone AMS Wrap Large, X-Large, XX-Large	B-715260004
Neoprene AMS Wrap Small, Medium	B-410000001
Neoprene AMS Wrap Large, X-Large	B-410000002

CTI®

ACCESSORIES

MX GEAR

Description	Item No
Patella Cup	B-711130011
Gear Guards	B-705000003
MX Kit S (patella cup, gear guards, undersleeve, stickers)	B-711140002
MX Kit M (patella cup, gear guards, undersleeve, stickers)	B-711140003
MX Kit L (patella cup, gear guards, undersleeve, stickers)	B-711140004
MX Kit XL (patella cup, gear guards, undersleeve, stickers)	B-711140005

SLEEVES

Circumference 6" (15cm) above mid patella	Softsleeve	Undersleeve	Oversleeve	Sportsleeve
<15" (38cm)	20643	20251	20051	20682
15-18" (38-45cm)	20644	20252	20052	20683
18-19.5" (45-49.5cm)	20644	20253	20053	20684
19.5-21.5" (49.5-54.5cm)	20645	20254	20054	20685
21.5-24.5" (54.5-62cm)	20645	20255	20055	N/A
24.5-28" (62-71cm)	20646	20256	20056	N/A
28-34" (71-86cm)	20647	20257	20057	N/A

OTHER ACCESSORIES

Description	Item No
CTi D-Ring Kit Standard 2x 1" / 2.5cm and 2x 1.5" / 3.8cm d-rings	B-704132600
CTi D-Ring Kit 2" / 5cm 2x 2" / 5cm d-rings	B-704132602
CTi Condyle Pad Kit	B-704130150
Extension Stop Kit	22130
Flexion Stop Kit	19254

REBOUND® DUAL

The Rebound DUAL knee brace delivers functional support for ligament instabilities, including ligament instability with osteoarthritis, in a single, versatile aluminum brace. Suitable for low-to-high impact activities as well as sliding sports, Rebound DUAL is your new, go-to knee brace.



INDICATIONS FOR USE

- Supports ACL, MCL, LCL, PCL or combined instabilities
- Ligament instabilities with osteoarthritis

FEATURES AND BENEFITS

- Polycentric hinge for stability, durability, and low-profile design
- Adjustable tibia strap captures tibia for optimal stability
- Optional modular SmartDosing® kit available when needed for fine-tuned pain management and offloading
- Breathable ActiveGrip™ calf liner with anti-migration capabilities

REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella

CUSTOM BRACE COLORS



CUSTOM LINER COLORS



REBOUND® DUAL

ORDERING INFORMATION

REBOUND DUAL

Size	Caliper	Circumference	Right	Left
XSmall	3-3.5" (7.6-8.9cm)	15-16.5" (38.1-41.9cm)	B-242510001	B-242610001
Small	3.5-4" (8.9-10.2cm)	16.5-18.5" (41.9-47cm)	B-242510002	B-242610002
Medium	4-4.5" (10.2-11.4cm)	18.5-20.5" (47-52.1cm)	B-242510003	B-242610003
Large	5-5.5" (12.7-14cm)	22.5-24.5" (57.2-62.2cm)	B-242510004	B-242610004
X-Large	5.5-6" (14-15.2cm)	24.5-27.5" (62.2-69.9cm)	B-242510005	B-242610005
XX-Large	6-6.5" (15.2-16.5cm)	27.5-30" (69.9-76cm)	B-242510006	B-242610006
XXX-Large	6-6.5" (15.2-16.5cm)	27.5-30" (69.9-76cm)	B-242510007	B-242610007

M/L at joint width in a standard, weight-bearing position 2 Thigh Circumference: 6" above mid-patella

REBOUND DUAL CUSTOM

Measuring	Right	Left
SmartMeasure, DMS, Coordinate Cast System, Computer Aided Design	B-155510001	B-155610001

REBOUND DUAL ACCESSORIES

Description	Item No
Tibial Wedge Kit	B-701090001
Flexion Stop Kit (0°, 30°, 45°, 60°, 75° and 90°)	B-701090004
XS-LG Strap Kit	B-701090005
XL-XXXL Strap Kit	B-701090006
D-Ring Kit (new version field serviceable rivets)	B-701090059
Extension Stop Kit (5°, 20°, 30° and 40°)	B-701090008
Condyle Pad Kit	B-701090009
Paradigm Strap Pad Kit	B-701090037
Acculign Screw Kit (10 pieces)	B-701090050
PCL Strap Kit XS-LG	B-701090057
PCL Strap Kit XL-XXXL	B-701090058
DUAL Flex Stop 0° (single stop)	DUA20075
DUAL Flex Stop 30° (single stop)	DUA20076
DUAL Flex Stop 0° Kit (2 stops in bag with label)	B-701090062
DUAL Flex Stop 30° Kit (2 stops in bag with label)	B-701090063



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace, including the CTi. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.



DIGITAL MEASUREMENT SYSTEM (DMS)

This system simply requires minimal physical measurements and two digital photos, which can then be sent to customer service. This system is a quick, easy and reliable method of measuring for a CTi without the need for training.

PARADIGM®

INDICATIONS FOR USE

- ACL, MCL, LCL, rotary and combined instabilities
- Protection and stabilization of surgical grafts
- Recommended for low to medium contact / impact

FEATURES AND BENEFITS

- Carbon fiber Triax® shells with flexible subshell allows for anatomical differences while providing firm control
- Hard-anodised, polycentric hinges
- Adjustable Tibial Strap captures the tibia and allows adjustments in hinge depth for a more customized fit
- Full ROM kit is available, enabling post surgical use
- Large, "open window" frame design for comfort and compliance



REQUIRED MEASUREMENTS

- Circumference: 15cm above mid-patella

CUSTOM STANDARD COLORS



ORDERING INFORMATION

Size	Circumference	Left	Right
X-Small	13.5-15.5" (34.3-39.4cm)	B-210600001	B-210500001
Small	15.5-17.5" (39.4-44.5cm)	B-210600002	B-210500002
Medium	17.5-20.5" (44.5-52.1cm)	B-210600003	B-210500003
Large	20.5-23.5" (52.1-59.7cm)	B-210600004	B-210500004
X-Large	23.5-26.5" (59.7-67.3cm)	B-210600005	B-210500005
XX-Large	26.5-29.5" (67.3-74.9cm)	B-210600006	B-210500006
Custom	N/A	B-110600020	B-110500020

ACCESSORIES

Size	Strap and Pad Kit	Liner Kit Left	Liner Kit Right
X-Small	B-701060001	B-703066001	B-703065001
Small	B-701060002	B-703066002	B-703065002
Medium	B-701060003	B-703066003	B-703065003
Large	B-701060004	B-703066004	B-703065004
X-Large	B-701060005	B-703066005	B-703065005
XX-Large	B-701060006	B-703066006	B-703065006

Description	Item No
Flexion / Extension Stop Kit	B-705000504

FORM FIT® KNEE MCL

INDICATIONS FOR USE

- Mild to moderate Medial Collateral Ligament (MCL) sprains (Grade I and Grade II)

FEATURES AND BENEFITS

- Single, steel ROM lateral hinge (0° or 20° extension stop, 70° flexion stop)
- Dynamic Force Strap resist valgus forces
- Hybrid wrap/sleeve for easy application
- Adjustable tibial strap for improved suspension

REQUIRED MEASUREMENTS

- Circumference: 6" (15cm) above mid-patella



ORDERING INFORMATION

Size	Circumference	Left	Right
X-Small	13-15.5" (34.3-39.4cm)	B-237600001	B-237500001
Small	15.5-17.5" (39.4-44.5cm)	B-237600002	B-237500002
Medium	17.5-20.5" (44.5-52.1cm)	B-237600003	B-237500003
Large	20.5-23.5" (52.1-59.7cm)	B-237600004	B-237500004
X-Large	23.5-26.5" (59.7-67.3cm)	B-237600005	B-237500005

C180™ ROCKET

INDICATIONS FOR USE

- ACL, MCL, LCL, PCL, rotary and combined instabilities
- Recommended for low to medium contact/impact
- Recommended for pediatric usage

FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with easy-to-use extension stops
- Heat-mouldable frame for customization
- Adjustable Tibial Strap captures the tibia and allows adjustments in hinge depth for a more customized fit
- Available with an optional patella cup
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Semi-rigid frame consisting of woven carbon fiber infused with acrylic resin
- Lightweight, streamlined design for patient comfort and compliance
- Non-corrosive materials, ideal for water sports

REQUIRED MEASUREMENTS

- Caliper measurement: Medial-lateral at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella



ORDERING INFORMATION

Size	Caliper	Circumference	Left	Right
Small	3-3.5" (7.6-8.9cm)	9.5-11" (24.1-27.9cm)	C180KSL	C180KSR
Large	3.3-3.75" (8.3-9.5cm)	10-12" (26-29.8cm)	C180KLL	C180KLR
X-Large	3.5-4" (8.9-10.2cm)	11-12.5" (27.9-32.1cm)	C1802XSL	C1802XSR

ACCESSORIES

Description	Item No	
Motocross kit Small / Large	10833	
Motocross kit X-Large	B-705130072	
Patella Cup Small / Large	25029	
Patella Cup X-Large	B-711130000	
Calf Liner Small	25061 (L)	26062 (R)
Calf Liner Large	25057 (L)	26058 (R)
Calf Liner X-Large	25106 (L)	25107 (R)
Thigh Liner Small	25060 (L)	26063 (R)
Thigh Liner Large	25056 (L)	26059 (R)
Thigh Liner X-Large	25104 (L)	25105 (R)

REBOUND® KNEE BRACE

With six variations to choose from, Össur's Rebound hinged knee brace product line is a comprehensive portfolio of high-quality braces designed to serve the needs of a broad range of patients.

INDICATIONS FOR USE

- Mild to moderate knee sprains and strains
- Mild to moderate Anterior Cruciate Ligament (ACL) and Posterior Cruciate Ligament (PCL) tears
- Medial Collateral Ligament (MCL) and Lateral Collateral Ligament (LCL) tears

FEATURES AND BENEFITS

- Featuring our unique, 3-layer Cooltech™ fabric, the Rebound Knee Brace is highly breathable, reduces heat retention and dries ultrafast.
- Low profile hinges with extension and flexion stop kit options (ROM version) and malleable arms to provide a custom fit and knee movement control.
- Extension Stops: 0, 5, 10 (pre-installed), 20, 30, 40
- Flexion stops: 45, 60, 75, 90
- Contoured pre-tensed inset fabric that reduces irritation and material bulging in the popliteal area as well as inventory
- Removable straps allow application from right or left side and strap pull tab allow for consistent donning by patient
- Universal buttresses are included with all sleeve versions for customizable patella support

REQUIRED MEASUREMENTS

- Thigh Circumference: 6" (15cm) above mid-patella

ORDERING INFORMATION

REBOUND ROM KNEE

Size	Circumference	Sleeve Long	Sleeve Short	Wrap Long	Wrap Short
X-Small	11.5-13.75" (29-35cm)	708152	708052	707152	707052
Small	13.75-16" (35-41cm)	708153	708053	707153	707053
Medium	16-18" (41-46cm)	708155	708055	707155	707055
Large	18-20.5" (46-52cm)	708157	708057	707157	707057
X-Large	20.5-22.5" (52 -57cm)	708158	708058	707158	707058
XX-Large	22.5-24.75" (57-63cm)	708159	708059	707159	707059
XXX-Large	24.75-28.5" (63-72cm)			707160	707060
XXXX-Large	28.5-32" (72-81cm)			707161	707061

REBOUND KNEE

Size	Circumference	Sleeve Short	Wrap Short
X-Small	11.5-13.75" (29-35cm)	704052	703052
Small	13.75-16" (35-41cm)	704053	703053
Medium	16-18" (41-46cm)	704055	703055
Large	18-20.5" (46-52cm)	704057	703057
X-Large	20.5-22.5" (52 -57cm)	704058	703058
XX-Large	22.5-24.75" (57-63cm)	704059	703059
XXX-Large	24.75-28.5" (63-72cm)	704060	703060
XXXX-Large	28.5-32" (72-81cm)	704061	



INNOVATOR®

INDICATIONS FOR USE

- Post-op/post-injury range of motion control or immobilization

FEATURES AND BENEFITS

- Flexion range can be quickly adjusted to 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105° and 120°. Extension range can be adjusted to 0°, 10°, 20°, 30°, 45°, 60°, 75°, and 90°
- Available in cool or full foam models
- One touch easy-to-adjust ROM hinge
- Breakaway struts for easy, tool-free length adjustments

REQUIRED MEASUREMENTS

- Circumference: 15cm / 6" above mid-patella



ORDERING INFORMATION

Size	Thigh circumference	Cool 61cm	Cool 71cm	Full foam 40cm	Full foam 61cm	Full foam 71cm
Regular	Up to 27" (69cm)	403100	403101	403300	403200	403201
X-Large	27-35" (69-89cm)	403102	403103	403301	403202	403203

ACCESSORIES

Description	Item No
Cool Length Reduction Kit	403150
Strap Length Reduction Kit (Cool)	403199
Strap Length Reduction Kit (Full Foam)	403299
Regular hinge set (Cool)	403140
Regular hinge set (Full)	403400

REPLACEMENT SOFTGOODS

Size	Cool	Full Foam
Regular	403175	403275
Regular Long		403276
X-Large	403176	403277
X-Large Long		403278

Note: will be revamp launched under Formfit brand in late of 2018.

REBOUND® POST-OP KNEE

INDICATIONS FOR USE

- ACL, PCL, MCL, LCL repairs/injuries
- Meniscal repairs
- Tibial plateau fractures
- Patellar tendon repairs
- Osteochondral repairs
- Condylar fractures
- Sprains/strains of the knee
- HTOs (High Tibial Osteotomy)

FEATURES AND BENEFITS

- Available in universal, left or right contour and pediatric
- Telescoping frame from 18.5"-27" (47cm-68.5cm) for Adult and 15"-21" (38cm-53cm) for Pediatric
- OneTouch™ easy to adjust ROM hinge with lock-out capability
- Color-coding and numbered frame for patient education
- Quick-fit, gripping buckles for easy application
- Tool-free design and multi-purpose Smart-Fit™ packaging

REQUIRED MEASUREMENTS

- Thigh circumference 6" (15cm) above mid-patella

ORDERING INFORMATION

Order Information	Brace	Order Information	Accessories
B-8381680	Rebound PO Knee Cool LFT	B-8381082	Kit Reb PO Knee Cool STD Liner
B-8381580	Rebound PO Knee Cool RGT	B-8382082	Kit Reb PO Knee Full STD Liner
B-8381480	Rebound PO Knee Cool UNIV	B-8382092	Kit Reb PO Knee Full XL Liner
B-8381490	Rebound PO Knee Cool UNIV XL	B-8381072	Kit Reb PO Knee Cool PED Liner
B-8382680	Rebound PO Knee Full LFT	B-8381092	Kit Reb PO Knee Cool XL Liner
B-8382580	Rebound PO Knee Full RGT	B-8380005	Kit Reb PO Knee ROM Stop Locks
B-8382480	Rebound PO Knee Full UNIV	B-8380006	Kit Alligator Velcro Kits
B-8382490	Rebound PO Knee Full UNIV XL	B-8380007	Kit Adhesive Velcro Adlt Kit
B-8381470	Rebound PO Knee Cool UNIV PED	B-8380008	Kit Adhesive Velcro Ped Kit
B-8385480	Rebound PO Knee Cool UNIVERSAL NO TRAY	B-703180005+	Incision Pad Kit (Pack of 10)

Note: Pediatric accommodates thigh circumferences 12"- 25" (30.5cm-63.5cm) and frame telescopes from 15"-21" (38cm-53cm)
 Universal, Contoured Left/Right accommodates thigh circumferences: 16"-32" (41cm-81cm) and frame telescopes from 18.5"-27" (47cm-68cm)
 XL accommodates thigh circumferences: 16"-38" (41cm-96.5cm) and frame telescopes from 18.5"-27" (47cm-68cm)



EXOFORM® KNEE IMMOBILIZER

INDICATIONS FOR USE

- Post-injury
- Post-trauma

FEATURES AND BENEFITS

- Slide-to-size adjustment enables you to provide the right fit for the largest variety of patient lengths and widths
- Popliteal pad provides superior fit and limits migration
- Cool foam version can be worn directly against the skin for greater comfort on warm days
- Full foam version provides warmth and compression

REQUIRED MEASUREMENTS

- Universal-fits leg inseams of 18-24" (46-61cm)



ORDERING INFORMATION

Description	Leg circumference	Item No
Cool	up to 28" (71.1cm)	222100
Full foam	up to 28" (71.1cm)	222000
Full foam XXXL	up to 38" (96.5cm)	222051

ACCESSORIES

Description	Item No
Replacement softgood	222001

Description	Item No
Cool	222101
Full foam	222001
3XL full foam	22051

Special Note: This product will be available Q4, 2018.

FORM FIT® KNEE IMMOBILIZER

INDICATIONS FOR USE

- Emergency immobilization
- Pre and post-surgical knee immobilization

FEATURES AND BENEFITS

- Adjustable medial/lateral stays and rigid posterior stay provide stability
- Soft breathable foam keeps patients cool
- Universal 3-piece design simplifies inventory needs

REQUIRED MEASUREMENTS

- Universal



ORDERING INFORMATION

Length	Item No
12" (30cm)	200009
14" (36cm)	200010
16" (41cm)	200016
18" (46cm)	200118
20" (51cm)	200020
22" (56cm)	200022
24" (61cm)	200030

ACCESSORIES

Description	Item No
Hot/Cold Therapy Gel Pad	GB7

FX™ PATELLA STABILIZER

INDICATIONS FOR USE

- Ideal for Patellofemoral Pain Syndrome
- Patellofemoral Stress Syndrome (Runner's Knee)
- Patellar Subluxation/Dislocation
- Chondromalacia
- Patellar Tendinopathy (Tendonitis or Jumper's Knee)
- Arthritis and Osgood-Schlatter disease.

FEATURES AND BENEFITS

- Offers superior patellar tracking through a combination of maximum coronal plane control and minimal flexion and extension resistance
- The detachable thigh strap provides a customizable fit. Additionally, the coronal adjuster allows for further specific fit
- The patented stabilizer offers superior control and assists in restoring normal kinematics to the joint
- Breathable Coolflex sleeve minimizes heat and perspiration

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



ORDERING INFORMATION

Size	Knee circumference	Item No
X-Small	10.5-12" (27-31cm)	B-218000001
Small	12-13.5" (31-34cm)	B-218000002
Medium	13.5-15" (34-38cm)	B-218000003
Large	15-17" (38-43cm)	B-218000004
X-Large	17-19" (43-48cm)	B-218000005
XX-Large	19-21" (48-53cm)	B-218000006

ACCESSORIES

Size	Knee circumference	Item No
XSmall	Osgood-Schlatter Buttress	B-304040037
Medium	Osgood-Schlatter Buttress	B-304040038
XLarge	Osgood-Schlatter Buttress	B-304040039

FORM FIT® HINGED LATERAL J

INDICATIONS FOR USE

- Lateral patellar subluxation and dislocations
- Lateral patellofemoral malalignment
- Patellofemoral maltracking and pain
- Chondromalacia due to lateral maltracking

FEATURES AND BENEFITS

- Breathable material
- Removable hinges
- Easy pull-on
- Lateral J buttress
- Open popliteal

REQUIRED MEASUREMENTS

- Thigh Circumference: 6" above mid-patella
- Calf Circumference: 6" below mid-patella



ORDERING INFORMATION

Size	Thigh circumference	Knee circumference	Calf circumference	Right leg	Left leg
X-Small	13-15.5" (33-39.4cm)	12-13" (30.5-33cm)	10-12" (25.4-30.5cm)	709152	709162
Small	15.5-18.5" (39.4-47cm)	13-14" (33-35.6cm)	12-14" (30.5-35.6cm)	709154	709164
Medium	18.5-21" (47-53.3cm)	14-15" (35.6-38.1cm)	14-16" (35.6-40.6cm)	709155	709165
Large	21-23.5" (53.3-59.7cm)	15-17" (38.1-43.2cm)	16-18" (40.6-45.7cm)	709157	709167
X-Large	23.5-26.5" (59.7-67.3cm)	17-19" (43.2-48.3cm)	18-20" (45.7-50.8cm)	709158	709168
XX-Large	26.5-29.5" (67.3-74.9cm)	19-21" (48.3-53.3cm)	20-22" (50.8-55.9cm)	709159	709169

FORM FIT® KNEE, HINGED

The Form Fit® Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.

The Form Fit Knee brace comes in either a wrap or a sleeve version.

INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Universal buttress
- Open popliteal

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee Circumference	Sleeve	Wrap
XSmall	13"-15.5" (33-39cm)	504252	
Small	15.5-18.5" (39-47cm)	504253	503253
Medium	18.5-21" (47-53cm)	504255	503255
Large	21-23.5" (53-60cm)	504257	503257
XLarge	23.5-26.5" (60-67cm)	504258	503258
XXLarge	26.5-29.5" (67-75cm)	504259	503259
XXXLarge	29.5-32" (75-81cm)		503260

Accessories

Description	Item No
Form Fit Knee Brace Buttress (Sleeve only)	503000



FORM FIT® KNEE ROM

The Form Fit® Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.



The Form Fit Knee brace comes in either a wrap or a sleeve version.

INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Universal buttress
- Open popliteal

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee Circumference	Sleeve		Wrap	
		Long ROM Sleeve	Short ROM Sleeve	Long ROM Wrap	Short ROM Wrap
XSmall	13"-15.5" (33-39cm)	508352	508252	507352	507252
Small	15.5-18.5" (39-47cm)	508353	508253	507353	507253
Medium	18.5-21" (47-53cm)	508355	508255	507355	507255
Large	21-23.5" (53-60cm)	508357	508257	507357	507257
XLarge	23.5-26.5" (60-67cm)	508358	508258	507358	507258
XXLarge	26.5-29.5" (67-75cm)	508359	508259	507359	507259
XXXLarge	29.5-32" (75-81cm)			507360	507260

NEOPRENE HINGED KNEE SUPPORT

INDICATIONS FOR USE

- Moderate ligament instabilities

FEATURES AND BENEFITS

- Accordion molding in the popliteal area reduces bunching behind the knee
- Contoured fit for additional comfort and protection
- The hinges are molded into the brace to maintain a low profile

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee circumference	Item No
Small	11.5-13" (29-33cm)	302153BLK
Medium	13-14.5" (33-37cm)	302155BLK
Large	14.5-16" (37-41cm)	302157BLK
X-Large	16-17.5" (41-44cm)	302158BLK
XX-Large	17.5-19" (44-48cm)	302159BLK



NEOPRENE WRAPAROUND HINGED KNEE SUPPORT

INDICATIONS FOR USE

- For patients with unstable or injured knees that may benefit from additional knee support, such as post-operation rehabilitation, tendonitis, bursitis, arthritis; and sprains and strains.

FEATURES AND BENEFITS

- Dual action hinges incorporate our patented three dimensional moulding technology
- The hinges are molded into the brace to maintain a low profile

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee circumference	Item No
Small	13-14" (33-35.6cm)	212003
Medium	14-15" (35.6-38cm)	212005
Large	15-16" (38-40.6cm)	212007
X-Large	16-17" (40.6-43.2cm)	212008
XX-Large	17-18" (43.2-45.7cm)	212009
3X-Large	18-23" (45.7-58.4cm)	212010
4X-Large	19-25" (48.3-63.5cm)	212011



NEOPRENE KNEE SUPPORT WITH STABILIZED PATELLA

INDICATIONS FOR USE

- For patients suffering from knee conditions that may benefit from additional knee and patella support such as patellar tendonitis, subluxation dislocation, post-operation rehabilitation, post arthroscopic surgery, condromalacia, bursitis, arthritis and tendonitis.

FEATURES AND BENEFITS

- Molded Hardware to maintain a low profile
- Accordion Molding in the popliteal area eliminates bunching
- Soft Cushioned Liner to offer unequaled comfort
- Varied areas of thickness give additional protection
- Lycra Liner offers a breathable contact area for comfortable and extensive use

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee circumference	Item No
Small	11.5-13" (29-33cm)	302563BLK
Medium	13-14.5" (29-33cm)	302565BLK
Large	14.5-16" (29-33cm)	302567BLK
X-Large	16-17.5" (29-33cm)	302568BLK
XX-Large	17.5-19" (29-33cm)	302569BLK



NEOPRENE KNEE SLEEVE

INDICATIONS FOR USE

- For patients suffering from weak and/or injured knees that may benefit from additional knee support.

FEATURES AND BENEFITS

- Unique accordion design enhances comfort by eliminating bunching in the popliteal area
- Open or closed patella

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Circumference	0.3cm Neoprene		0.6cm Neoprene	
		Open Patella	Closed Patella	Open Patella	Closed Patella
Small	11.5-13" (29-33cm)	402013BLK	412013BLK	302013BLK	312013BLK
Medium	13-14.5" (33-37cm)	402015BLK	412015BLK	302015BLK	312015BLK
Large	14.5-16" (37-41cm)	402017BLK	412017BLK	302017BLK	312017BLK
X-Large	16-17.5" (41-44cm)	402018BLK	412018BLK	302018BLK	312018BLK
XX-Large	17.5-19" (44-48cm)	402019BLK	412019BLK	302019BLK	312019BLK



REBOUND® DIABETIC WALKER

Offloading has proven to be a key factor in foot ulcer treatment and entails reducing pressure of the plantar area in general and the ulcer site in particular. The Rebound Diabetic Walker is uniquely designed for effective treatment of plantar foot ulcers, with easy access and re-fitting during wound care, while improving the patient's mobility during the healing process.

The Rebound Diabetic Walker has major advantages over diabetic shoes and shoe inserts in its healing effectiveness and over total contact casting due to its ease of use, speed of application and low cost.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimise healing while maintaining function and mobility.



FUNCTIONAL HEALING-WHERE MOBILITY MEANS RECOVERY

INDICATIONS FOR USE

- Plantar ulcers needing offloading

FEATURES AND BENEFITS

- Effective tri-laminate insole with removable shock-absorbing hexagons that promote ulcer offloading and a Plastazote, heat-moldable top layer that evenly distributes the pressure and minimizes friction
- Unique pressure relief valve for preventing over-inflation, while balancing compression with customized fit
- Positioning Portals allow the clinician to check the fit of the walker at the heel bottom, in order to ensure appropriate positioning of the foot in the walker
- Anatomically designed shell that immobilizes the foot and ankle at 90 degrees while providing protection
- Flex Edge® technology for relief of pressure points
- Rocker Bottom outsole for a smooth roll over while allowing for individual natural toe out
- Soft and breathable, low friction, machine washable liner for hygiene promotion and prevention of infection
- Tamper evident strapping for compliance monitoring

REQUIRED MEASUREMENTS

- Shoe size



REBOUND® DIABETIC WALKER

ORDERING INFORMATION

REBOUND DIABETIC WALKER

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe size UK	Shoe size EU	Item No
Small	4.5-7	5.5-8.5	4-6.5	36-39	B-242500002
Medium	7-10.5	8.5-12	7-10	39.5-44	B-242500003
Large	10.5-12.5	12-14	10.5-12	44.5-46.5	B-242500004
X-Large	12.5-15	14-16.5	12.5-14.5	47-50	B-242500005

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Liners	Insoles
Small	B-418500102	B-418500202
Medium	B-418500103	B-418500203
Large	B-418500104	B-418500204
X-Large	B-418500105	B-418500205

WALKER COVERS

Size	Hygiene covers (2 pcs)	Durable weather cover	Disposable weather cover (20 pcs)
Small	B-418500602	B-418500502	B-418500512
Medium	B-418500603	B-418500503	
Large	B-418500604	B-418500504	B-418500514
X-Large	B-418500605	B-418500505	

DIABETIC SOCKS (2 PCS)

Size	Item No
Small, Medium	B-418500703
Large, X-Large	B-418500704

COMPLIANCE STRAPS (5 PCS)

Size	Item No
Universal	B-418500300



HYGIENE COVER

DURABLE WEATHER COVER

DISPOSABLE WEATHER COVER

ACHILLES WEDGE KIT

INSOLE

REBOUND® AIR WALKER

The Rebound Air Walker features a host of details designed to provide added protection, superior compression, simplified application and usage, improved comfort and, ultimately, successful patient outcomes.



INDICATIONS FOR USE

REBOUND AIR WALKER

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase. If used with the Achilles Wedge Kit, Achilles ruptures can be treated as well.

REBOUND AIR WALKER LOW TOP

Foot and ankle conditions that may benefit from motion restriction, such as:

- Stable foot fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase.

FEATURES AND BENEFITS

- Full circumferential shell with additional semi-soft toe bumper for frontal toe protection
- Soft Flex Edge™ reduces pressure points
- Full rocker bottom is especially engineered to be low and wide to help promote a natural gait
- Stretch-to-fit shell for optimal fit
- Shell can be trimmed for additional stretch in the calf
- Cushioning and breathable softgood provides comfort
- Linked medial/lateral bladders for easy inflation and equal compression
- Integrated pump and air release valve for easy access

REQUIRED MEASUREMENTS

- Shoe size



REBOUND® AIR WALKER

ORDERING INFORMATION

REBOUND AIR WALKER

Size	Men's Shoe Size US	Women's Shoe Size US	Standard	Low Top
X-Small	2-4	3.5-5.5	B-242900001	B-242900061
Small	4-7	5.5-8.5	B-242900002	B-242900062
Medium	7-10.5	8.5-12	B-242900003	B-242900063
Large	10.5-12.5	12-14	B-242900004	B-242900064
X-Large	12.5-15	14-16.5	B-242900005	B-242900065

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Slim Liner / High Top	Slim Liner / Low Top
X-Small	B-418030101	B-418030161
Small	B-418030102	B-418030162
Medium	B-418030103	B-418030163
Large	B-418030104	B-418030164
X-Large	B-418030105	B-418030165

REBOUND ACHILLES WEDGE KIT

Size	Item No
Small	AWK0001
Medium	AWK0005
Large	AWK0009
Universal	AWK00015

WALKER COVERS

Size	Hygiene covers (2 pcs)	Durable weather cover	Disposable weather cover (20 pcs)
Small	B-418500602	B-418500502	B-418500512
Medium	B-418500603	B-418500503	
Large	B-418500604	B-418500504	B-418500514
X-Large	B-418500605	B-418500505	

OTHER ACCESSORIES

Description	Item No
Adhesive Insole for Wedge	WHW-W0107
Strap Extension	B-502100101
Rebound Hightop Socks SM/MD	B-418500703
Rebound Hightop Sock LG/XL	B-418500704
EvenUP Small*	10350
EvenUP Medium*	10351
EvenUP Large*	10352

*Available for all walkers except the ROM walkers

AFO DYNAMIC

The AFO Dynamic provides a strong, lightweight solution for patients with varying degrees of drop foot secondary to a variety of neurological conditions. Össur has taken the carbon fiber technology from its Flex-Foot® product and incorporated it into the new AFO Dynamic to provide increased support over traditional ankle foot orthoses. The energy storing properties of carbon fiber make it an ideal material for this dynamic ankle foot orthosis because it provides real energy return. Best for those people leading busy lives, the AFO Dynamic will ensure a more efficient walk for the user.



INDICATIONS FOR USE

Drop foot secondary to:

- Cerebro Vascular Accident (CVA)
- Other neurological conditions
- Multiple Sclerosis
- Mild knee instability
- Charcot Marie Tooth disease (HMSN)

Contraindications:

- Severe ankle-foot deformities
- Fluctuating edema
- Severe spasticity

FEATURES AND BENEFITS

- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- The tapered carbon heel stores energy and manages forces at heel strike
- Full-length toe lever provides full support while foot is loaded, as well as reduced shock to the contralateral limb
- Creates medial-lateral stability at mid-stance
- Reduces energy consumption by reducing hip hike
- Medial strut makes it less visible in shoe
- Open anterior tibia shell provides superior anatomical fit, and reduces heat and perspiration
- Weight limit: 120kg (265lbs)

REQUIRED MEASUREMENTS

- Shoe size
- Calf height at thickest part of the calf



AFO DYNAMIC

ORDERING INFORMATION

AFO DYNAMIC

Size	Calf Height	Shoe Size US		Shoe Size UK		Shoe Size EUcm		Left	Right
		Male	Female	Male	Female	Male	Female		
X-Small	30-34cm	3-5	4-6.5	2.5-4.5	< 4	34-37	< 37	O-102021	O-102011
Small	32-36cm	5.5-7.5	7-9	5-7	4.5-6.5	37.5-40	37.5-40	O-102022	O-102012
Medium	34-38cm	8-10.5	9.5-12	7.5-10	7-9.5	40.5-44.5	40.5-44	O-102024	O-102014
Large	36-40cm	11-14		10.5-13.5		45-48.5		O-102026	O-102016

ACCESSORIES

ACCESSORY KITS

Size	Item No
X-Small	O-102101
Small	O-102102
Medium	O-102104
Large	O-102106

Note: Accessory kits come with foam padding, upper and lower Velcro straps, and Velcro dots.

AFO LIGHT

AFO Light provides support only where it is needed, eliminating unnecessary bulk and providing a cosmetically-appealing look. Using the best of Flex-Foot® technology and the energy storing properties of carbon fiber, AFO Light offers a strong and lightweight solution for people with mild to moderate Drop Foot.



INDICATIONS FOR USE

- Mild to moderate drop foot
- Contraindications:
 - Severe ankle-foot deformities
 - Medial-lateral instability
 - Severe spasticity
 - Extreme activity
 - Fluctuating edema

FEATURES AND BENEFITS

- Energy return and propulsion assistance of the limb at toe off
- Supports the foot for ground clearance during swing phase
- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- Gradual plantar flexion at initial contact
- Stiffness gradient at heel and the toe correspond to ground reaction forces as they are applied to the foot
- Easy to don and doff
- Weight limit: 100kg (220lbs)

REQUIRED MEASUREMENTS

- Shoe size
- Calf height at thickest part of the calf

ORDERING INFORMATION

AFO LIGHT

Size	Calf height	Shoe Size US		Shoe Size UK		Shoe Size EU		Left	Right
		Male	Female	Male	Female	Male	Female		
Small	32-36cm	Up to 7.5	7-9	< 7	4.5-6.5	< 40	37.5-40	O-121022	O-121012
Medium	34-38cm	8-10.5	9.5-12	7.5-10	7-9.5	40.5-44.5	40.5-44	O-121024	O-121014
Large	36-40cm	11-14		10.5-13.5		45-48.5		O-121026	O-121016
X-Large	38-42cm	12-15		11.5-14.5		46-50		O-121028	O-121018

ACCESSORIES

Description	Item No
Replacement Kit	O-122101

REBOUND FOOT-UP®

Rebound Foot-Up is a lightweight ankle/foot orthosis (AFO) that offers dynamic support for drop-foot or related disorders that require dorsiflexion support. The Rebound Foot-Up supports the ankle dorsiflexion during swing phase. Dynamic support is provided by an elastic strap connected to the shoe by either a transparent plastic inlay placed between shoelaces, or by a nut and bolt fastened through the eyelets. The positioning of the fastening allows for individualizing the support by providing the option of positioning the dorsiflexion support on the medial or lateral side. When using the plastic inlay, the dorsiflexion support becomes neutral.

The Rebound Foot-Up Ankle Cuff has a smart clamp that holds the cuff in place and facilitates one-handed donning and doffing.

Optional: Rebound Foot-Up Foot Wrap can be used when barefoot or in open shoes or sandals. The adjustable elastic strap is attached to the Foot Wrap using a nut and bolt fastened through the eyelets.



INDICATIONS FOR USE

- Drop foot

FEATURES AND BENEFITS

- Dynamic and discreet support for drop-foot or related disorders that require dorsiflexion support
- Improved function, ease of use, comfort, and quality
- Allows for voluntary plantar flexion
- The elastic strap is adjustable and easily replaceable.
- The Rebound Foot-Up Ankle Cuff and Rebound Foot-Up Foot Wrap are made from breathable material that can be worn comfortably for long periods of time.
- The Rebound foot-Up is hand-washable
- Positioned medial or lateral support
- Contoured to the shape of leg
- Adjustable strap length
- Smart clamp facilitates one-handed donning / doffing
- Optional Rebound Foot-Up® Foot Wrap can be used without shoes or in open shoes or sandals.

REQUIRED MEASUREMENTS

- Circumference of ankle, measured around the malleoli

ORDERING INFORMATION

ANKLE CUFF

Prod. nr.	Description	BeigeAnkle Circumference
RFU50001	Rebound Foot-Up Ankle Cuff S/M	18-25cm
RFU50002	Rebound Foot-Up Ankle Cuff L/XL	26-33cm

FOOT WRAP

Prod. nr.	Description	BeigeAnkle Circumference
RFU50101	Rebound Foot-Up Foot Wrap	20-28cm

ACCESSORIES

Prod. nr.	Description
RFU50201	Rebound Foot-Up Attachment Kit (includes plastic inlay, elastic strap, and 1 x nut and screw)
RFU50202	Rebound Foot-Up replacement Nuts And Bolts (5 nuts included: 2 short and 3 long screws)

FOOT-UP®

Foot-Up is a lightweight ankle orthosis that offers dynamic support for Drop Foot or similar complaints. The revolutionary orthosis provides visible improvement in gait by supporting the foot the moment it is raised.



INDICATIONS FOR USE

- Drop foot
- Support of dorsiflexion

FEATURES AND BENEFITS

- Streamlined
- Padded
- Lightweight
- Easy to fit
- Invisible when worn
- Separate inlays available
- Comfortable to wear
- Available in black and beige
- Breathable material

REQUIRED MEASUREMENTS

- Basic version: Circumference of ankle, measured around the malleoli
- Shoeless version: Circumference of foot immediately behind the ball of the foot

ORDERING INFORMATION

FOOT-UP

Size	Ankle Circumference	Beige	Black
Medium	7-8.5" (18-21cm)	07811 -1	07810-1
Large	8.5-10" (22-26cm)	07811 -2	07810-2
X-Large	10-13" (27-33cm)	07811 -3	07810-3

SHOELESS VERSION*

Size	Foot Circumference	Beige	Black
Small	< 8" (20cm)	27811-1	27810-1
Medium	8-9.5" (20-24cm)	27811-2	27810-2
Large	9.5-11" (24-28cm)	27811-3	27810-3
X-Large	> 11" (28cm)	27811-4	27810-4

*Shoeless version item comprises of foot wrap only and must be ordered with the Foot-Up unless the patient already has one.

ACCESSORIES

Description	Beige	Black
Foot-Up plastic insert	17811	17810

FORM FIT® WALKER, AIR

INDICATIONS FOR USE

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use

FEATURES AND BENEFITS

- High and low top versions available
- Integrated pneumatic system
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Standard	Low Top
X-Small	2-4	3.5-5.5	1.5-3.5	33-35.5	A-W0200BLK	A-W0300BLK
Small	4.5-7	6-8	4-6.5	36-39	A-W0400BLK	A-W0500BLK
Medium	7.5-10.5	8.5-11.5	7-10	39.5-44	A-W0600BLK	A-W0700BLK
Large	10.5-12.5	11.5-13.5	10.5-12	44.5-46.5	A-W0800BLK	A-W0900BLK
X-Large	12.5 +	13.5 +	12.5-14.5	47-50	A-W1000BLK	A-W1100BLK

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Standard	Low Top
X-Small	12349	12356
Small	10675	12357
Medium	10678	12358
Large	10681	12359
X-Large	12347	12360

OTHER ACCESSORIES

Description	Item No
10cm softgood width extension	WKSG-1BLK
Tibia and calf plates for use when additional splinting is required	WKSP100
Hot/Cold therapy insert for safe and reusable cryotherapy	WKGEL100
Heel wedge	AWK0005
Adhesive insole for heel wedge	WHW-W0107

FORM FIT® WALKER

INDICATIONS FOR USE

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use

FEATURES AND BENEFITS

- High and low top versions available
- Shock absorbing sole reduces impact during heel strike
- Low, wide rocker bottom promotes a stable, natural gait
- Breathable, seamless liner

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Standard	Low Top
X-Small	2-4	3.5-5.5	1.5-3.5	33-35.5	W0200BLK	W0300BLK
Small	4.5-7	6-8	4-6.5	36-39	W0400BLK	W0500BLK
Medium	7.5-10.5	8.5-11.5	7-10	39.5-44	W0600BLK	W0700BLK
Large	10.5-12.5	11.5-13.5	10.5-12	44.5-46.5	W0800BLK	W0900BLK
X-Large	12.5 +	13.5 +	12.5-14.5	47-50	W1000BLK	W1100BLK

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Standard	Low Top
X-Small	WKSG200BLK	WKSG0300BLK
Small	WKSG400BLK	WKSG500BLK
Medium	WKSG600BLK	WKSG700BLK
Large	WKSG800BLK	WKSG900BLK
X-Large	WKSG1000BLK	WKSG1100BLK

OTHER ACCESSORIES

Description	Item No
10cm softgood width extension	WKSG-1BLK
Tibia and calf plates for use when additional splinting is required	WKSP100
Hot/Cold therapy insert for safe and reusable cryotherapy	WKGEL100
Heel wedge	AWK0005
Adhesive insole for heel wedge	WHW-W0107

RANGE OF MOTION WALKER

INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction or controlled range of motion, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Achilles ruptures
- During postoperative or rehab phase.

FEATURES AND BENEFITS

- Malleable metal struts for adjustable fit
- Integrated ROM setting adjustability from 0-40° of dorsiflexion and plantar flexion in 10° increments
- Ability to lock adjustment tabs to prevent unintended changes

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	With Regular Liner	With Pneumatic Liner
X-Small	up to 4	up to 5.5	<3.5	<35.5	ROM1200	ROM0200
Small	4.5-7	6-8	4-6.5	36-39	ROM1400	ROM0400
Medium	7.5-10.5	8.5-11.5	7-9	39.5-43	ROM1600	ROM0600
Large	10.5-12.5	11.5-13.5	9.5-10.5	43.5-45.5	ROM1800	ROM0800
X-Large	12.5 and up	13.5 and up	11+	46+	ROM2000	ROM1000

DH OFFLOADING WALKER®

INDICATIONS FOR USE

- Treatment of plantar ulcers

FEATURES AND BENEFITS

- Shear force reduction by stabilizing foot and ankle at 90°
- Wide low profile rocker bottom provides a smooth roll-through and increased stability for altered gaits
- The softgood's patented inner bootie prevents foot from migrating on the insole
- Designed around Össur's unique Active Offloading Insole
- Provides a safe environment for the simultaneous use of wound care material
- Easy to apply and modify

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

DH OFFLOADING WALKER

Size	Shoe Size US (Men)	Shoe Size US (Women)	Shoe Size EU	Item No
X-Small	2-4	3.5-5.5	33-35.5	DH0200BLK
Small	4.5-7	6-8	36-39	DH0400BLK
Medium	7.5-10.5	8.5-11.5	39.5-44	DH0600BLK
Large	10.5-12.5	11.5-13.5	44.5-46.5	DH0800BLK
X-Large	12.5+	13.5+	47-50	DH1000BLK

ACCESSORIES

REPLACEMENT PARTS

Size	Softgoods	Insoles
X-Small		11742C
Small	10958	PP04C
Medium	10959	PP06C
Large	10960	11740C
X-Large		11741C

DH OFFLOADING POST-OP STYLE SHOE

INDICATIONS FOR USE

- Plantar wounds

FEATURES AND BENEFITS

- Designed to accommodate either the left or the right foot
- Designed around Össur's patented Active Offloading Insole
- Designed to accommodate bulky dressings
- Super-soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe piece aides in wound accessibility for daily wound site inspection

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

DH POST-OP SHOE

Size	Shoe Size US (Men)	Shoe Size US (Women)	Shoe Size EU	Item No
Small	4-6	6.5-8.5	36-39	10340
Medium	6.5-8.5	9-11	39.5-44	10341
Large	9-11	12+	44.5-46.5	10342
X-Large	11.5-14	13+	47-50	10343

ACCESSORIES

REPLACEMENT PARTS

Size	Insoles
Small	12365
Medium	12366
Large	12367
X-Large	12368

PEDIATRIC WALKER

INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- During postoperative or rehab phase

FEATURES AND BENEFITS

- High quality rigid outer skeleton for protection
- Dense foam lining within the sole to provide comfort and support
- Shock absorbing sole to improve comfort and protects the user from slips, providing additional security in each step
- Uprights give additional stability by restricting motion of the leg
- Lightweight design helps in keeping kids active

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

PEDIATRIC WALKER

Size	Shoe Size US	Item No
Small	Child 6.5-9.5	PEW0400
Medium	Child 10-Youth 1	PEW0600
Large	Youth 1.5-2.5	PEW0800

REPLACEMENT LINERS

Size	Shoe Size US	Item No
Small	Child 6.5-9.5	PEW0001
Medium	Child 10-Youth	PEW0002
Large	Youth 1.5-2.5	PEW0003

METAL UPRIGHT AIR WALKER

INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- During postoperative or rehab phase

FEATURES AND BENEFITS

- Malleable struts allow you to customized fit
- Rocker bottom sole promotes a more natural gait
- Pneumatic pump allows air compression to be adjusted as needed for a more customized fit.



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's shoe size US	Women's shoe size US	Item No
X-Small	up to 4	up to 5.5	MAW0200
Small	4.5-7	6-8	MAW0400
Medium	7.5-10.5	8.5-11.5	MAW0600
Large	10.5-12.5	11.5-13.5	MAW0800
X-Large	12.5 & up	13.5 & up	MAW1000

REBOUND® ANKLE BRACE

INDICATIONS FOR USE

- Grade 1, 2, 3 and high ankle sprains
- Transition out of a walker and/or cast after fracture or surgery
- Chronic ankle instability
- Prophylactic use

FEATURES AND BENEFITS

- Adjustable calf band accommodates a broad range of leg circumferences and creates a more custom fit
- Open front, crossover strap design allows simple and secure application
- Version with stability strap offers numerous strapping configurations for a range of stability levels. Can be stepped down or removed as the patient progresses through rehabilitation



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe size US	Shoe Size UK	Shoe Size EU	Left	Right
Small	6-7.5	7.5-9	2-7.5	35-39.5	B-231602512	B-231502512
Medium	8-12	9.5-13.5	8-10.5	39.5-42	B-231602513	B-231502513
Large	12.5-16	N/A	11+	42+	B-231602514	B-231502514

ACCESSORIES

Description	Item No
Width Extender Kit	B-231900000
Stability Strap Kit	B-701160080

FORM FIT® ANKLE BRACE

The Form Fit® Ankle brace provides stability and protection for the ankle-whether used prophylactically, for chronic instability or following an injury. A hybrid between a rigid ankle stirrup and a soft ankle support, the Form Fit Ankle ensures a close fit which limits inversion, eversion, flexion and extension movements. The low profile design fits easily into shoes and the quick lace up feature makes it easy to don and doff. Breathable material ensures it can be comfortably worn all day.

The Form Fit Ankle with Figure 8 strap version, offers an even greater level of support, by further limiting the range of motion of the ankle.



INDICATIONS FOR USE

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

FEATURES AND BENEFITS

- Breathable material
- Quick lace Up
- Inversion/Eversion Control
- Available with Figure 8 strap

REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone.

ORDERING INFORMATION

Size	Ankle Circ.*	Form Fit Ankle Brace	Form Fit Ankle with Figure 8
X-Small	11-12" (28-30.5cm)	B-212010001	B-212000001
Small	12-13" (30.5-33cm)	B-212010002	B-212000002
Medium	13-14" (33-35.5cm)	B-212010003	B-212000003
Large	14-15" (35.5-38cm)	B-212010004	B-212000004
X-Large	15-16" (38-40.5cm)	B-212010005	B-212000005

*Measure ankle circumference diagonally over the base of the heel and across the ankle bone.

FORM FIT® ANKLE WITH SPEEDLACE

The Form Fit® Ankle Brace with Speedlace provides stability and protection for the ankle-whether used prophylactically, for chronic instability or following an injury. Made from durable materials, its low profile, quick-lace design is comfortable and fits easily into shoes. Removable semi-rigid stays provide additional support and the Figure 8 straps are designed for inversion/eversion control.

INDICATIONS FOR USE

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

FEATURES AND BENEFITS

- Durable materials
- Quick-Pull closure straps
- Figure 8 straps
- Removable stays
- Circumferential cuff


REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone

ORDERING INFORMATION

Size	Ankle Circ.*	Item No
XXSmall	9-10" (23cm-25.5cm)	W-10621
XSmall	10-11" (25.5cm-28cm)	W-10622
Small	11-12" (28cm-30.5cm)	W-10623
Medium	12-13" (30.5cm-33cm)	W-10625
Large	13-14" (33cm-35.5cm)	W-10627
XLarge	14-15" (35.5cm-38cm)	W-10628
XXLarge	15-16" (38cm-40.5cm)	W-10629

*Measure ankle circumference diagonally over the base of the heel and across the ankle bone.



AIRFORM® UNIVERSAL INFLATABLE ANKLE STIRRUP

INDICATIONS FOR USE

- Acute Ankle Injuries
- Post-cast healing and post-walker support
- Chronic Instability

FEATURES AND BENEFITS

- Universal design fits both left and right ankles, reducing inventory needs
- Soft Flex-Edge molding provides a comfort barrier between the shell and ankle
- Inflatable air bladders provides adjustable pneumatic compression
- Adjustable heel strap allow for easy width adjustments



ORDERING INFORMATION

Size	White Single	White 6 pack	Black Single
Youth/Pony	80251	802516	
Adult	80250	802506	80350

AIRFORM® INFLATABLE ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Inflatable bladder provides adjustable pneumatic compression
- Side specific right or left sizing



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	70200	70201
Adult	80200	80201

AIRFORM® PRE-INFLATED ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Pre-inflated, pulsating air bladder with regulated chamber reduces bottoming out and decreases edema
- No straws or over inflation
- Adjustable pivoting heel strap increases comfort and allows for full range of motion
- Side specific right or left sizing



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	70100	70101
Adult	80100	80101

FORM FIT® HONEYCOMB ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Patented Honeycomb pads
- Flex-Edges and adjustable pivoting heel straps



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	40100	40101
Adult	30100	30101

FORM FIT® FOAM ANKLE BRACE

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Patented Honeycomb pads
- Flex-Edges and adjustable pivoting heel straps



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	20100	20101
Adult	50100	50101

EXOFORM® DORSAL NIGHT SPLINT

INDICATIONS FOR USE

- Plantar fasciitis
- Post-static pain
- General heel pain
- Drop foot
- Achilles tendonitis
- Lower extremity overuse injuries

FEATURES AND BENEFITS

- Adjustable aluminum stays can be positioned to maintain flexion between 85°-90°
- Soft, Flex edge over molding improves patient comfort
- Full calf and heel strapping prevent brace from rotating
- Slip resistant padding

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Item No
Medium	5-9	6 - 10	4.5-8	37.5-42.5	W-50085
Large	9.5-14	10.5-14	8.5-13.5	43-48.5	W-50087

FORM FIT® NIGHT SPLINT

INDICATIONS FOR USE

- Plantar fasciitis
- Lower extremity overuse injuries
- Achilles tendinopathy

FEATURES AND BENEFITS

- Straps adjust from 0°-10° of dorsiflexion
- Contoured for a comfortable fit
- Soft foam interior is washable and breathable
- quick release buckles for easy application
- Slip resistant tread

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Standard
	Male	Female	Male	Female	Male	Female	
Small	Up to 7	Up to 7.5	Up to 6.5	Up to 5	Up to 39	Up to 38	50023
Medium	7.5-10	8-10.5	7-9.5	5.5-8	39.5-44	38.5-43	50025
Large	10.5 and up	11 and up	10 and up	8.5 and up	44.5 and up	43.5 and up	50027

ACCESSORIES

Description	Size	Item No
Foam Wedge		10700
Replacement Liner	Small	11760
Replacement Liner	Medium	11761
Replacement Liner	Large	11762

AIRFORM® NIGHT SPLINT

INDICATIONS FOR USE

- Plantar Fasciitis
- Achilles Tendonitis
- Drop foot
- Post-static pain

FEATURES AND BENEFITS

- Light weight, low profile design
- Easy to adjust
- Available in singles or 6 packs

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Standard	6 Pack
	Male	Female	Male	Female	Male	Female		
Medium	5-9	6-10	4.5-8.5	3.5-7.5	37-43	36-42	50045	500456
Large	9.5-14	10.5-15	9-13.5	8-12	43.5-48.5	42.5-48	50047	500476

AFO LEAF SPRING

AFO Leaf Spring is a pre-fabricated polypropylene ankle-foot orthosis designed to support flaccid drop foot. It provides a semirigid section for toe clearance and support, yet has a thin, flexible foot plate for easy trimming. The absence of a heel section makes the Leaf Spring more comfortable to wear and provides a better fit in shoes.



INDICATIONS FOR USE

- Drop foot
- Mild drop foot secondary to other neurological pathologies

FEATURES AND BENEFITS

- Lightweight injection-molded polypropylene
- Thin, flexible foot part which is easy to trim with scissors
- Absence of heel section for improved comfort

REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size US	Shoe Size UK	Shoe Size EU	Height	Left	Right	Foot Plate Length
Small-Medium	7.5-10.5	3.5-6.5	36-39	10.75" (27cm)	O-110227	O-110127	24cm
Large	11-14	7-10.5	39.5-45	12" (30cm)	O-110230	O-110130	28cm

MESH TOP POST-OP SHOE

INDICATIONS FOR USE

- For recovery after soft tissue procedures and post trauma

FEATURES AND BENEFITS

- Breathable mesh upper with reinforced heel counter and padded heel collar, more stability overall to the rear of the foot while reducing slippage
- Semi-rigid sole provides superior protection and control after osteotomies
- Extended strap and dual side closure system accommodates extra bulky dressings



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Male	Female
	Male	Female	Male	Female	Male	Female		
Small	4-6	5-6	3.5-5.5	2.5-3.5	35.5-38	35-36	17013	17003
Medium	6.5-8.5	7-8	6-8	4-5.5	38.5-42	36.5-38.5	17015	17005
Large	9-11	9 +	8.5-10.5	6 +	42.5-45	39 +	17017	17007
X-Large	11.5-14		11-13.5		45.5-48.5		17018	

SOFT TOP POST-OP SHOE

INDICATIONS FOR USE

- For recovery after soft tissue procedures and post trauma

FEATURES AND BENEFITS

- The cotton-poly blend upper aids in breathability and comfort
- The contact closure system is designed for ease of use



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Male	Female
	Male	Female	Male	Female	Male	Female		
Pediatric	Preschool 10-Youth 1						18002	
Small	4-6	5-6	3.5-5.5	2.5-3.5	35.5-38	35-36	18013	18003
Medium	6.5-8.5	7-8	6-8	4-5.5	38.5-42	36.5-38.5	18015	18005
Large	9-11	9 +	8.5-10.5	6 +	42.5-45	39 +	18017	18007
X-Large	11.5-14		11-13.5		45.5-48.5		18018	

CANVAS ROCKER BOTTOM CAST SHOE

INDICATIONS FOR USE

- For use with a short leg cast

FEATURES AND BENEFITS

- Cushioned rocker bottom cast sandal
- Contact closure, open toe



REQUIRED MEASUREMENTS

- Insole length

ORDERING INFORMATION

Size	Insole length	Blue	Black
X-Small	7" (18cm)	308XS	
Small	7.5-9" (19-23cm)	308SM	B308SM
S/Medium	9-9.75" (23-25cm)	308SM/MED	B308S/M
Medium	10-10.5" (25.5-26.5cm)	308MED	B308MED
M/Large	10.5-11" (26.5-28cm)	308MLG	B308MLG
Large	11-12" (28-30.5cm)	308LG	B308LG
X-Large	12-13.75" (30.5-35cm)	308XLG	B308XLG



SPECIALITY CARE,
CASTING & TRACTION



STOCKINETTE

- 100% synthetic stockinette
- Tubular stockinette for use with fiberglass, 25 yard dispenser box

ORDERING INFORMATION

Width	U/M	Synthetic	Cotton
1" (2.5cm)	1 ea	2101	2000A
2" (5cm)	1 ea	2102	2001
3" (7.6cm)	1 ea	2103	2002
4" (10cm)	1 ea	2104	2003
6" (15.2cm)	1 ea	2106	2005

TECHFORM® PREMIUM

INDICATIONS FOR USE

- Immobilization

FEATURES AND BENEFITS


- Glides off the roll for easier application
- Soft, conformable rolls allow for easy handling
- "Rub-responsive" resin formula adapts to the technician's technique-the more it is rubbed, the smoother the cast gets
- Excellent lamination ensures a strong and durable cast
- Smoother finished cast is non-abrasive and will not snag on clothing
- Precision welded ends do not fray and lay down easily
- Finger holes in the pouch allow for easy opening



ORDERING INFORMATION

WIDTH	WH	BK	OR	RD	PK	HP	LB	DB	PR	GN	RB
1.5" (3.8cm)	154WH										
2" (5cm)	204WH	204BK	204OR	204RD	204PK	204HP	204LB	204DB	204PR	204GN	204RB
3" (7.6cm)	304WH	304BK	304OR	304RD	304PK	304HP	304LB	304DB	304PR	304GN	304RB
4" (10cm)	404WH	404BK	404OR	404RD	404PK	404HP	404LB	404DB	404PR	404GN	404RB
5" (12.7cm)	504WH										

Note: Packed in boxes of 10



TECHFORM® ORTHO-ROLL®

Developed in collaboration with orthopaedic professionals, Techform Ortho-Roll represents the first major breakthrough in orthopaedic splinting technology in decades.

In focus groups, Techform Ortho-Roll rated from "Good" to "Excellent" in every key category, including handling, conformability, padding, water retention, dispensing method, ease-of-use, overall appearance, and overall satisfaction. Techform Ortho-Roll simply performs better.

FEATURES AND BENEFITS

- Quickdry is a patented, revolutionary orthopaedic padding material. Quickdry padding never really stays wet. Water simply passes straight through it. The result is less risk of skin maceration. Plus, it saves you application time
- Instructions are right on the box
- Box opens easily along perforation
- Closure rod creates an air tight seal



ORDERING INFORMATION

TECHFORM ORTHO-ROLL

Item No	Size	Item No	Size
TR0115	1"x15ft.	TR0415	4"x15ft.
TR0215	2"x15ft.	TR0515	5"x15ft.
TR0315	3"x15ft.	TR0615	6"x15ft.

ACCESSORIES

Item No	Size
13302	Positioning Strap (12 pack)
13332	Closure Rod
P-550606	Splinting Guide Poster

TECHFORM® PRE-CUT SPLINTS

Techform Pre-Cut Splints have the same key features as our patented Ortho-Roll splinting material, including our patented Quickdry padding. Each splint comes pre-packaged in an individual foil pouch, perfect for those who prefer to eliminate waste of splint material.



FEATURES AND BENEFITS

- Our selection includes 5 sizes, a full variety of lengths and widths to accommodate all your splinting needs.
- Pre-cut splints come in various lengths so you can choose the appropriate size to match your specific needs

ORDERING INFORMATION

Item No	Size	Quantity
TC0312	3"x12" (8x30cm)	10/Box
TC0335	3"x35" (8x90cm)	5/Box
TC0415	4"x15" (10x40cm)	10/Box
TC0430	4"x30" (10x76cm)	8/Box
TC0530	5"x30" (13x76cm)	10/Box