

live with pride

www.arrexcanada.com

www.arrexcanada.com

HEAD OFFICE
TORONTO, CANADA

BRANCH OFFICE
DEIRA, DUBAI,
KUWAIT CITY, KUWAIT

SOLE IMPORTERS, DISTRIBUTORS & WHOLESALERS FOR ASIA

TRUSTED SINCE 1920

BOMBAY SURGICAL Co.

HOSPITAL, SURGICAL & REHABILITATION EQUIPMENTS

AN ISO 9001:2008 CERTIFIED COMPANY



designs@hyangraphics.com

arrex BROCHURE® 2016

arrexTM

DEALER STAMP



HOSPITAL AIDS
WHEEL CHAIRS
WALKING AIDS
COMMODE
BODY SUPPORTS

BROCHURE© 2016 GLOBAL RELEASE



ABOUT US

ARREX Canada was established in 1982 by an experienced group of medical professionals who have dedicated their entire professional life in the pursuit of creating innovative healthcare products.

All our products have been designed in consultation with orthopaedic doctors to ensure that patient comfort and rehabilitation is given top priority. Through the distribution of high end healthcare products, ARREX has become a trusted international brand that never compromises on quality.

Since our products are high quality and affordable we have become the #1 trusted source for wheelchairs, walking aids, hospital furniture, body support, medical silicon products, home care and patient rehabilitation products.



MISSION

To provide affordable, high quality healthcare products.
Give top priority to patient care and comfort.
Provide top notch customer support throughout the life cycle of each product.



live with pride
www.arrexcanada.com



INDEX



ARREX SERIES	CATEGORY	PAGE
WHEEL CHAIR		
• VOYAGER SERIES	COMPACT ALUMINIUM WHEELCHAIR	01
• EFFICIENCY SERIES	BASIC ALUMINIUM WHEELCHAIR	03
• ZEN SERIES	PREMIUM ALUMINUM WHEELCHAIRS	04
• COMFORT SERIES	RECLINING WHEELCHAIR	08
• ELEMENT SERIES	BASIC STEEL WHEELCHAIR	11
• PRIMO SERIES	PREMIUM STEEL WHEELCHAIR	14
• TITAN SERIES	HEAVY DUTY STEEL WHEELCHAIR	16
• SPROUT SERIES	WHEELCHAIRS FOR KIDS & YOUNG ADULTS	19
• IMPERIAL SERIES	POWER WHEELCHAIR	21
• XENIA SERIES	COMMODO WHEELCHAIR	25
• CUSTOM SERIES	CUSTOM WHEELCHAIR	30
COMMODO		
• VITA SERIES	BASIC COMMODO CHAIR	31
• VITA PLUS SERIES	PREMIUM COMMODO CHAIR	33
• VITA RISE SERIES	COMMODO RAISER	37
BATH BENCHES		
• LAGOON SERIES	BATH BENCHES	37
REHAB & WALKING AIDS		
• MARATHON SERIES	LIGHT WEIGHT WALKER	41
• MARATHON CHILD SERIES	WALKERS FOR KIDS & YOUNG ADULTS	45
• MARATHON ROLLATOR SERIES	PREMIUM ROLLATOR	47
• MARATHON PHYSIO SERIES	REHAB WALKING AIDS	49
• MARATHON SUPPORT SERIES	WALKING AIDS	51
HOSPITAL BEDS		
• MONOREST SERIES	MANUAL BEDS	53
• TRITRON SERIES	ELECTRIC BED	55
• FYSIO SERIES	PHYSIO-THERAPY COUCHES	58
HOSPITAL FURNITURE & APPARATUS		
• NOAH SERIES	HOSPITAL FURNITURE	59
• ALEVA SERIES	SUCTION MACHINE	63
• PURE SERIES	OXYGEN CONCENTRATORS	65
AIR BEDS & STRETCHERS		
• STRAND SERIES	STRETCHERS	66
• AIR BEDS		72
BODY SUPPORTS & SILICONE PRODUCT		
• HAND & ELBOW		75
• LEG & ANKLE SUPPORT		78
• WAIST & BACK SUPPORT		82
• NECK SUPPORT		83
• BANDAGE		84
• SILICONE PRODUCT		85
• EXERCISER		87
• MISCELLANEOUS		88
• REHABILITATION AIDS		89
• SPARE PARTS		90
• MISCELLANEOUS		94

COMPACT ALUMINIUM WHEELCHAIRS

These FEATHERWEIGHT Aluminium wheelchairs are especially designed to help you navigate through the tight spaces of an airplane or railway compartment. All VOYAGER Series wheelchairs are highly compact and fold-up or collapse to a fraction of their size making it an ideal companion for travelling. Although lightweight, the bearing capacity of these wheelchairs is not compromised and it has the strength and durability of any of our regular models.



NANUK

Frame Type	Flip up Armrest	Wheel Type	
		Rear	Front
Aluminium	✓	Solid Mag Wheel	Solid Castor

51 (20)	38 (15)	15 (6)	12.5 (5)	45 (18)	88 (34.6)

DIM	kg	kg
46x34x66	6.9	100

AIRLIFT

Frame Type	Flip up Armrest	Wheel Type	
		Rear	Front
Aluminium	✓	Solid Mag Wheel	Solid Castor

38 (15)	32 (12.6)	15 (6)	10 (4)	45 (18)	84 (33)

DIM	kg	kg
39x29x70	6	100



TINO-X

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid PU Mag Wheel	Solid Castor

45 (18)	37 (14.6)	20 (8)	12.5 (5)	45 (18)	89 (35)

DIM	kg	kg
75x25x68	7	100



BASIC ALUMINIUM WHEELCHAIRS

ARREX designed these LIGHTWEIGHT Aluminium wheelchairs for patients who need the basic functionality of a wheelchair. These wheelchairs are full-size when open and are easily folded when not in use or travelling. Since the body of the wheelchair is made entirely of Aluminum these are lighter than steel wheelchairs.



JOEL

Frame Type	Wheel Type	
	Rear	Front
Aluminium	Solid Wheel	Solid Castor

61 (24)	41 (16)	59 (23)	20 (8)	49 (19)	67 (34)

DIM	kg	mm
103x25x89	11.7	100



FELIX

Frame Type	Detachable \ Flip up Footrest	Wheel Type	
		Rear	Front
Aluminium	✓	PU Mag Wheel	Solid Castor

56 (22)	41 (16)	30 (12)	20 (8)	49 (19)	97 (38)

DIM	kg	mm
70x26x79	14	100

For patients who need a LIGHTWEIGHT Aluminium wheelchair with a lot more functionality. The wheelchairs in these series have the premium features like flip-up armrests, detachable footrests and some have cushion seats as well.

PREMIUM ALUMINIUM WHEELCHAIRS

KONRAD-46 / KONRAD-41



Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

63 (25)	46 (18)	40 (16)	15 (6)	51 (20)	91 (36)	78x28x78	13.7	100

Frame Type	Detachable \ Flip up Armrest	Footrest	Wheel Type	
			Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

61 (24)	41 (16.1)	40 (15.7)	15 (5.9)	51 (20)	91 (35.8)	78x28x79	13.7	150

HUGO-46 / HUGO-41



Frame Type	Flip up Armrest	Wheel Type	
		Rear	Front
Aluminium	✓	Pneumatic Wheel	Solid Castor

64 (25)	46 (18)	61 (24)	15 (6)	51 (20)	94 (37)	94x28x73	14.6	100

Frame Type	Flip up Armrest	Wheel Type	
		Rear	Front
Aluminium	✓	Pneumatic Wheel	Solid Castor

59 (25.2)	41 (18.1)	61 (24)	15 (5.9)		94 (37)	94x28x73	14.6	150

ROSS



Frame Type	Flip up / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Wheel	Solid Castor

68 (26.8)	46 (18.2)	62 (24.4)	15 (5.9)	51 (20)	90 (35)

DIM	kg	mm
84x33x93	17.6	120

ZANE



Frame Type	Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

62 (24.4)	46 (18.2)	58 (22.8)	15 (5.9)		89 (35)

DIM	kg	mm
80x28x90	13.6	120

OMEGA



Frame Type	Flip Up / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

62 (24.4)	46 (18.2)	38 (14.6)	15 (5.94)	51 (20)	89 (35)

DIM	kg	mm
78x28x76	13.9	120

ARROW



Frame Type	Fixed / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

57 (22.4)	43 (17)	30 (11.8)	20 (8.4)	52 (20.4)	93 (36.6)

DIM	kg	mm
74x26x83	12	120

FABER



Frame Type	Fixed / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

57 (22.4)	43 (17)	30 (11.8)	20 (8.4)	52 (20.4)	93 (36.6)

DIM	kg	mm
74x26x83	12	120

DELUX



Frame Type	Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Mag Wheel	Solid Castor

64 (25.2)	46 (18.1)	62 (24.4)	15 (5.9)	52 (20.4)	90 (35.4)

DIM	kg	mm
80x28x91	14.2	120



LANZO

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Pneumatic Mag Wheel	Solid Castor

64 (24)	46 (18)	40 (16)	15 (6)	50 (20)	92 (36)

78x33x93	14.4	100

DETACHABLE FOOTREST
CALF STRAP



MAXIMA

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Quick release PU Mag Wheel	Solid Castor

67 (26)	46 (18)	61 (24)	15 (6)	51 (20)	91 (36)

83x33x93	16.7	150

When a patient spends a lot of time on their wheelchair, they need to have the ability to rest or even sleep in it. Our COMFORT+ wheelchairs are designed like a first class seat in an airplane- where a patient can easily take a nap and even flatten the chair completely to sleep comfortably.

RECLINING WHEELCHAIRS



JANA

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Powder Coating Steel	✓	✓	✓	PU Mag Wheel	Tpr Castor

53 (21)	46 (18)	30 (12)	20 (8)	53 (21)	118 (46)

82x61x84	35.5	150



FABIO

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Powder Coating Steel	✓	✓	✓	Pneumatic Wheel	Tpr Castor

68 (27)	46 (18)	62 (24)	20 (8)	58 (23)	121-133 (47-52)

88x67x84	37.7	150





ALMOS-46

Frame Type	Detachable \ Flip up Armrest	Footrest	Reclining	Wheel Type Rear	Front
Aluminium	✓	✓	✓	Pneumatic Wheel	Solid Castor

68 (27)	46 (18)	62 (24)	20 (8)	51 (20)	128 (50)

DIM	18	100

ALMOS-41

Frame Type	Detachable \ Flip up Armrest	Footrest	Reclining	Wheel Type Rear	Front
Aluminium	✓	✓	✓	Pneumatic Wheel	Solid Castor

63 (24.8)	41 (16.1)	62 (24)	20 (8)	51 (20)	128 (50)

DIM	18	100



DORIS

Frame Type	Detachable \ Flip up Armrest	Footrest	Reclining	Wheel Type Rear	Front
Powder Coating Steel	✓	✓	✓	Solid Wheel	Solid Castor

65 (26)	41 (16)	58 (23)	20 (8)	50 (20)	119 (47)

DIM	24.3	100



PEARCE

Frame Type	Detachable \ Flip up Armrest	Footrest	Reclining	Wheel Type Rear	Front
Aluminium	✓	✓	✓	Solid Mag Wheel	Solid Castor

66 (26)	44 (17.3)	58 (22.8)	20 (7.9)		86 (33.6)

DIM	18	120

BASIC STEEL WHEELCHAIRS

Our most affordable wheelchairs till date. Although low-cost we have maintained the integrity of our ELEMENT series wheelchairs, ensuring that it can easily bear a full-size grown adult and is rugged enough to endure everyday abuses.



TOMMY

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Wheel	Solid Castor

64 (25)	46 (18)	58 (23)	20 (8)	49 (19)	87 (34)

DIM	17.1	100



TOMMY LIFT

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Wheel	Solid Castor

64 (25)	46 (18)	58 (23)	20 (8)	49 (19)	87 (34)

DIM	17.5	100



TOMMY MAG

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Mag Wheel	Solid Castor

65 (26)	46 (18)	62 (24)	20 (8)	49 (19)	89 (35)

DIM	17.8	100



TOMMY PLUS

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Mag Wheel	Solid Castor

65 (26)	46 (18)	62 (24)	20 (8)	49 (19)	89 (35)

DIM	17.8	100

ROSA



Frame Type	Wheel Type	
	Rear	Front
Powder Coating Steel	Solid Wheel	Solid Castor

62 (24)	46 (18)	59 (23)	20 (8)	49 (19)	87 (34)

DIM	16.2	100

LITT



Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Wheel	Solid Castor

64 (25)	46 (18)	58 (23)	20 (8)	49 (19)	87 (34)

DIM	17.1	100

CRUISER



Frame Type	Detachable \ Flip up Footrest	Wheel Type	
		Rear	Front
Chromed Steel	✓	Solid Wheel	Solid Castor

56 (22)	43 (17)	20 (8)	20 (8)	50 (20)	96 (38)

DIM	12.9	100

Generous amount of premium features like detachable armrest, elevating footrest and heavy-duty castors, make our PRIMO series wheelchairs the best choice for patients who want luxury at an affordable cost.

PREMIUM STEEL WHEELCHAIRS

MARTIN



Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Solid Mag Wheel	PU Castor

64 (25)	46 (18)	61 (24)	20 (8)	50 (20)	88.5 (35)

DIM	17.9	100

SAVINA



Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Chromed Steel	✓	✓	Solid Wheel	PU Castor

65 (26)	46 (18)	62 (24)	20 (8)	51 (20)	89 (35)

DIM	19.6	100

PERLA LIFT



Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Pneumatic Mag Wheel	Solid Castor

62 (24)	46 (18)	31 (12)	20 (8)	50 (20)	88 (35)

DIM	kg	kg
71x25x90	14.3	75

PERLA



Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Chromed Steel	✓	✓	Pneumatic Wheel	PU Castor

65 (26)	46 (18)	62 (24)	20 (8)	51 (20)	89 (35)

DIM	kg	kg
80x28x91	18.4	100

CERSEI



Frame Type	Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front

62 (24)	46 (18)	8 (3.14)	8 (3.14)	50 (20)	88 (35)

DIM	kg	kg
71x25x90	14.3	75

Our double reinforced TITAN series wheelchairs are the only choice for patients who need a little bit of more room. The heavy duty castors and wheels can bear loads up to 125 kgs and the robust frame can also endure its share of stress.

HEAVY DUTY STEEL WHEELCHAIRS



VICTOR 51

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Pneumatic Mag Wheel	Rubber Castor

64 (25.1)	51 (20)	62 (24)	20 (8)	54 (21)	93 (37)

DIM	kg	kg
83x33x93	21.5	150

VICTOR 56

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Pneumatic Mag Wheel	Rubber Castor

69 (27.1)	56 (22)	62 (24)	20 (8)	54 (21)	93 (37)

DIM	kg	kg
83x33x93	22	150

VICTOR 60

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Pneumatic Mag Wheel	Rubber Castor

73 (29)	60 (24)	62 (24)	20 (8)	54 (21)	93 (37)

DIM	kg	kg
83x33x93	22	150



HARVEY

Frame Type	Flip up / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Wheel	Solid Castor

68 (26.7)	51 (20)	62 (25.4)	15 (5.9)	50 (19.6)	89 (35)

DIM	kg	mm
83x30x90	18	180



MOLLY

Frame Type	Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Mag Wheel	Solid Castor

78 (31)	61 (24)	30 (12)	20 (8)	50 (20)	96 (38)

DIM	kg	mm
80x32x83	19.3	180



TERZO

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Pneumatic Mag Wheel	Solid Castor

78 (31)	61 (24)	30 (12)	20 (8)	50 (20)	96 (38)

DIM	kg	mm
80x32x83	19.3	180



SKY

Frame Type	Flip up / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Wheel	Mag Wheel

68 (26.7)	51 (20)	62 (25.4)	15 (5.9)	50 (19.6)	89 (35)

DIM	kg	mm
83x30x90	18	180



GUSTAF - 56 / GUSTAF -51

Frame Type	Flip up / Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Chromed Steel	✓	✓	Solid Wheel	Solid Castor

73 (28.7)	56 (22)	62 (25.4)	15 (5.9)	50 (19.6)	89 (35)

DIM	kg	mm
83x30x90	18	180

68 (26.7)	51 (20)	62 (25.4)	15 (5.9)	50 (19.6)	89 (35)

DIM	kg	mm
83x30x90	18	180



KELSO

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Pneumatic Wheel	Solid Castor

81 (32)	61 (24)	62 (24)	20 (8)	51 (20)	93 (37)

DIM	kg	mm
86x37x83	21.8	125



WHEELCHAIRS FOR KIDS & YOUNG ADULTS

Just like a child is not comfortable wearing an adult's shirt, in the same way a child is never at home using an adult wheelchair. Since our SPROUT series wheelchairs are specifically designed for young adults and kids, it is far more easy to maneuver and handle. Special children have special needs and ARREX knows how to cater to them.



ENU

Frame Type	Flip up Armrest	Wheel Type	
		Rear	Front
Aluminium	✓	Solid Wheel	Solid Castor

52 (20)	35 (14)	58 (23)	15 (6)	48 (19)	90 (35)

DIM	13	75



ELMO

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Wheel	Solid Castor

55 (22)	35 (14)	55 (22)	15 (6)	47 (18.5)	88 (35)

DIM	16	100



NONA-X

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Aluminium	✓	✓	✓	PU Mag Wheel	Solid Castor

52 (20)	38 (15)	40 (16)	15 (6)	49 (19)	97-114 (38-45)

DIM	22	75



TINI

Frame Type	Reclining	Wheel Type	
		Rear	Front
Aluminium	✓	Solid Mag Wheel	Solid Castor

42 (16.5)	30 (11.8)	30 (11.8)	15 (6)	49 (19)

DIM	21	50



CLARA

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Chromed Steel	✓	✓	✓	Solid Mag Wheel	Solid Castor

55 (21.6)	41 (16)	40 (16)	15 (6)	49 (19)	107 (42)

DIM	20.2	100

POWER WHEELCHAIRS

The freedom and independence that a patient experiences when they can get from one place to another without the need of an attendant or the fear of injuring their own hands is unparalleled. IMPERIAL series power/electric wheelchairs come standard with a high charge battery that can operate continuously for 20 km.

QADIR

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Rubber Foam	Rubber Castor

--	--	--	--	--	--

--	--	--



TAMAR

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Rubber Foam	Rubber Castor

--	--	--	--	--	--

--	--	--



ZULA

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Aluminium & Powder Coating	✓	✓	✓	Rubber Foam	Rubber Castor

--	--	--	--	--	--

--	--	--



HERMOS

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Resin Coating Steel	✓	✓	Rubber Foam	Rubber Foam

--	--	--	--	--	--

--	--	--



THOR

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓		Strong Rubber Foam	Strong Rubber Castor

--	--	--	--	--	--

--	--	--



SNOW

Frame Type	Flip up Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Wheels	Solid Castor

63 (24.8)	41 (16.1)	62 (24.4)	17 (6.6)	49 (19.2)	95 (37.4)

DIM	kg	mm
91x44x86	33	125

SPECTER

Frame Type	Flip up Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Wheels	Solid Castor

63 (24.8)	41 (16.1)	62 (24.4)	17 (6.6)	49 (19.2)	95 (37.4)

DIM	kg	mm
91x44x86	33	125

JABBAR

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Rubber Foam	Rubber Castor

64 (25)	43 (17)	32 (13)	20 (8)	56 (22)	97 (38)

DIM	kg	mm
81x69x72	62.8	110



HECTOR

Frame Type	Flip up Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Solid Wheels	Solid Castor

63 (24.8)	41 (16.1)	62 (24.4)	17 (6.6)	49 (19.2)	95 (37.4)

DIM	kg	mm
91x44x86	33	125

63 (24.8)	41 (16.1)	62 (24.4)	17 (6.6)	49 (19.2)	95 (37.4)

DIM	kg	mm
91x44x86	33	125

YALE

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Cushioned Seat & Metal body	✓		Rubber Foam	Rubber Foam

63 (25)	46 (18)	26 (10)	26 (10)	64-70 (25-27.8)	135 (53)

DIM	kg	mm
140x71x79	105.2	110



COMMODE WHEELCHAIRS

A growing number of our customers were demanding a solution where a commode and wheelchair are combined so as to save on storage space and for other ergonomic reasons. Our XENIA series wheelchairs bridges the gap between a fully functional wheelchair, replete with features, and a portable commode.



RENZO

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Chromed Steel	✓	✓	✓	Pneumatic Wheel	PU Castor

63 (25)	46 (18)	61 (24)	20 (8)	53 (21)	101 (40)	81x64x65	26	130	



MORSO

Frame Type	Fixed		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Pneumatic Wheel	Solid Castor

64 (25)	44 (17)	58 (23)	20 (8)	50 (20)	87 (34)	92x29x89	19	100	



MORSO BASIC

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	Solid Castor With Brakes	Solid Castor

64 (25)	44 (17)	58 (23)	20 (8)	50 (20)	87 (34)	92x29x89	19	100	



MORSO PLUS

Frame Type	Wheel Type	
	Rear	Front
Aluminium	Pneumatic Wheel	Solid Castor

68.5 (27)	46 (18)	58 (23)	15 (6)	51 (20)	90 (35)	94x29x92	15.2	100	



MORSO COMFORT

Frame Type	Detachable \ Flip up		Reclining	Wheel Type	
	Armrest	Footrest		Rear	Front
Chromed Steel	✓	✓	✓	Solid Wheel	Solid Castor

64 (25)	44 (17)	58 (23)	20 (8)	53 (21)	129 (51)

107x33x78	25.6	100



EDDA / EDDA-51

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Powder Coating Steel	✓	✓	Solid Wheel With Lock	Solid Castor

55 (22)	44 (17)	12.5 (5)	12.5 (5)	51 (20)	81 (32)	54x53x53	12	100

62 (25.4)	51 (20)	12.5 (5)	12.5 (5)	51 (20)	81 (32)	54x53x53	12	100



EDDA-S

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Stainless Steel	✓	✓	Castor Wheel With Lock	Castor Wheel

55 (22)	44 (17)	12.5 (5)	12.5 (5)	51 (20)	88 (35)

54x53x56	10.5	100



EZIO

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Chromed Steel	✓	✓	Solid Wheel	Solid Wheel

60 (24)	45 (18)	58 (23)	15 (6)	53 (21)	90 (35)

75x28x91	19.6	100



ZURI

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Aluminium	✓	✓	Pneumatic Wheel	Solid Castors

61 (24)	46 (18.1)	58 (22.8)	15 (5.9)	48 (18.8)	93 (36.6)

80x35x71	15	100

HEAVY DUTY COMMODE WHEEL- CHAIRS



GOLIATH

Frame Type	Wheel Type	
	Rear	Front
Chromed Steel	PU Wheel	PU Castor

89 (35)	70 (28)	59 (23)	20 (8)	70 (28)	105 (41)

DIM	kg	in
112x35x115	34	150

HEAVY DUTY COMMODE WHEEL- CHAIRS



SAMSON

Frame Type	Detachable \ Flip up		Wheel Type	
	Armrest	Footrest	Rear	Front
Chromed Steel	✓	✓	Solid Wheel	PVC Castor

70 (28)	51 (20)	59 (23)	20 (8)	52 (20)	88 (35)

DIM	kg	in
78x36x89	25	100

From time to time we have customers whose specific needs are not met with our standard wheelchairs then we go back to our drawing board to make a customized product. Our CUSTOM series showcases a couple of our past projects. Please contact us with your specific needs and we will customize a wheelchair specially for you.

CUSTOM WHEELCHAIRS



TRANSPORTER

Frame Type	Wheel Type	
	Rear	Front
Powder Coating	Pneumatic Mag Wheel	Solid Castor

65 (26)	46 (18)	40 (16)	20 (8)	49 (19)	89 (35)

DIM	kg	in
80x30x76	15	150



SCANNER- MRI WHEELCHAIR

Frame Type	Detachable		Wheel Type	
	Armrest	Footrest	Rear	Front
Virgin Plastic	✓	✓	Mag Wheel	Solid Castor

64 (25)	46 (18)	58 (23)	20 (8)	49 (19)	87 (34)

DIM	kg	in
94x23x89	17.1	100

BASIC COMMODE CHAIRS

Our VITA series commode chairs provide the basic functionality of portability and convenience. All the chairs that you see in this series are foldable and lightweight, thus making them perfect for travelling.



V10

Frame Type	FOLDABLE
Powder Coating Steel	✓

41(16)	38(15)	45(18)	40(16)	67x20x39	(CM)	N.WT(kg)
						4.7



V30

Frame Type	FOLDABLE	POT ATTACHED
Powder Coating Steel	✓	✓

60(24)	48(19)	45(18)	62(32)	41(16)	40(16)	

(CM)	N.WT(kg)
88x30x60	12.8



V20

Frame Type	FOLDABLE	POT ATTACHED
Powder Coating Steel	✓	✓

48(19)	43(17)	44(17)	78(31)	40(16)	35(14)	

(CM)	N.WT(kg)
79x43x55	16.7



V50

Frame Type	FOLDABLE	POT ATTACHED
Powder Coating Steel	✓	✓

45(17.7)	43(16.9)	42(16.5)	77(30.3)	40(15.7)	28(11)	

(CM)	N.WT(kg)
63x32x44	8.9

PREMIUM COMMODE CHAIRS

Our VITA PLUS series premium commode chairs would be perfect for patients who need features like height adjustability, armrests and the option of using their commode chairs as bath benches.



VP10

Frame Type	HEIGHT ADJUSTABLE	POT ATTACHED
Powder Coated	✓	✓

48(19)	48(19)	37-47(15-18.5)	71-80(28-31.5)	44(17)	33(13)
--------	--------	----------------	----------------	--------	--------

(CM)	N.WT(kg)
56x40x36	6



VP20 STEEL

Frame Type	FOLDABLE	HEIGHT ADJUSTABLE	POT ATTACHED
Chromed Steel	✓	✓	✓

48(19)	45(18)	45-55(18-22)	78-88(31-35)	43(17)	33(13)
--------	--------	--------------	--------------	--------	--------

(CM)	N.WT(kg)
50x21x76	6



VP20 ALUMINIUM

Frame Type	FOLDABLE	HEIGHT ADJUSTABLE	POT ATTACHED
Aluminium	✓	✓	✓

48(19)	45(18)	45-55(18-22)	78-88(31-35)	43(17)	33(13)	50x23x76	3.5
--------	--------	--------------	--------------	--------	--------	----------	-----



VP30

Frame Type	FOLDABLE	HEIGHT ADJUSTABLE	POT ATTACHED
Chromed Steel	✓	✓	✓

52(20)	44(17)	45-55(18-22)	77-87(30-34)	43(17)	32(13)	52x23x79	6.2
--------	--------	--------------	--------------	--------	--------	----------	-----



VP40 STEEL

Frame Type	Wheel Type		HEIGHT ADJUSTABLE	POT ATTACHED
	Rear	Front		
Chromed Steel	Solid Castor With Brakes	Solid Castor	✓	✓

55(22)	45(18)	7.5(3)	7.5(3)	45-55(18-22)	78-88(31-35)	50x21x76	7	100
--------	--------	--------	--------	--------------	--------------	----------	---	-----



VP40 ALUMINIUM

Frame Type	Wheel Type		HEIGHT ADJUSTABLE	POT ATTACHED
	Rear	Front		
Aluminium	Solid Castor With Brakes	Solid Castor	✓	✓

55 (22)	45 (18)	7.5 (3)	7.5 (3)	45-55 (18-22)	75-88 (31-35)

50x21x76	7	150



VP60

Frame Type	FOLDABLE	POT ATTACHED
Aluminium	✓	✓

60 (23.6)	48 (18.8)	45 (17.7)	82 (32.2)	41 (16.1)	41 (16.1)	80x30x60	6.4



VP50

Frame Type	FOLDABLE	HEIGHT ADJUSTABLE	POT ATTACHED
Aluminium	✓	✓	✓

60 (23.6)	48 (18.8)	45 (17.7)	82 (32.2)	41 (16.1)	41 (16.1)

80x30x60	6.4



VP70

Frame Type	Wheel Type		HEIGHT ADJUSTABLE	POT ATTACHED
	Rear	Front		
Powder Coated	Solid Castor With Brakes	Solid Castor	✓	

55 (22)	45 (18)	12	12	45-55 (18-22)	82-90	50x21x76	9	180



VP90 / VP90-51

Frame Type	FOLDABLE	POT ATTACHED
Powder Coated	✓	✓

51 ()	44	43-58	76-91	40	33	41x58x43	12.3

58 ()	51	43-58	76-91	40	33	41x58x43	12.3



New
Outer Pack



CLAMP ATTACHES
TO THE CERAMIC POT

VA20
2 INCH COMMODE RAISERS

VA40
4 INCH COMMODE RAISERS

VA60
6 INCH COMMODE RAISERS

Commode Raiser are used on to raise height western commode

BATH BENCHES

Slip and falls are one of the most hazardous form of accidents for recuperating patients; our LAGOON series bath benches and chairs provide security for the patient and peace of mind for the caregiver.



LA10

42 (17)	50 (20)	35-45 (14-18)	74-84 (29-33)	30 (12)	40 (16)	68x19x39	6		



LA30

59 (23)	48 (19)	50 (20)	68 (35)	37 (15)	39 (15)	102x20x54	6.5		



LA40

39 (15)	51 (20)	47 (19)	52 (20)	30 (12)	74x19x52	6.4			



LA50

38 (15)	48 (19)	51-63 (20-25)	69-81 (27-32)	28 (11)	55x41x78	5.4			



LA60

42 (17)	50 (20)	35-45 (14-18)	39-49 (15-19)	30 (12)	60x17x36.5	3.5			



LA80

41 (16)	51 (20)	39-49 (15-19)	76-86 (30-34)	37 (15)	37 (15)	70x56x41.5	8.5		



LA90

27 (11)	41 (16)	14 (6)	23 (9)	14 (6)	45x27x19 (CM)	N.WT(kg)	



LA100

51 (20)	41 (16)	29-14 (11-6)	29 (11)	23 (9)	50x42x46 (CM)	N.WT(kg)	6.4



LA130

42 (17)	50 (20)	42-52 (16.5-20.5)	45-55 (17.3-21.7)	30 (12)	50x43x45 (CM)	N.WT(kg)	2.6



LA150

51 (20)	41 (16.1)	34-44 (13.4-17.3)	73-83 (28.7-32.7)	37 (14.6)	41 (16.1)	53x50x25 (CM)	N.WT(kg) 6.6



LA110

42 (16.5)	52 (20.4)	42-52 (16.5-20.4)	75-85 (29.5-33.4)	30 (11.8)	40 (15.7)	50x46x45 (CM)	N.WT(kg) 3.5



LA120

41 (16.1)	50 (20)	47 (18.5)	80 (31.5)	30 (12)	32 (12.6)	90x22x51 (CM)	N.WT(kg) 8.8



LA160

42 (16.5)	52 (20.4)	42-52 (16.5-20.4)	75-85 (29.5-33.4)	30 (11.8)	40 (15.7)	50x46x45 (CM)	N.WT(kg) 3.5



LA170

23 (9)	45 (17.7)	18 (7)	23 (9)	83x25x40 (CM)	N.WT(kg)	1.6	

LIGHT WEIGHT WALKERS



M10

50 (20)	47 (19)	52-62 (20-24)	79-89 (31-35)	60x34x89.5 (CM)	N.WT(kg)	7



M20

49 (19)	44 (17)	43-61 (17-24)	74-92 (29-36)	57x19x79 (CM)	N.WT(kg)	6.7



M30

50 (20)	44 (17)	78-96 (31-38)	58x18x82 (CM)	N.WT(kg)	5.5	



M40

48 (19)	44 (17)	82-100 (32-39)	57x19x86 (CM)	N.WT(kg)	5.1	



M50

47 (19)	46 (18)	78-96 (31-38)	53x17x81 (CM)	N.WT(kg)	4.8	



M60

47 (19)	46 (18)	78-96 (31-38)	53x17x81 (CM)		4.8 N.WT(kg)



M70

47 (19)	46 (18.1)	77-95 (30.3-37.4)	58x18x81 (CM)		5.2 N.WT(kg)



M80-S

43 (16.9)	38 (14.9)	65-75 (25.5-29.5)	76x68x35 (CM)		14.9 N.WT(kg)



M80-M

44 (17.3)	38 (14.9)	75-85 (29.5-33.4)	86x68x35 (CM)		15.2 N.WT(kg)



M80-L

45 (17.7)	38 (14.9)	85-95 (33.4-37.4)	96x68x35 (CM)		16 N.WT(kg)



M90


82 (32.2)	41 (16.1)	78-98 (30.7-38.5)	76x19x96 (CM)		9.4 N.WT(kg)

WALKERS FOR KIDS & YOUNG ADULTS



MC10

56 (22)	35 (14)	45-55 (18-22)	66x18x65 (CM)	6 N.WT(kg)



MC40

55 (21.6)	35 (13.7)	41-51 (16.1-20)	59x30x52 (CM)	5.2 N.WT(kg)



MC20

54 (21)	45 (18)	52-87 (20-34)	62x22x67 (CM)	8 N.WT(kg)



MC50

39	62	51-61	52.5x22x63 (CM)	3.9 N.WT(kg)

PREMIUM ROLLATORS

MR10

68 (27)	46 (18)	53 (21)	78-90 (31-35)	62x24x84 (CM)	7.4 N.WT(kg)



MR20

59 (23)	38 (15)	82-97 (32-38)	63x25x70 (CM)	7.7 N.WT(kg)	



MR30

68 (27)	46 (18)	53 (21)	78-90 (31-35)	62x24x84 (CM)	7.4 N.WT(kg)

MR40

60 (23.6)	47 (18.5)	47 (18.5)	76-86 (29.9-33.8)	70x20x75 (CM)	10.3 N.WT(kg)





MR50

68 (27)	46 (18)	53 (21)	78-90 (31-35)	62x24x84 (CM)	N.WT(kg)	7.4



MP20

74 (29)	45 (18)	55 (22)	110-155 (43-61)	74x15x85 (CM)	N.WT(kg)	11.9

REHAB WALKING AIDS



MP10

- Length: 65cm, Weight: 18kg, Max load: 120 kg
- Adjustable- Grip angle
- Armrest width adjusted from 53cm-76cm
- Armrest height adjusted from 71cm-110cm
- Leg width adjusted from 64cm-70cm
- Strong supported by gas spring
- Castors with brake.



MP30

35 (13.7)	34 (13.4)	75-82 (29.5-32.2)	43x14.5x76.5 (CM)	N.WT(kg)	3.2	



WALKING AIDS



MS10 MS20 MS30 MS40 MS50 (BROWN) MS70/MS70-X/Y FOLDING STICK MS210 MS240 MS60 MS250

HEIGHT ADJUSTABLE HEIGHT ADJUSTABLE HEIGHT ADJUSTABLE HEIGHT ADJUSTABLE



MS140 ELBOW CRUTCH MS150 (L/M/S) UNDERARM CRUTCH MS160 (L/M/S) UNDERARM CRUTCH MS270 (L/M/S) UNDERARM CRUTCH MS-VI STICK FOR THE VISUALLY IMPAIRED MS-PU PICK UP STICK MS230



MS80 MS90 AL MS90 STEEL MS110 MS120 MS200 AL



MS180




MS190



MS220

MANUAL BEDS


Monorest Series beds are the most economical hospital beds offered by Arrex. We have designed the manual hospital beds with heavy-duty steel and height, back and leg adjustability. All the features of the bed can be easily operated with a simple turn of a lever, placed at the feet of the patient. The bed offers you hospital like comfort in your house.



UNO 2 FUNCTION MANUAL BED

- H: 24 Inch with mattress
- Weight: 90kg
- Manual Bed With Two Function, Side Railings, Saline Rod, Mattress & Urine Hook.

Max. 175 kg	0~80°	0~40°




AMADO 3 FUNCTION MANUAL BED

- Weight: 100kg
- Manual Bed With Three Function, Side Railings, Saline Rod, Mattress & Urine Hook.

Max. 125 kg	45-78 cm	0~80°	0~40°


Tritron Series are automatic hospital beds designed with height, back and leg adjustability that can be controlled just by a click on the remote, making it easier for the patient and caretaker to operate the bed. With TRITRON Series Electronic Beds, ARREX offers you 5 star nursing comforts at an affordable price.



THEO 2 FUNCTION ELECTRONIC BED

- Weight: 110kg
- Manual Bed With Three Function, Side Railings, Saline Rod, Mattress & Urine Hook.


Max. 175 kg	0~80°	0~40°



RYERSON 3 FUNCTION ELECTRONIC BED

- Weight: 110kg
- Manual Bed With Three Function, Side Railings, Saline Rod, Mattress & Urine Hook.

Max. 125 kg	45-78 cm	0~80°	0~40°



EXCELSIOR 5 FUNCTION ELECTRONIC BED

Size: 2160x1000x480/740mm

Standard configuration:


- Advanced medical motor system (4 pcs motors, 1 pc control box, 1pc handsets)
- Low stable support system • Double side castor
- Mattress, Saline Stand, Urine hook.

N.W:120KG



Fysio Series offers physiotherapy couches that serves your patient with all the comfort, safety and helps enhance your patients experience during treatment. We have designed the couches with height, back, leg and armrest adjustability features. A perfect companion for a patient going through rehab.

PHYSIOTHERAPY COUCHES



FW MASSAGE TABLE

- Length: 6ft
- Self Weight: 25kg
- Dynamic max load: 180 kg
- Wooden Frame
- Folding & easy to carry



FX

- Length: 210cm
- Dynamic max load: 200 kg
- Static max load: 250kg
- Adjustable height from 46 cm to 95 cm raised by actuator system
- Two main pillows supported and angle-adjusted by gas springs
- Wheel lifting (Central wheel lift, Dual wheel lift)



FY

- Length: 200cm
- Dynamic max load: 200 kg • Static max load: 250kg
- Adjustable height from 46 cm to 95 cm raised by actuator system
- Three main pillows supported and angle-adjusted by gas springs
- Armrests height and direction adjustable (optional parts)
- Wheel lifting (Central wheel lift, Dual wheel lift)



FZ



- Length: 200cm
- Dynamic max load: 200 kg
- Static max load: 250 kg
- Pillow configurations
- 50cm-3cm-section; 70cm-3-section,
- 50cm-3cm-section; 50cm-4-section,
- 70cm-5-section; 70cm-2 section(Physiotherapy)
- Adjustable height from 46cm to 95cm raised by actuator system
- Three main pillows supported and angle adjusted by gas springs
- Armrests height and direction adjustable (optional parts)
- Wheel lifting (Central wheel lift, dual wheel lift)



FOOD TRAY

Attaches on Side Railings of any Hospital bed
Size: L: 99cm W: 30cm



Arrex brings you high-quality medical and hospital furniture, designed to meet and exceed the facilities provided at a hospital. Our furniture enhances patient care and creates a comfortable healing environment. Also, our home care furniture maximizes functionality without sacrificing comfort and aesthetics.

HOSPITAL FURNITURE



KINGSTON

Examination Table

Length: 6.5ft
Width: 3.5ft
Height: 2ft



NOAH TR

- Treatment Trolley



NOAH BL

- Bedside Locker
- Comes with Rotating Wheels



ORIGINS CHAIR

- Transfusion Chair



NOAH ATTEND CHAIR

- Patient Attend Chair



NOAH OBT 10

- Over bed table with Height adjustable



NOAH SS 10

- Saline Stand
- Size: H 136-230cm



NOAH SS 20

- Saline Stand
- Size: H 133-223cm



NOAH CHAIR

- Suitable for home care



NOAH STOOL

- Doctor Stool



NOAH X1

- Premium Quality.
- Includes concealed bright lights
- Vertical Mounting
- Single X-Ray View Box
- Size 50cmx15cmx67cm
- N.Wt 6kg, G.Wt 7.5kg.



NOAH X2

- Premium Quality.
- Includes concealed bright lights
- Horizontal Mounting
- Double X-Ray View Box
- Size 86cmx15cmx67cm
- N.Wt 8.5kg, G.Wt 11.5kg.



SUCTION MACHINE

Aleva Series offer compact medical suctioning devices, designed for reliable, portable operations. Our suction machines are equipped with high performance features, easy to read gauges, adjustable flow regulators and heavy-duty compressors. These machines are also light weight, easy to store and transport.



AS - HOSPITAL

- Heavy duty compressor.
- Double bottle with pedal switch.
- Suitable for operation theater purpose.
- Volume capacity – 5 liters.



AS - HOME

- Voltage: 220V+22V, 50Hz+1Hz
- Power: 120VA
- Pumping rate: 30L/min
- Continuous Working Time: 0-30min
- Bottle volume: 2500ml, Glass Bottle- 2pcs
- Ideal for use at home or Hospital
- Lightweight Compressor & Portable with wheels.
- N. wt.: 11.5kg



AS - TABLE TOP

- Voltage: AC220V+10%, 50Hz+2%
- Input power: 90VA
- Pumping rate: 18L/min
- Continuous Working Time: 0-30min
- Rest Time> 30min
- Bottle volume: 1000ml
- Ideal for use at home
- N. wt.: 4.32kg



AS - PORTABLE

- Voltage: AC220V+10%, 50Hz+2%, 12VDC
- Input power: 90VA
- Pumping rate: 18L/min
- Continuous Working Time: 0-30min
- Rest Time> 30min
- Bottle volume: 1000ml
- Comes with Car Charger
- Ideal for use in ambulance
- N. wt.: 6.7kg



AS - MANUAL

- Manually operable.
- Volume capacity – 1 liters.

OXYGEN CONCENTRATORS

Pure Series offers features home care providers want and patients deserve. Oxygen concentrators are commonly used for patients with respiratory disorders and are the best replacements for oxygen cylinders. The concentrators are light weight, low maintenance, quiet and ideal for hospital and personal use.



POX - HOME

- Flow range from 1 – 5 liter / min
- Oxygen concentration: 95.5% - 87%.
- New improved noise less compressor.
- Can be operated for 24 hours nonstop.
- Comes with LOW PURITY ALARM
- SD Card compatible
- Real Time Oxygen concentration display
- Multi-figures Display
- Large color LCD Screen
- Build in Clinical Software for nursing home & Doctors




POX - PORTABLE

- Portable: 3kg with carry bag
- 2x Batteries: Up to 6 hours working time
- Accurate pulse oxygen supply
- 12V DC power
- Rated equivalent flow- 5l/min

Strand Series displays a wide range of stretchers offered by Arrex. We have designed the stretchers with heavy-duty steel frame and patients both safety and comfort has been taken into special consideration. The stretchers are appropriate for hospital and personal use.

STRETCHERS SERIES



STR 10 HOSPITAL TROLLEY

- Size (L 6ft x W 2 ft x H 22in) :
- High position :190 x 55 x 86cm
- Low position :190 x 62 x 54cm
- Size of the castore : Ø 160mm
- Comes with Mattress
- load Bearing : 159kg
- Weight : 85kg




STR 20 ALUMINIUM ALLOY STRETCHER FOR AMBULANCE CAR

- Weight Bearing: 159Kg
- Size (LxWxH): 197 x 55 x 93
- High position: 190 x 55 x 90cm
- Low position: 190 x 60 x 24cm
- Maximum Angle of the Back : 60°
- Weight: 40kg



STR 30 ALUMINIUM ALLOY FOLDAWAY STRETCHER

- Comes with a strong Safety Bag
- Unfolded Size (LxWxH): 215 x 55 x 15cm
- Folded Size (LxWxH): 108 x 17 x 9cm
- Self-weight: 6kg
- Load Bearing: 159kg



STR 40 ALUMINIUM ALLOY FOLDAWAY STRETCHER WITH WHEELS

- Unfolded Size (LxWxH): 185 x 52 x 18cm
- Folded Size (LxWxH): 94 x 52 x 9cm
- Self-weight: 6kg
- Load Bearing: 159kg




STR 60 FOLDABLE CLOTH STRETCHER

- Load Bearing: 120 kg
- Completely Foldable
- Comes with a bag



STR 70 ALUMINIUM ALLOY SCOOP STRETCHER

- Max. unfolded Size (LxWxH) : 202 x 42 x 6cm
- Min. unfolded Size (LxWxH) : 168 x 42 x 6cm
- Can be dismantled into two pieces
- Self-weight : 10kg
- Load Bearing : 159kg



STR 80 SPINE BOARD

- High strength engineering plastics.
- It can easily float on water.
- Immobile patients can be safely transported with 3 safety belts.
- Easy to use
- MRI Compatible




STR 90 BASKET STRETCHER

- Unfolded Size (L 7 ft x W 2 ft):
- Self-weight: 20 kg
- Load Bearing: 180 kg




STR 100 NEIL ROBERTSON STRETCHER

- Self-weight : 10 kg
- Load Bearing : 155 kg



STR 110 FIBER SCOOP STRETCHER

- Unfolded Size: (L 7ft x W 1.4ft)
- Self-weight: 9 kg
- Load Bearing: 250kg
- Material: Aluminum alloy
- Uses: Transporting wounded persons,
- Apply to hospitals, stadium, ambulance, battlefield & other medical units




STR 120 MRI SCOOP STRETCHER

- Unfolded Size: (L 6ft x W 1.4ft)
- Self-weight: 8.5 kg
- Load Bearing: 250kg
- Shell: High Density Polyethylene
- Latches: Aluminum, Stainless Steel
- MRI Compatible



STR 130 CPR BOARD

- Size (60cm x 41cm x 7cm)
- Self Weight: 1.6kg
- The CPR Board provides a firm surface elasticity
- Helps to fix position and ensure the patient open airway
- The head cup maintains proper alignments with the patients head
- The board adopts PE material with no discharge contaminator.



STR 140 HEAD IMMOBILIZE

- Size (40cm x 26cm x 16cm)
- Self Weight: 1kg
- It is made of high density plastic materials, large ear hole for monitoring the patients ear canal.
- It can be used with scoop stretcher and spine board at the same time.
- Water proof plastic coating makes it easy to clean & prevents bacterial growth on it.



SMARTY EVACK CHAIR

- Size (71cm x 53cm x 94cm)
- Self Weight: 10kg
- Load Bearing: 159kg
- It is made of high stretcher aluminium alloy material
- The stretcher is foldable
- The back of stretcher comes with two foldable handles
- At the front, the stretcher comes with two flexible uplift rod.



New
Outer Pack

AB 10 BUBBLE MATTRESS

- Size (L 6ft x W 3ft x H 3in) :
- Bubble Type 3 Inch Thick Mattress
- Load Bearing: 100kg
- Comes with Air Electric Compressor



AB 20 TUBULAR MATTRESS

- Size (L 6ft x W 3ft x H 6in) :
- Tube Type 6 Inch Thick Mattress
- Load Bearing: 150kg
- Comes with Air Electric Compressor



AB 30 22" TUBULAR MATTRESS

- Size: (L 6ft x W 3ft x H 22in)
- Tube Type: 22 Inch Thick Mattress
- Load Bearing: 200kg
- Comes with Air Electric Compressor
- Foam Mattress is not needed below the airbed

HAND & ELBOW SUPPORT



RX713

Sleeve type elbow support. It helps to protect elbow. Reduce the pain caused by muscle injury. Reduce the hurts caused by friction and collision.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX723

Adjustable elbow support. It helps to protect elbow. Reduce the pain caused by muscle injury. Reduce the hurts caused by friction and collision.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX732

Helps to protect elbow. Reduces the pain caused by muscle injury and hurts caused by friction and collision.

- Good flexibility & great comfort.
- Sizes: S / M / L / XL / XXL



RX722

This product could help to protect elbow. Reduce the pain caused by muscle injury. Reduce the hurt caused by friction and collision.

- Made from superior neoprene material.
- Good flexibility & great comfort.
- Wind proof, water proof.
- Good thermal preservation & good elasticity.
- Sizes: S / M / L / XL / XXL



RX721

Adjustable Hand support. Bindingly designed. It is very helpful for protecting wrist when performing exercises. Prevent dislocation of wrist joint. Reduces the hurt caused by friction and collision.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL




RX741

En lace designed. It is very helpful for protecting wrist when performing exercises. Prevent dislocation of wrist joint. Reduce the hurts caused by friction and collision.

- Good flexibility & great comfort.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL

LEG & ANKLE SUPPORT



RX701

Avoid hurt of wrist and dislocation of wrist joint when doing strenuous exercises. Reduce the hurts caused by friction and collision.

- Good flexibility & great comfort.
- Sizes: S / M / L / XL / XXL



RX702

Avoid hurt of wrist and dislocation of wrist joint when doing strenuous exercises. Reduce the hurts caused by friction and collision.

- Good flexibility & great comfort.
- Sizes: S / M / L / XL / XXL



RX711

Bindingly designed. It is very helpful for protecting wrist when performing exercises. Prevent dislocation of wrist joint. Reduce hurt caused by friction and collision.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX743

Prevent the hurt of thighbone and the dislocation of the thigh joint. Help to prevent and quicken the speed of full recovery from thigh wrench and rheum arthritis.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX700

Protect the muscles around Knees. Prevent dislocation of knee joint. Reduce hurts caused by friction and collision. Support effect alleviate the pressure beared by tendon and joint from outside force. Uses aluminium inches inside it for better support. Sleeve type model.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Windproof & waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX728

Sleeve type knee support. Protect the muscles around Knees. Prevent dislocation of knee joint. Reduce hurts caused by friction and collision. Support effect alleviate the pressure beared by tendon and joint from outside force.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX744

Protect the muscles around knees. Prevent dislocation of knee joint. Reduce hurts caused by friction and collision.

- Good flexibility & great comfort.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX746

Helps to prevent and cure the drawing hurt and also contusion of shank muscles. Help to keep body temperature, accelerate blood circulation and speed up the full recovery.

- Good flexibility & great comfort.
- Good thermal preservation & good elasticity.
- Sizes: S / M / L / XL / XXL



RX728A

Protect the muscles around knees. Prevent dislocation of knee joint. Reduce hurts caused by friction and collision. Support effect alleviate the pressures beared by tendon and joint from outside force.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Windproof & waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX717

Prevent the hurt of anklebone and the dislocation of the ankle joint. Help to prevent and quicken the speed of full recovery from ankle wrench and rheum arthritis.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL



RX715

Helps to prevent and cure the drawing hurt and contusion of shank muscles. Help to keep body temperature, accelerate blood circulation and speed up the full recovery.

- Made from superior neoprene material.
- Good flexibility & great comfort.
- Wind proof, water proof.
- Good thermal preservation & good elasticity.
- Sizes: S / M / L / XL / XXL



RX729

Adjustable Ankle support. Openly designed. Prevent the hurt of anklebone and the dislocation of the ankle joint. Help to prevent and quicken the speed of full recovery from ankle wrench and rheum arthritis.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes: S / M / L / XL / XXL

WAIST & BACK SUPPORT



RX747

Prevent the hurt of anklebone and the dislocation of the ankle joint. It help to prevent and quicken the speed of full recovery from ankle wrench and rheum arthritis.

- Good flexibility & great comfort.
- Good thermal preservation and good elasticity.
- Sizes S/M/L/XL/XXL



RX767

Adjustable ankle support. Prevent the hurt of anklebone and the dislocation of the ankle joint. Help to prevent and quicken the speed of full recovery from ankle with the rheum arthritis.

- Made from superior Neoprene material.
- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes S/M/L/XL/XXL



RX727

Openly designed, this product's degree of tightness could be adjusted very easily. It is effective to protect the loin muscles from hurt. Used as binder in gymnastic body building exercises.

- Good flexibility & great comfort.
- Wind proof, waterproof.
- Good thermal preservation and good elasticity.
- Sizes S/M/L/XL/XXL



RX727A

Openly designed, this product's degree of tightness could be adjusted very easily. It is effective to protect the loin muscles from hurt. Used as a binder in gymnastic body building exercises.

- Made from superior neoprene material.
- Good flexibility & great comfort.
- Wind proof, water proof.
- Good thermal preservation & good elasticity.
- Sizes S/M/L/XL/XXL



RX216

Openly designed. This product's degree of tightness could be adjusted very easily. It is effective to protect the loin muscles from hurt. Used as binder in gymnastic body building exercises. Adjustable belt type model and has a slim cloth.

- Good flexibility & great comfort.
- This product has strong tensile force and elasticity in longitudinal and latitudinal way.
- To protect your waist.
- Extra length for more protection area.
- Easy to wear, suitable for any sports.
- Sizes S/M/L/XL/XXL



ARREX TRUSTED SILICONE PRODUCTS FOR FAST HEALING



BANDAGE



AR615

Designed in four steps, the product's degree of tightness could be adjusted very easily. It is effective to protect the loin muscles from hurt. Used as binder in gymnastic body building exercises.

- Made from superior neoprene material.
- Good flexibility & great comfort.
- Wind proof, water proof.
- Good thermal preservation & good elasticity.
- Sizes: S / M / L / XL / XXL



RX705

Enlace designed. This product is very helpful for protecting joints and muscles when performing exercises. Prevent dislocation of joint. Reduces the pain caused by collision.

- Good thermal preservation and good elasticity.
- Good flexibility & great comfort.

NECK SUPPORT



NS10-TRACTIONFORNECK ACE

As we know Shoulder and neck support the head. Neck tractor can support the head in suitable space, at the same time it stretches the neck slightly under the fixed situation to ease the burden for neck, removing the sore.

- Ease shoulder sore & neck painful distension.
- Strength of action : Adjusted by the user.
- Operate easily



RX704

Enlace designed. It is very helpful for protecting joints and muscles when performing exercises. Prevent dislocation of joint. Reduces the hurt caused by friction and collision.

- Good thermal preservation and good elasticity.
- Good flexibility & great comfort.

SILICONE PRODUCT



SS10-BREAST PROSTHESIS

Our silicone breast has superior quality, perfect seating and ideal shape which are resistant to chlorinated water & seawater. It is provided with special flush fitting edges all around which ensures smooth and elegant adaptation to the body.

- Sheathed in a thin, flexible very durable and polyurethane film.
- Ideal for swimming, sports and gymnastics and doesn't produce any pressure marks or sores in the skin.
- Handmade, same softness as in original tissue and doesn't cause any inconvenience.



SS20

Gel pad inside provides extra support. Walking on hard surface cause shock waves to be transmitted to the heel, Foot, Knee, Hips and the Back. Heel Cushion is ideal for treatment of PLANTAR FASCITIS and can be worn to prevent injury to the Heel bone.

- Recommended for diabetic feet, rheumatoid feet & metatarsalgia.
- Heel cushion is made up of soft non-irritable silicone.
- Sizes: I/II/III



SS30

Walking on hard surface cause shock waves to be transmitted to the heel, Foot, Knee, Hips and the Back. Heel Cushion is ideal for treatment of PLANTAR FASCITIS and can be worn to prevent injury to the Heel bone.

- Recommended for diabetic feet, rheumatoid feet & metatarsalgia.
- Heel cushion is made up of soft non-irritable silicone.
- Sizes: I/II/III



SS40

Relief for heel, ankle, knee and back. Formulated to absorb shock, vibration and provide even pressure distribution. Soft silicone with dual density dots for Full-foot comfort and protection. Multiple vents enhance air circulation.

- Ideal as replacement insole in casual or dress shoes.
- Hand Washable; physician recommended.
- Sizes: I/II/III/IV/V



SS50

Relief for heel, ankle, knee and back. Formulated to absorb shock, vibration and provide even pressure distribution. Soft silicone with dual density dots for full-foot comfort and protection. Multiple vents enhance air circulation.

- Ideal as replacement insole in casual or dress shoes.
- Hand Washable; physician recommended.
- Sizes: I/II/III/IV/V

EXERCISER



RX678

Grip Exercisers steel body with padded hand grips to increase strength & mobility in patients hand packed in pairs.

- Improves blood circulation
- Builds strong muscles in the forearm and in the body



RX772

Waist Exercisers twisting plastic is an ideal exercise for sliming & strengthening of stomach and hip region. It is small & compact which is easier to use.

- It stimulates the activeness of sympathetic nerve.
- Extends body height. It massages the meridian points on the bottom of foot.
- It improves blood circulation to all parts of human body through nerves.



RX679

A useful aid for physiotherapy and exercise of the hands, palm and fingers. Increases the strength of the patient and improves blood circulation.

- Good flexibility & great comfort.
- Suitable for post paralysis rehabilitation patients.



BREAST DEVELOPER

Used in enlarging of breasts.



RX673

Wrist and arm exercisers and ideal physiotherapy aid for exercising the wrist and forearm. Dynamic but can also be used as a static finger splint indicated for patient suffering from after effects of paralysis attack. Made of flexible still and padded with rubber. Adjustable to any size of hand. Can be used for right and left hand. Very durable and easy to use.

- Provides muscle strength to wrist & forearm.
- Prevents frequent wrist injuries.



POCKETBRA

Comes with flexible material & lace, with special design of the inner pocket for silicone breast. Sizes: 75, 80, 85, 90cm

MISCELLANEOUS



BR10

Back Rest



BR10 - X

Back Rest



BR20

Wheelchair Food Table



BR30

Oxygen Cylinder trolley for wheelchair



GR10 - GRIP BAR
(L/M/S)



GR20
BATHROOM
SIDE RAIL
(L/M/S)

GR30



UR10
URINAL FOR MALE



UR20
URINAL FOR FEMALE



UR30
URINAL WITH LID



SE10
SOFT CUSHION SEAT



SE20

WHEEL CHAIR SPARE PARTS



HOSPITAL SPARE PARTS



ABS PANEL FRONT



ABS PANEL BACK



ABS SIDE RAILING

ABS SIDE RAILINGS USED FOR EXCELSIOR BED



REMOTE



HOSPITAL BED CASTOR WHEELS



STAINLESS STEEL SIDE RAILING

MISCELLANEOUS



PATIENT HOIST



DOCTORS VISITING BAG



DEAD BODY BAG



MC30

58 (23)	45 (18)	31 (12)	54x41x34 (CM)	2.7 N.WT(kg)