

Table of contents

VENOSAN® Classes, Styles and Materials	ř
VENOSAN® Measurement Size Guide	ii
VENOSAN® Selection Chart	iv
VENOSAN® Made to Measure Stockings Fax Order Form	V
VENOSAN® Made to Measure	vi
VENOSAN® MEDICAL GRADE	
VENOSAN® 4000	2
VENOSAN® 4002 Armsleeves	6
VENOSAN® 5000	3
VENOSAN® 6000	12
VENOSAN® 7000	16
VENOSAN® 7002 Armsleeves	20
VENOSAN® Lymphoedema	22
Made to Measure Armsleeves	24
VENOSAN® 2002 Gloves	25
VENOSAN® 8000 Ulcerfit	27
VENOSAN® SUPPORT SOCKS AND STOCKINGS	
VENOSAN® MicroFiberline	31
VENOSAN® Silverline	35
VENOSAN® Legline	
VENOSAN® Athletic	43
VENOSAN® Sportsline	45
RESOXX by VENOSAN®	48
VENOSAN® AES	50
CAROSKIN™ MEDICAL GRADE FLAT KNIT STOCKINGS	
CAROSKIN™ Made to Measure Flat Knit Stockings	52
ACCESSORIES	57

VENOSAN® Classes, Styles & Materials



MILD COMPRESSION

18-21 mmHg

CLASS II

MODERATE COMPRESSION

23-32 mmHg

CLASS III

FIRM COMPRESSION

34-46 mmHg

styles



AD Below Knee



ADH Below Knee with Self Supporting Top



AGH Thigh High with Self Supporting Plain Top



AGH Thigh High with Self Supporting Lace Top



AGG Thigh Length with Waist Attachment



AT Pantyhose



ATM Maternity Pantyhose

material composition

VENOSAN® 4000: 18–21 mmHg 23–32 mmHg Custom M2M

71% Tactel® Polyamide, 29% Elastane

VENOSAN® 5000: 18–21 mmHg 23–32 mmHg Custom M2M

40% Tactel® Polyamide, 30% Elastane, 20% Cotton, 10% X-Static® – The Silver Fiber™

VENOSAN® 6000: 18–21 mmHg 23–32 mmHg

70% Nylon Polyamide, 30% Elastane

VENOSAN® 7000: 18–21 mmHg 23–32 mmHg Custom M2M (CII, CIII ONLY)

62% Tactel® Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure


VENOSAN® 8000 Ulcerfit: 23–32 mmHg 34–46 mmHg

OVERSTOCKING: 62% Tactel® Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure

LINING: 90% Polyamide, 10% Elastane

VENOSAN® Measurement Size Guide





VENOSAN® Selection Chart

Orthostatic Stress	CVI	Lipodema	PTS	Lymphoedema
The term orthostasis literally means standing upright. Orthostatic stress may be defined as the development of symptoms during upright standing relieved by recumbency.	Chronic Venous Insufficiency	Uneven distribution of fat cells in the subcutaneous regions generally in the legs, arms or abdomen	Post Thrombotic Syndrome (long term effects than can occur after venous thrombosis and their consequences)	Localized fluid retention caused by a compromised lymphatic system

Characteristics	Lipoedema	Lymphoedema
Sex	Females only	Males & Females
Age of onset	Puberty	Usually Puberty
Obesity	Yes	Variable
Symmetry	Usually symmetrical	Always symmetrical
Swollen feet	Never	Usually
Skin consistency	Normal or thinner	Thicker & firmer
Pitting of skin	Usually absent	Always pits
Easy bruising	Often	No increase
Pain & Discomfort	Often	Infrequent
Tenderness	Often	Infrequent

VENOSAN® Selection Chart



Article	VENOSAN® 7000	VENOSAN® 5000	VENOSAN® 4000
Compression Class	Class II (18–21 mmHg) Class III (23–32 mmHg) Custom M2M	Class I (18–21 mmHg) Class II (23–32 mmHg) Class III (34–46 mmHg) Custom M2M	Class I (18–21 mmHg) Class II (23–32 mmHg) Custom M2M
Material	62%TACTEL® 26% LYCRA® 8% Cotton 4% SeaCell® Pure	40%TACTEL® 30% LYCRA® 20% Cotton 10% X-STATIC® Silver	71% TACTEL® 29% LYCRA®
Characteristics	 Short Stretch High Stiffness Antibacterial Fungicidal Skin-friendly SeaCell® Pure is effective against, neurodermatitis Re-mineralizes the skin Natural fibres Open / Closed toe The feel good fibre 	 Long / Medium Stretch Antibacterial / Antifungal Silver disperses heat Enhances microcirculation Deodorant Anti-Static Temperature regulating Skin friendly Natural fibres Open / Closed toe 	 Long / Medium Stretch Climate effect/ wicks away moisture Breathable / cooler wearing Fine, almost transparent Comfortable Extremely elastic / ease of donning Skin friendly Open / Closed toe
Indication	 Moderate to severe indications Secondary varicosities CVI I, II, III Lipodema grade 2-3 	 Mild to Moderate Indications Primary varicosities CVI I, II Sweaty feet 	Mild to moderate indicationsPrimary varicositiesPost sclerotherapyCVI I

- Post-op recovery (Ulceration)
- After acute Thrombosis (after bandaging)
- PTS (Post Thrombotic Syndrome)
- Light lymphoedema
- Pronounced orthostatic stress
- Inflammatory skin problems
- Hot or cold weather
- Post-op recovery
- Light skin problems (related to veins)
- Lipodema grade 1
- During pregnancy
- Light varicose vein symptoms
- Post-op

Wearer/Target Group

- Men & Women
- People whose professions involve intensive periods of standing & demand much from their stockings
- Inflammatory skin problems
- Long term stocking wearers
- Men & Women
- Young & older long tern stocking wearers
- People that sweat
- Clients complaining of excessive heat
- People with recurring ulcers
- People with athlete's foot or foot odour problems
- Clients using walkers that complain of static electricity due to foot dragging
- Men & Women of all ages
- Slimmer legs
- Arthritic Hands
- Clients complaining of excessive heat
- Clients that sweat in nylon stockings

VENOSAN® Made to Measure Stockings Fax Order Form

MEDICAL ACCESSORIES OF AUSTRALIA 11/45 Lang Parade, Milton QLD 4064		1300 52 07 3870						ries.co es.com			
Name of patient:	Mrs	sM	r	Sender	(stam	p) Cu	stome	r #			
Date of birth:											
ast order no.:											
Ooctor:											
Customer no.:				Date:							
		Compi	ressio	n class			S	tyle			
Quality		I	II	III	AD	ADH	AG	AGH	AGG	AT	ATM Maternity
VENOSAN® 4000 VENOSAN® 5000 VENOSAN® 7000											
Colour Special requirements: Mexico Beige Black						ntity Pair Pieces Pieces		la	lain top lace top pen toe losed toe		ne dots
Length (cm) Left leg IT IK IG CF Left leg IF CF Left leg IE CB ID CD IC CC IB1 CB1 CB1 CB1 CB1 CC CA IA IA IZ		back	ents shou ith the pand if poss g has bee swelling n. All mea uld be ma g. ade stockie e returned	annot for			ht leg	Right	Length Right I IT IK IG IF IE ID IC IB1 IB IA	eg	

All measurements should be made on the bare leg. The manufacturer is not responsible for incorrect measurements. Custom made stockings cannot be returned.

VENOSAN® Made to Measure Stockings



VENOSAN® 4000

BREATHABLE

Tactel® climate effect wicks away moisture, so ensuring excellent comfort in wear.

EXTREMELY ELASTIC

Easy to put on because of the flat stress/strain curve of the knitted fabric.

COMFORTABLE

Material and visual elegance ensure optimum comfort in wear.

Material (CCL. I, II)

71% Tactel® Polyamide 29% Elastane

Colour

CCL I (Black, Mexico) CCL II (Black, Mexico, Marokko)

VENOSAN® 5000

COUNTERING ODOUR

X-Static® removes bacteria, so preventing unpleasant odours.

BALANCING FUNCTION

Optimum comfort in wear is achieved by combining X-Static® silver-coated fibre (the high conductivity of silver regulates the temperature) with Tactel® microfibre (rapidly wicks moisture away from the skin).

ANTISTATIC

Because of the high conductivity of silver, X-Static®, keeps you static-free.

Material (CCL. I, II, III)

40% Tactel® Polyamide 30% Elastane 20% Cotton 10% X-Static®-The Silver Fiber™

Colour

CCL I & III (Mexico) CCL II (Black, Mexico, Marokko)

VENOSAN® 7000

Once again we are able to present a genuine innovation: for the first time SeaCell® Pure fibre is being incorporated in VENOSAN® 7000. This cellulose fibre contains brown algae and has very special effects on human skin.

Benefits

The soft material improves comfort in wear for all those who prefer to feel natural materials against their skin.

- Antibacterial and fungicidal action (inhibits the growth of bacteria and fungi). This effect of SeaCell® Pure is long lasting, even after 150 washing cycles.
- Algae have an anti-inflammatory effect and relieve itching.
- The minerals contained in the algae replenish mineral levels in the skin.

Material (CCL, II, III)

62% Tactel® Polyamide 26% Elastane 8% Cotton 4% SeaCell® Pure

Colour

CCL I & II Black, Beige CCL III Beige

VENOSAN® Made to Measure

Stockings & Armsleeves

VENOSAN® 4000 + 7000 MTM

	min.	max.
а	18	50
У	24	52
b	16	50
b1	19	60
С	22	66
d	21	65
е	26	70
f	32	85
g	33	90
h	66	170
t	52	170

VENOSAN® 5000 MTM

	min.	max.
а	18	37
У	24	45
b	16	40
b1	19	50
С	23	56
d	21	55
е	26	62
f	32	82
g	34	88
h	68	150
t	60	150

VENOSAN® MTM ARMSLEEVES

	min.	max.
b	16	40
d	22	59
е	22	
f	22	70
g	23	80

Measurements

Measurements should be taken with the patient standing and if possible, after the leg has been rested and swelling is at a minimum. All measurements should be made on the bare leg.

Take care

- Custom made stockings may not be returned for credit.
- The manufacturer cannot be held responsible for incorrect measurements.

Guarantee

- We guarantee the premium quality of our VENOSAN® stockings!
- In the event of manufacturing or material faults the stocking will be replaced free of charge within six months.

Measuring points

- T waist
- **K** centre point of the crotch
- **H** hips (over the buttocks)
- **G** thigh (5 cm below K)
- **F** mid thigh (between K and E)
- **E** mid knee (over the kneecap)
- **D** 4 cm below the knee
- C calf at its maximum girth
- **B1** beginning of calf muscle
- **B** ankle (narrowest part above the ankle bone)
- Y heel (foot position 90°)
- A beginning of small toe

Compression Classes

I (18-21 mmHg)

II (23-32 mmHg)

III (34–46 mmHg)





Quality of Life

VENOSAN® 4000 is manufactured in a refined sandwich construction and is the first compression stocking in the world to be awarded the Tactel® climate effect label.

RREATHARI E


Tactel® climate effect wicks away moisture, so ensuring excellent comfort in wear

EXTREMELY ELASTIC

Easy to put on because of the flat stress/strain curve of the knitted fabric.

COMFORTABLE

Material and visual elegance ensure optimum comfort in wear







VENOSAN® 4000 Size Chart (cm) For an accurate fit to your exact measurements, use the chart below.



Circumfere	nce	Small	Medium	Large	X-Large
cB Ankle		20-23 cm	23-26 cm	26-29 cm	29-31 cm
cC Wides	t part of calf	29-39 cm	33-42 cm	36-45 cm	39-47 cm
cD 4 cm b	elow the knee	28-38 cm	32-41 cm	35-44 cm	38-46 cm
cG Thigh		48-63 cm	54-67 cm	58-72 cm	63-74 cm
Length					
ID AD/AD	Н	short	35-40 cm		
		long	37-42 cm		
IG AGH		short	63-72 cm		
		long	68-78 cm		
IK AGG/A	T	short	67-74 cm		
		long	74-85 cm		

Compression CLASS I		e at the ankle		rial: 71% Tactel® Polyamide Elastane	,	Tactel	OUTSTANDARD VESTINGNER VESTINGNER VESTINGNER VESTINGNER VESTINGNER VESTINGNER VESTINGNER VESTINGNER	C E
Colour	Article	Toe	Length		Small	Medium	Large	X-Large
	AD Below Knee	Open	long		41001	41002	41003	41004
SN SG Black Mexico		Closed	long		41011	41012	41013	41014
		Open	short		41021	41022	41023	41024
SM Marokko		Closed	short		41031	41032	41033	41034
	ADH Below Knee Silicone top	Open	long	self supp. plain top	41041	41042	41043	41044
SN SG Black Mexico		Closed	long	self supp. plain top	41051	41052	41053	41054
		Open	short	self supp. plain top	41061	41062	41063	41064
	\$	Closed	short	self supp. plain top	41071	41072	41073	41074
	AGH Thigh High Silicone top	Closed	long	self supp. lace top	41601	41602	41603	41604
SN SG Black Mexico		Open	long	self supp. plain top	41611	41612	41613	41614
		Closed	long	self supp. plain top	41621*	41622*	41623*	41624*
SM Marokko	/ [Closed	short	self supp. lace top	41631	41632	41633	41634
*AGH Closed long and Closed short self supporting plain		Open	short	self supp. plain top	41641	41642	41643	41644
top not available in Mexico.		Closed	short	self supp. plain top	41651*	41652*	41653*	41654*
	AT Pantyhose	Closed	long		41221	41222	41223	41224
SN SG Black Mexico		Closed	short		41231	41232	41233	41234
SM Marokko								
	ATM Maternity	Closed	long		41301	41302	41303	41304

SN

Black

SG

Mexico





Pressure at the ankle 23–32 mmHg

Material: 71% Tactel® Polyamide, 29% Elastane









02/100 22			-9	20 /0 Eldotano				
Colour	Article	Toe	Length		Small	Medium	Large	X-Large
	AD Below Knee	Open	long		42001	42002	42003	42004
SG SM Mexico Marokko		Closed	long		42011	42012	42013	42014
		Open	short		42021	42022	42023	42024
SN Black		Closed	short		42031	42032	42033	42034
	ADH Below Knee	Open	long	self supporting plain top	42041	42042	42043	42044
SG SM Marokko	Silicone Top	Closed	long	self supporting plain top	42051	42052	42053	42054
CN		Open	short	self supporting plain top	42061	42062	42063	42064
SN Black		Closed	short	self supporting plain top	42071	42072	42073	42074
	AGH Thigh High Silicone Top	Closed	long	self supporting lace top	42601	42602	42603	42604
SG SM Mexico Marokko		Open	long	self supporting plain top	42611	42612	42613	42614
		Closed	long	self supporting plain top	42621*	42622*	42623*	42624*
SN Black		Closed	short	self supporting lace top	42631	42632	42633	42634
*AGH Closed long and Closed short		Open	short	self supporting plain top	42641	42642	42643	42644
self supporting plain top not available in Mexico.		Closed	short	self supporting plain top	42651*	42652*	42653*	42654*
	AGG 🖫	Open	long		42151	42152	42153	42154
SG Mexico	RIGHT THE	Open	long		42171	42172	42173	42174
	THORN THORN							
	AT Pantyhose	Closed	long		42221	42222	42223	42224
SN SG Black Mexico		Closed	short		42231	42232	42233	42234
SM Marokko								
	ATM Maternity	Closed	long		42301	42302	42303	42304
SN SG Black Mexico		Closed	short		42311	42312	42313	42314
Black Mexico								




VENOSAN® 4002 Armsleeves

Quality of Life

world to be awarded the Tactel® climate effect label.

COMFORTABLE



Evaporation

Tactel® micro touch

Inside Tactel® multisoft

Moisture

Skin



VENOSAN® 4002 Armsleeves Size Chart (cm)

Circum	ference	Small	Medium	Large	X-Large
сC	Wrist	15–17 cm	17–20 cm	20-23 cm	23-26 cm
cD	Forearm (widest part)	22-26 cm	26-30 cm	30-34 cm	34-38 cm
cG	Upper Arm (Bicep)	24-30 cm	30-36 cm	36–42 cm	42–48 cm

VENOSAN® 4002 Armsleeves

Compression CLASS II		at the wrist mmHg	Material : 71% Tactel® Polyamide, 29% Elastane			Tactel Control of the		
Colour	Article	Length	Hand extension	Self supporting top	Small	Medium	Large	X-Large
SG Beige	CG	long	without	without	42501	42502	42503	42504





The Magic of Silver

VENOSAN® 5000 combines the natural characteristics of cotton with the antibacterial, regulating and antistatic properties of X-Static® and Tactel® high-tech yarns. A top quality product, developed for active women and men.

COUNTERING ODOUR

X-Static® removes bacteria, so preventing unpleasant odours.

BALANCING FUNCTION

Optimum comfort in wear is achieved by combining X-Static® silver-coated fibre (the high conductivity of silver regulates the temperature) with Tactel® microfibre (rapidly wicks moisture away from the skin).

ANTISTATIC

Because of the high conductivity of silver, X-Static®, keeps you static-free.







VENOSAN® 5000 Size Chart (cm) For an accurate fit to your exact measurements, use the chart below.



Circumference	Small	Medium	Large	X-Large
cB Ankle	20-23 cm	23-26 cm	26–29 cm	29-31 cm
cC Widest part of calf	29–39 cm	33–42 cm	36–45 cm	39-47 cm
cD 4 cm below the knee	28-38 cm	32-41 cm	35–44 cm	38-46 cm
cG Thigh	48–63 cm	54–67 cm	58-72 cm	63-74 cm
Length				
ID AD	long	37–42 cm		
IG AGH	long	68-78 cm		
IK AGG/AT	long	74–85 cm		



Pressure at the ankle 18–21 mmHg

Material: 40% Tactel® Polyamide, 30% Elastane, 20% Cotton, 10% X-Static®−The Silver Fiber™









Colour	Article	Toe	Length		Small	Medium	Large	X-Large
	AD Below Knee	Open	long		45001	45002	45003	45004
SN SG Black Mexico		Closed	long		45011	45012	45013	45014
SM Marokko								
	AGH Thigh High Silicone top	Closed	long	self supp. lace top	45601	45602	45603	45604
SN SG Black Mexico		Open	long	self supp. plain top	45611	45612	45613	45614
SM Marokko								
		Open	long		45151	45152	45153	45154
SG Mexico	THOM	Open	long		45171	45172	45173	45174
SN SG Black Mexico SM Marokko	AT Pantyhose	Closed	long		45221	45222	45223	45224

VENOSAN® 5002 The Magic of Silver





Pressure at the ankle 23–32 mmHg

Material: 40% Tactel® Polyamide, 30% Elastane, 20% Cotton, 10% X-Static®−The Silver Fiber™







Colour	Article	Тое	Length		Small	Medium	Large	X-Large
	AD Below Knee	Open	long		46001	46002	46003	46004
SG SM Mexico Marokko		Closed	long		46011	46012	46013	46014
SN Black								
AGF Thigh H	AGH Thigh High Silicone Top	Open	long	self supporting plain top	46611	46612	46613	46614
SG SM Mexico Marokko	Silicone Top	Closed	long	self supporting plain top	46621	46622	46623	46624
SN Black								
	AGG	Open	long		46151	46152	46153	46154
SG Mexico		Open	long		46171	46172	46173	46174
Mexico	THOM							
THE PROPERTY OF THE PROPERTY O								
	AT Pantyhose	Closed	long		46221	46222	46223	46224
SG SM Mexico Marokko								
SN								
Black	₩							

VENOSAN® 5003



Pressure at the ankle 34–46 mmHg

Material: 40% Tactel® Polyamide, 30% Elastane, 20% Cotton, 10% X-Static®−The Silver Fiber™









Colour	Article	Toe	Length	Small	Medium	Large	X-Large
	AD Below Knee	Open	long	47001	47002	47003	47004
SG SM Mexico Marokko		Closed	long	47011	47012	47013	47014
SN Black							





Affordable & Compliant

FEATURES

- Breathable
- Sits perfectly
- Complies with medical requirements
- Ease of care

BENEFITS

- Our development team took up the challenge to prove that medical compression stockings do not have to be designed only with the elderly in mind. The result is a stocking which is sheer and elegant, an indispensable companion for the modern person.
- Thanks to its sheerness and the special yarn combination, VENOSAN® 6000 offers a feeling of well-being no matter what the temperature.
- The 4 sizes and 2 lengths, combined with the highly elastic yarns, make for a stocking which sits perfectly on the leg and allows full freedom of movement in all situations.
- VENOSAN® 6000 can be machine washed, like all modern textiles today.





VENOSAN® 6000 Size Chart (cm) For an accurate fit to your exact measurements, use the chart below.



Circumference	Small	Medium	Large	X-Large
cB Ankle	20-22 cm	23-25 cm	26-28 cm	29-31 cm
cC Widest part of calf	30-39 cm	33-41 cm	36-44 cm	39–47 cm
cD 4 cm below the knee	29-38 cm	32-40 cm	35–43 cm	38-46 cm
cG Thigh	49-63 cm	54-66 cm	58-70 cm	63-74 cm
Length				
ID AD/ADH	short	34-39 cm		
	long	39-44 cm		
IG AGH	short	59-67 cm		
	long	68-78 cm		
IK AGG/AT	short	67-74 cm		
	long	74-85 cm		

VENOSAN® 6001

Compression CLASS I		e at the ankle		e rial: 70% Nylon Polyamide, Elastane			GOTE/SERVING GOTE	· (€
Colour	Article	Тое	Length		Small	Medium	Large	X-Large
	AD Below Knee	Open	long		61001	61002	61003	61004
CG CN Beige Black		Closed	long		61011	61012	61013	61014
		Open	short		61021	61022	61023	61024
		Closed	short		61031	61032	61033	61034
	AGH Thigh High Silicone top	Open	long	self supp. plain top	61601	61602	61603	61604
CG CN Beige Black		Closed	long	self supp. plain top	61611	61612	61613	61614
		Open	short	self supp. plain top	61621	61622	61623	61624
		Closed	short	self supp. plain top	61631	61632	61633	61634
	AT Pantyhose	Closed	long		61211	61212	61213	61214
CG CN Beige Black		Closed	short		61231	61232	61233	61234

Compression CLASS II		ire at the an		Material: 70% Nylon Polyamide, 30% Elastane			MOTERNICHIN METERNICHI METERNICHI METERNICHI METERNICHI METERNICHI	interest of the second of the	
Colour	Article	Toe	Length		Small	Medium	Large	X-Large	
	AD Below Knee	Open	long		62001	62002	62003	62004	
CG CN Beige Black		Open	short		62021	62022	62023	62024	
		Closed	long		62011	62012	62013	62014	
		Closed	short		62031	62032	62033	62034	
	AGH Thigh High Plain Top	Open	long	self supporting plain top	62601	62602	62603	62604	
CG CN Beige Black	Plain Top	Closed	long	self supporting plain top	62611	62612	62613	62614	
		Open	short	self supporting plain top	62621	62622	62623	62624	
		Closed	short	self supporting plain top	62631	62632	62633	62634	
	AT Pantyhose	Closed	long		62211	62212	62213	62214	
CG CN Beige Black		Closed	short		62231	62232	62233	62234	
Black Black									
	AGG	Open	long		62151	62152	62153	62154	
CG Beige	THE BEAT	Open	long		62171	62172	62173	62174	
		Open	short		62161	62162	62163	62164	
) Find a	Open	short		62181	62182	62183	62184	





Pressure at the ankle 34–46 mmHa

Material: 70% Nylon Polyamide, 30% Flastane



<u>::</u> .	GÜTEZÜSCHEN MEDESMISCHE KOMPTESSIONS- STRÜMPYE	Company ®	(F	
•	BAD	amening to Data for Standard CO Security 2000 code(1) 2013 Colors	-		

CLASS III	34-4	imm oi	1g 70% Nylon Polyamide, 30% Ela	70% Nyion Polyamide, 30% Elastane		STORY Control of the story con	
Colour	Article	Toe	Length	Small	Medium	Large	X-Large
	AD Below Knee	Open	long	63001	63002	63003	63004
CG Beige		Open	short	63021	63022	63023	63024
	AGH Thigh High Plain Top	Open	long	63601	63602	63603	63604
CG Beige	Plain Top	Open	short	63621	63622	63623	63624
20.90							
	AT Pantyhose	Closed	long	63201	63202	63203	63204
CG Beige		Closed	short	63221	63222	63223	63224
23.gc							





The Vitality of the Sea

Once again we are able to present a genuine innovation: for the first time SeaCell® Pure fibre is being incorporated in VENOSAN® 7000. This cellulose fibre contains brown algae and has very special effects on human skin.

BENEFITS

- Antibacterial and fungicidal action (inhibits the growth of bacteria and fungi). This effect of SeaCell® Pure is long lasting, even after 150 washing cycles.
- Algae have an anti-inflammatory effect and relieve itching.
- The minerals contained in the algae replenish mineral levels in the skin.
- Short stretch.





VENOSAN® 7000 Size Chart (cm) For an accurate fit to your exact measurements, use the chart below.



Circumference	Small	Medium	Large	X-Large
cB Ankle	20-23 cm	23–26 cm	26-29 cm	29–31 cm
cC Widest part of calf	29-39 cm	33–42 cm	36-45 cm	39–47 cm
cD 4 cm below the kn	ee 28-38 cm	32-41 cm	35-44 cm	38-46 cm
cG Thigh	48–63 cm	54–67 cm	58–72 cm	63-74 cm
Length				
ID AD/ADH	short	35–40 cm		
	long	37–42 cm		
IG AGH	short	63–72 cm		
	long	68-78 cm		
IK AGG/AT	short	67–74 cm		
	long	74–85 cm		

VENOSAN® 7001

Compression CLASS I	Pressure at the ankle 18–21 mmHg		Material: 62% Tactel Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure		re SeaCell Nove Tactel		CE
Colour	Article	Тое	Length	Small	Medium	Large	X-Large
AD Below Knee	Open	long	71001	71002	71003	71004	
		Closed	long	71011	71012	71013	71014
		Open	short	71021	71022	71023	71024
		Closed	short	71031	71032	71033	71034



Pressure at the ankle
23–32 mmHg

Material: 62% Tactel® Polyamide,
26% Elastane, 8% Cotton, 4% SeaCell® Pure









	_	_			_	_	_	
Colour	Article	Toe	Length		Small	Medium	Large	X-Large
	AD Below Knee	Open	long		72001	72002	72003	72004
SG SN Black		Closed	long		72011	72012	72013	72014
		Open	short		72021	72022	72023	72024
		Closed	short		72031	72032	72033	72034
	ADH Below Knee	Open	long		72041	72042	72043	72044
SG SN Black	Silicone Band	Closed	long		72051	72052	72053	72054
		Open	short		72061	72062	72063	72064
		Closed	short		72071	72072	72073	72074
	AGH Thigh High Silicone Top	Open	long	self supporting plain top	72601	72602	72603	72604
SG SN Black		Closed	long	self supporting plain top	72611	72612	72613	72614
		Open	short	self supporting plain top	72621	72622	72623	72624
		Closed	short	self supporting plain top	72631	72632	72633	72634
	AGG	Open	long		72151	72152	72153	72154
SG Beige	SIGHT LEF	Open	short		72161	72162	72163	72164
	RIGHT	Open	long		72171	72172	72173	72174
		Open	short		72181	72182	72183	72184





Pressure at the ankle
34–46 mmHg

Material: 62% Tactel® Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure









Colour	Article	Тое	Length	Small	Medium	Large	X-Large
	AD Below Knee	Open	long	73001	73002	73003	73004
SG SN Black		Closed	long	73011	73012	73013	73014
		Open	short	73021	73022	73023	73024
		Closed	short	73031	73032	73033	73034
	ADH Below Knee Silicone Band	Open	long	73041	73042	73043	73044
SG SN Black		Closed	long	73051	73052	73053	73054
		Open	short	73061	73062	73063	73064
		Closed	short	73071	73072	73073	73074
	AGH Thigh High Silicone Top	Open	long	73601	73602	73603	73604
SG SN Black		Closed	long	73611	73612	73613	73614
		Open	short	73621	73622	73623	73624
		Closed	short	73631	73632	73633	73634



VENOSAN® 7002 Armsleeves

The Vitality of the Sea

Once again we are able to present a genuine innovation: for the first time SeaCell® Pure fibre is being incorporated in VENOSAN® 7000. This cellulose fibre contains brown algae and has very special effects on human skin.

BENEFITS

- Antibacterial and fungicidal action (inhibits the growth of bacteria and fungi). This effect of SeaCell® Pure is long lasting, even after 150 washing cycles.
- Algae have an anti-inflammatory effect and relieve itching.
- The minerals contained in the algae replenish mineral levels in the skin.
- Short stretch.









Circumf	ference	Small	Medium	Large	X-Large
сС	Wrist	14–17 cm	17-20 cm	20-23 cm	23-26 cm
cD	Forearm (widest part)	22-26 cm	26-30 cm	30-34 cm	34-38 cm
cG	Upper Arm (Bicep)	24-30 cm	30-36 cm	36-42 cm	42-48 cm
Length					
C to G	Wrist to Upper Arm	short	35–45 cm		
		long	40-50 cm		

VENOSAN® 7002 Armsleeves

Compression CLASS II	Pressure at the wrist 23–32 mmHg		Material: 62% Tactel® Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure		SeaCe sea and	Pure Tactel	OUTEZEICHEN MICHTSICHE O O OCESSTÜNDIGE INTERPRE	(E
Colour	Article	Length	Hand extension	Self supporting top	Small	Medium	Large	X-Large
	CG	long	without	without	72501	72502	72503	72504
SG Beige		short	without	without	72511	72512	72513	72514
	CG-H	long	without	plain	72521	72522	72523	72524
SG Beige		short	without	plain	72531	72532	72533	72534
	ZG	Long	with	without	72541	72542	72543	72544
SG Beige		Short	with	without	72551	72552	72553	72554
SG Beige	ZG-H	Long	with	plain	72561	72562	72563	72564
		Short	with	plain	72571	72572	72573	72574



VENOSAN® Lymphoedema Care

Comfort and Care

The Lymphoedema Care range offers reflief from post-mastectomy swelling, phlebitis, radiation therapy or other edema swelling of the arm. Our Lymphoedema Armsleeve provides tissue compression to encourage lymph fluid removal and maintain limb size Gauntlets also available.

BENEFITS

- Assists with the management of Lymnphoedema of the arm.
- Helps maintain the limb size by compressing the tissue and encouraging the removal of lymph fluid.
- Gradient pressure indicated for post trauma, surgery or radiation to help maintain limb size.
- May be worn in conjunction with compression gauntlet in cases of Lymphoedema of the hand.





VENOSAN® Lymphoedema Care Size Chart (cm)



M.	Circ	cumference	Small	Medium	Large	X-Large
	cA	Palm	18–20 cm	20–22 cm	22–25 cm	25–28 cm
	- cC	Wrist	15–17 cm	17–19 cm	19-22 cm	22–26 cm
	- cD	Forearm (widest part)	23-27 cm	27-31 cm	31–35 cm	35–41 cm
	cG	Upper Arm (Bicep)	26-30 cm	30-35 cm	35-40 cm	40–45 cm

VENOSAN® Lymphoedema Care Armsleeves

Compres		Pressure at 20–30 n		rial: 75% Nylon, 25% S x free)	pandex			
Colour	Article	Length	Hand extension	Self supporting top	Small	Medium	Large	X-Large
Beige	CG	Long	without	without	LM89100	LM89101	LM89102	LM89103
Beige	ZG	Long	with	without	LM89200	LM89201	LM89202	LM89203

VENOSAN® Lymphoedema Care Gauntlet

Compress	SII	Pressure at the wrist 20–30 mmHg	Material : 80% Nylon, 20% Spand (Latex free)	ex			
Colour	Article		s	mall	Medium	Large	X-Large
Beige			LN	189300	LM89301	LM89302	LM89303

Made to Measure Armsleeves Fax Order Form



MEDICAL ACCESSORIES OF AUSTRALIA P. 1300 527 127 sales@medicalaccessories.com.au 11/45 Lang Parade, Milton QLD 4064 F. 07 3870 5944 www.medicalaccessories.com.au Name of patient: Mrs Mr Sender (stamp) Customer # Date of birth: Last order no.: Doctor: Date: Customer no.: Style Compression class Quality CG Ш Ш CG-H ZG-H СН ZG ZH w/ shoulder strap w/ hand extension VENOSAN® 7000 VENOSAN® 4000 (Beige only) Special requirements Colour left right black Quantity Right arm Left arm Length (cm) Circumference (cm) Length (cm) ΙH ΙH IG IG сG сG IF IF сF cF cЕ cЕ ΙE ΙE cD cD ID ID cС сC cA cA cZ cZ C C IA IA ΙZ ΙZ

Measure the arm when it is in the relaxed position. The manufacturer is not responsible for incorrect measurements. Custom made armsleeves cannot be exchanged.



VENOSAN® 2002 Gloves

Comfort and Ease

The VENOSAN® 2002 compression gloves are recommended for patients with Lymphoedema of the hand.

BENEFITS

- Therapeutic compression
- Comfort
- Ease of application



VENOSAN® 2002 Gloves

Compression CLASS II	Pressure at 23–32 r		Material : 77% Nylon Polyamid 23% Lycra® Elastane	e		(E
Colour	Article	Colour	Finger Extension	Small	Medium	Large
	YB	Nude	with 1/4	94541	94542	94543
VE Nude	ХВ	Nude	with 3/4	94521	94522	94523

VENOSAN® 2002 Gloves Size Chart (cm)

Palm Circumference		Small	Medium	Large
XB	Palm	16-19 cm	19-22 cm	22-25 cm



VENOSAN® 8000 Ulcerfit

Successful Therapy

VENOSAN® 8000 ULCERFIT is a compression stockings system for the treatment of leg ulcers (crural ulcer). The treatment system consists of two specially knitted stockings. They complement each other in their medical effect and generate a special pressure configuration.

BENEFITS

- The liner is normally worn day and night. This stops wound dressings from slipping.
- The overstocking supports treatment during the day by guaranteeing the defined pressure.

When treating a leg ulcer (wound) your doctor integrates the stocking system into the correct phase of wound management, which always includes a suitable dressing.





VENOSAN® 8000 Size Chart (cm) For an accurate fit to use the chart below.

For an accurate fit to your exact measurements, use the chart below.

Circumference	Small	Medium	Large	X-Large
cB Ankle	20-23 cm	23-26 cm	26-29 cm	29-31 cm
cC Widest part of calf	29-39 cm	33-42 cm	36-45 cm	39-47 cm
cD 4 cm below the knee	28-38 cm	32-41 cm	35–44 cm	38-46 cm
Length				
ID AD	short	35–40 cm		
	long	37-42 cm		

VENOSAN® 8002 Ulcerfit

*One overstocking open toe 18–21 mmHg and two white liners closed toe 10 mmHg

CLASS II

Pressure at the ankle w/Liner 18–21 mmHg

Material: 62% Tactel® Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure











Colour	Article	Length		Small	Medium	Large	X-Large
SG Beige	AD Below Knee	long 37–42 cm	14½–16½ in.	82001	82002	82003	82004
		short 35–40 cm	13¾-15¾ in.	82021	82022	82023	82024

VENOSAN® 8003 Ulcerfit

*One overstocking open toe (23–32 mmHg) and two white liners closed toe (10 mmHg)



Pressure at the ankle w/Liner 34–46 mmHg

Material: 62% Tactel® Polyamide, 26% Elastane, 8% Cotton, 4% SeaCell® Pure









Colour	Article	Length	Small	Medium	Large	X-Large
	AD Below Knee	long 37–42 cm	83001	83002	83003	83004
SG Beige		short 35–40 cm	83021	83022	83023	83024

VENOSAN® 8000 Liners



Compression 10 mmHg	Material: 90% Polyamide, 10% Elastane	(€			
Colour	Article Length	Small	Medium	Large	X-Large
SG White	AD Below Knee 35-42 cm / 13¾-16½ in.	80001	80002	80003	80004

closed toe for liners Pressure at the ankle 10 mmHg







VENOSAN® MicroFiberline

Stylish comfort, fashionably soft

VENOSAN® MicroFiberline is a graduated compression support sock with a closed toe. Micro Fiber softness and excellent breathability ensure wearer comfort. For active men and women in work and travel.

- Graduated compression
- Stylish comfortable support
- Fashionably soft comfor
- Reinforced heel and toe
- Micro Fiber softness
- Comfortable non-restricting top band



VENOSAN® MicroFiberline Men

SUPPORT	Pressure at the ankle 15–20 mmHg	Material: 82% Nylon, 10% Spandex (Latex Free), 8% Cotton	C€

Swatch	Article	Colour	Тое	Small	Medium	Large	X-Large
MC Black	AD Below Knee	Black	closed	79002	79003	79004	79005

VENOSAN® MicroFiberline Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large
Ankle	19–21.5 cm	21.5–24 cm	24–26.5 cm	26.5-32 cm
Calf	28-38 cm	33-43.5 cm	35.5–46 cm	38.5–51.5 cm

VENOSAN® MicroFiberline Shoe Size Chart

	Small	Medium	Large	X-Large
US	7 & under	7½–10	10½–12	12½+

VENOSAN® MicroFiberline Lady



SUPP	ORT		re at the ankle	Material: 90% Nylon, 10		C€		
Swate	ch	Article	Colour	Тое	Small	Medium	Large	X-Large
	MC MC	AD Black Beige	Black	closed	79221	79222	79223	79224
223			Beige	closed	79205	79206	79207	79208

VENOSAN® MicroFiberline Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large
Calf	26.5-37 cm	29-39.5 cm	32-43 cm	34.5–47 cm
Ankle	16.5-20.5 cm	20.5-24.5 cm	24.5-28 cm	28-32 cm

VENOSAN® MicroFiberline Shoe Size Chart

	Small	Medium	Large	X-Large
US	4–5	5½–7½	8–10½	11–13

VENOSAN® MicroFiberline Men

SUPPORT		re at the ankle	Material: 82% Nylon, 10	Material: 82% Nylon, 10% Spandex (Latex Free), 8% Cotton			
Swatch	Article	Colour	Тое	Small	Medium	Large	X-Large
MC Black	AD Below Knee	Black	closed	89002	89003	89004	89005

VENOSAN® MicroFiberline Men Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large
Calf	28-38 cm	33-43.5 cm	35.5–46 cm	38.5–51.5 cm
Ankle	19-21.5 cm	21.5-24 cm	24-26.5 cm	26.5-32 cm

VENOSAN® MicroFiberline Men Shoe Size Chart

	Small	Medium	Large	X-Large
US	7 & under	7½–10	10½–12	12½+







Everyday Comfort

SILVERLINE® offers long-lasting protection from foot odour and athlete's foot. X-Static®, the silver fibre, promotes hygiene and anti-microbial properties. Skin remains deodourized and free of bacteria and fungus—particularily effective in protecting against foot odour and athlete's foot. Cool and dry in Summer and warm in Winter, SILVERLINE® fibres wick moisture away from the skin for comfort in warm weather. X-Static®, the silver fibre, conducts heat in Winter to retain warmth. The silver in X-Static® yarns is an excellent conductor and eliminates static electricity.

SPECIFICALLY DESIGNED

- Enhanced comfort and fit with knitted top band
- Elimination of bacteria and fungus associated with foot odour
- Excellent regulator of heat







VENOSAN® Silverline Gentleman







Swatch	Article	Colour	Тое	Small	Medium	Large	X-Large	XX-Large
	AD Below Knee	Black	closed	49405	49406	49407	49408	49409
VS Black								

VENOSAN® Silverline Gentlemen Size Chart (cm) For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large	XX-Large
Calf	28-38 cm	33-43.5 cm	35.5-46 cm	38.5–51.5 cm	41–53 cm
Ankle	19-21.5 cm	21.5-24 cm	23-26.5 cm	26.5-32 cm	32–36 cm

VENOSAN® Silverline Gentlemen Shoe Size Chart For persons with large calf muscles, it is recommended that the next larger size be used.

	Small	Medium	Large	X-Large	XX-Large
US	7 & under	7½–10	10½–12	12½–14	13–15

VENOSAN® Silverline Lady



SUPPORT		ure at the ankle 30 mmHg	Materia 9% X-S1	l: 80% Nylon, FATIC® Silver F	11% Spandex, Plated Nylon (I	Latex Free)	J-Sta	TIC CE
Swatch	Article	Colour	Тое	Small	Medium	Large	X-Large	XX-Large
AD Below Knee	Black	closed	49505	49506	49507	49508	49509	
VS VS Black Beige		Beige	closed	49511	49512	49513	49514	49515

VENOSAN® Silverline Lady Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large	XX-Large
Ankle	16.5–20.5 cm	20.5–24.5 cm	24.5–28 cm	28-32 cm	32–36 cm
Calf	26.5–37 cm	29-39.5 cm	32-43 cm	34.5–47 cm	37–49.5 cm

VENOSAN® Silverline Lady Shoe Size Chart For persons with large calf muscles, it is recommended that the next larger size be used.

	Small	Medium	Large	X-Large	XX-Large
US	4–5	5½-7½	8–10½	11–13	14–15







Elegant Luxury

Elegant appearance for your legs, tangible prophylaxis for your health! That is VENOSAN® LEGLINE®, the effective compression stockings and support socks with an extremely fine mesh.

SPECIFICALLY DESIGNED

- Therapeutically correct, graduated compression
- Sheer good looks and superior comfort fit
- Fashionable colours






VENOSAN® Legline 20 Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.




	X-Small	Small	Medium	Large	X-Large
Ankle	20-21 cm	22-23 cm	24-25 cm	26-27 cm	28-29 cm
Calf	31–33 cm	34–36 cm	37–39 cm	40–42 cm	43–45 cm
Below Knee	29-31 cm	32–34 cm	35-37 cm	38-40 cm	41–43 cm
Mid-thigh	42–44 cm	45–47 cm	48-50 cm	51-53 cm	54-60 cm

AD & AGH



XS X-SMALL S SMALL M MEDIUM L LARGE

AT & ATM



VENOSAN® Legline 20

SUPPORT		e at the ankle mmHg	Material: 35% Lycr	65% Nylon, a Elastine				
Colour	Article	Colour	Toe	X-Small	Small	Medium	Large	X-Large
	AD Below Knee	Black	closed	303BL	304BL	305BL	306BL	
VL VL Nude		Nude	closed	303NU	304NU	305NU	306NU	
VL		Sahara	closed	303SA	304SA	305SA	306SA	
Sahara								
	AGH Thigh High Lace top	Black	closed	323BL	324BL	325BL	326BL	
VL VL Nude		Nude	closed	323NU	324NU	325NU	326NU	
VL		Sahara	closed	323SA	324SA	325SA	326SA	
Sahara								
	AT Pantyhose	Black	closed	333BL	334BL	335BL	336BL	337BL
VL VL Nude		Nude	closed	333NU	334NU	335NU	336NU	337NU
VL		Sahara	closed	333SA	334SA	335SA	336SA	337SA
Sahara								
	ATM Maternity	Black	closed	343BL	344BL	345BL	346BL	347BL
VL VL Black Sahara		Sahara	closed	343SA	344SA	345SA	346SA	347SA


VENOSAN® Legline 30 Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.




	X-Small	Small	Medium	Large	X-Large
Ankle	20-21 cm	22-23 cm	24–25 cm	26-27 cm	28-29 cm
Calf	31–33 cm	34–36 cm	37-39 cm	40-42 cm	43–45 cm
Below Knee	29-31 cm	32-34 cm	35–37 cm	38-40 cm	41–43 cm
Mid-thigh	42-44 cm	45–47 cm	48–50 cm	51–53 cm	54–60 cm

AGH



XS X-SMALL S SMALL M MEDIUM L LARGE

ΑT



 XS
 X-SMALL

 S
 SMALL

 M
 MEDIUM

 L
 LARGE

 XL
 X-LARGE (AT ONLY)

VENOSAN® Legline 30

SUPPORT		ure at the ankle			% Nylon, 35%l Polyamid (La		x Free)	C€
Colour	Article	Colour	Тое	X-Small	Small	Medium	Large	X-Large
	AD Below Knee	Black	closed	403BL	404BL	405BL	406BL	
VL Black Nude		Nude	closed	403NU	404NU	405NU	406NU	
		Sahara	closed	403SA	404SA	405SA	406SA	
VL Sahara								
	AGH Thigh High Lace Top	Black	closed	423BL	424BL	425BL	426BL	
VL VL Nude		Nude	closed	423NU	424NU	425NU	426NU	
		Sahara	closed	423SA	424SA	425SA	426SA	
VL Sahara								
	AT Pantyhose	Black	closed	433BL	434BL	435BL	436BL	437BL
VL VL Nude		Nude	closed	433NU	434NU	435NU	436NU	437NU
		Sahara	closed	433SA	434SA	435SA	436SA	437SA
VL Sahara								



VENOSAN® Athletic

Performance Socks for Men and Women

Designed for everyday wear and broad gender appeal. Whether the need for athletic, travel or daily living, you will be excited these compression therapy options exist. Cushioned foot beds and balloon toe pockets, provide day – long benefits of comfort and relief from the effects of venous disease and / or daily living.

MATERIAL

75% Nylon, 12% Polyester, 13% Spandex. Made in the USA.

INDICATIONS

Everyday wear, travel, athletic need, mild oedema, prevention and treatment of mild venous symptoms, heaviness and fatigue in the legs.

BENEFITS

- Unisex compression socks
- Reciprocated Heel and toe for outstanding fit and comfort
- Seamless toe, ventilated foot for breathability and arch support
- Cushioned achilles tendon area for better protection and comfort



VENOSAN® Athletic Unisex Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.



SUPPORT

Pressure at the ankle 20-30 mmHg

Material: 75% Nylon, 12% Polyester, 13% Spandex







Swatch	Article	Colour	Toe	Small	Medium	Large	X-Large
W B	AD Below Knee	White	closed	A59202W	A59203W	A59204W	A59205W
White Black		Black	closed	A59202B	A59203B	A59204B	A59205B
F P		Pink	closed	A59202F	A59203F	A59204F	A59205F
Pink Purple		Purple	closed	A59202P	A59203P	A59204P	A59205P
S Blue		Blue	closed	A59202S	A59203S	A59204S	A59205S

VENOSAN® Athletic Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large
Calf	29-38 cm	30-42 cm	32-46 cm	34–50 cm
Ankle	18-21 cm	22-25 cm	25-29 cm	29-33 cm

VENOSAN® Athletic Shoe Size Chart

	Small	Medium	Large	X-Large
Ladies' US	6-8	8½-11	11½–13	13+
Men's US	7 & under	7½-10	10½-12	12½+



VENOSAN® SportsLine

Performance Support

VENOSAN® SportsLine is a graduated compression support sock for active men and women. Coolmax® draws moisture away from the skin and helps regulate your body temperature to keep you cool.


WHY COOLMAX?

Coolmax® is the original performance fabric that's been helping top athletes push harder and longer but is now being discovered for the cool comfort it brings to everyday apparel.

Coolmax® has proven physiological benefits that, depending on your needs, can be applied to or found in a whole new range of clothing styles – from moderate and high-intensity performance apparel to everyday casual wear.

So consumers get the level of performance they need to stay cool and comfortable in any situation.







VENOSAN® SportsLine Unisex

SUPPORT		Pressure at the ankle Material: 65% Coolmax® 9% Spandex (Latex Free)			6 Nylon,		
Swatch	Article	Colour	Тое	Size III	Size IV	Size V	Size VI
	AD Below Knee	White	closed	69202	69203	69204	69205
SP SP Black		Black	closed	89202	89203	89204	89205

VENOSAN® SportsLine Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large
Calf	28.6-38 cm	30.5-42 cm	31.8-46 cm	34-50cm
Ankle	17.8–21 cm	22.2–25.5 cm	25-29 cm	29-33 cm

VENOSAN® SportsLine Shoe Size Chart

	Small	Medium	Large	X-Large
Ladies' US	6-8	8½-11	11½-13	13+
Men's US	7 & under	7½-10	10½-12	12½+

VENOSAN® SportsLine Unisex



SUPPORT	Pressure at 20–30 n		Material: 65% (Nylon, 9% Spar			COOL MAX	□ □ (€
Swatch	Article	Colour	Тое	Small	Medium	Large	X-Large
	AD Below Knee	White	closed	79202	79203	79204	79205
SP SP White		Black	closed	99202	99203	99204	99205

VENOSAN® SportsLine Size Chart (cm)

For an accurate fit to your exact measurements, use the chart below.

	Small	Medium	Large	X-Large
Ankle	18–21 cm	22–25 cm	25-29 cm	29-33 cm
Calf	29-38 cm	30-42 cm	32-46 cm	34-50 cm

VENOSAN® SportsLine Shoe Size Chart

	Small	Medium	Large	X-Large
Ladies' US	6-8	8½-11	11½–13	13+
Men's US	7 & under	7½-10	10½-12	12½+



RESOXX® by VENOSAN®

Pure Performance

Sport socks are the product of many decades of experience and competence in the manufacture of medical compression stockings together with a patented and tested energy pressure point called the Support Energy Point (S.E.P.). This means the specialist expertise of the world market leader **VENOSAN®** combined with S.E.P. is now available for sports.

The **RESOXX® Socks** combine their very latest sock construction made of ultramodern high tech fibres, a special knitting process, and a simple design, thereby best meeting the needs of athletes.

RESOXX® have a degressive pressure gradient from the ankle upwards that has been specially adapted for sports. The blood circulation in the legs is accelerated to provide optimum support for the calf pump. This is ensured by a maximum pressure on the ankle (approx. 25 mmHg) which decreases by controlled reduction to 50% at the calf, thereby helping athletes break down metabolic substances (e.g. lactate) more quickly and the muscle fibres are supplied more rapidly with new oxygen-enriched blood.

- Regenerative compression
- With silver ions, antibacterial
- Friction free movement
- «Red» energy pressure point (S.E.P.)



RESOXX Socks by VENOSAN® Size Chart



US Shoe Size	6–7½		8-	8–9½ 10–		-12 13-		-14
Europe Shoe Size	36–38		39–41		42–44		45–47	
Calf Circumference (cm)	31–37	38–44	32–38	39–45	33–39	40–46	34-40	41–48
	S	S-Plus	M	M-Plus	L	L-Plus	XL	XL-Plus

RESOXX Tubes by **VENOSAN®** Size Chart

Circumference	Small	Medium Large		X-Large
Ankle (cm)	20–23	23–26	26–29	29–31
Calf (cm)	29–39	33–42	36–45	39–47

RESOXX® by VENOSAN® Socks



RESOXX® by VENOSAN® Tubes

SUPPORT		e at the ankle 2 mmHg	Material: 40% Tactel® Poly 20% Cotton, 10% X-Static®	cane, Tactel	
Colour	Article	Small	Medium	Large	X-Large
BL WH Black White	Tubes	02S	02M	02L	02X



VENOSAN® AES

Thrombosis and embolism prevention Anti-embolism stockings

The anti-embolism stockings are very easy to put on and offer a high level of comfort in wear due to their breathability.

PERFECT FIT

Anatomical shaping and comprehensive choice of sizes.

COLOUR-CODED

Allows easy and reliable size recognition after washing.

SAFE

Because of the top inspection window.



VENOSAN® AES Thrombosis prevention and Anti-embolism stockings



SUPPORT

Pressure at the ankle 18 mmHg

Material: 90% Nylon Polyamide, 10% Spandex (Latex Free) CE

rticle & watch	Ankle Circumference	Calf Circumference	Thigh Circumference	Lengt	h	Size		Size colour code (toe)	Length colour code (top)	Art. No
AD Below knee	18-20 cm	to 32 cm		to	40.5 cm	S	regular	yellow	natural	SG57000
				over	40.5 cm	S	long	yellow	blue	SG57001
	20-22 cm	32-37 cm		to	43 cm	M	regular	natural	natural	SG57002
				over	43 cm	M	long	natural	blue	SG57003
	23-25 cm	37-42 cm		to	45.5 cm	L	regular	blue	natural	SG57004
				over	45.5 cm	L	long	blue	blue	SG57005
	25-28 cm	42-47 cm		to	45.5 cm	XL	regular	green	natural	SG57006
G				over	45.5 cm	XL	long	green	blue	SG57007
hite	28-30 cm	47-52 cm		to	45.5 cm	XXL	regular	purple	natural	SG57008
				over	45.5 cm	XXL	long	purple	blue	SG57009
G	18-20 cm	to 32 cm	45-64 cm	to	73.5 cm	S	short	yellow	yellow	SG57020
Length					73.5-78.5 cm	S	regular	yellow	natural	SG5702
				over	79 cm	S	long	yellow	blue	SG5702
	20-22 cm	32-37 cm	53-71 cm	to	73.5 cm	M	short	natural	yellow	SG5702
					73.5-78.5 cm	M	regular	natural	natural	SG5702
B				over	79 cm	M	long	natural	blue	SG5702
	23-25 cm	37-42 cm	61-80 cm	to	73.5 cm	L	short	blue	yellow	SG5702
G .					73.5-78.5 cm	L	regular	blue	natural	SG5702
hite				over	79 cm	L	long	blue	blue	SG5702
	25-28 cm	42-47 cm	70-84 cm	to	78.5 cm	XL	regular	green	natural	SG5702
				over	79 cm	XL	long	green	blue	SG5703
GG	18-20 cm	to 32 cm	45-64 cm	to	73.5 cm	S	short	yellow	yellow	SG5705
h Belt					73.5-78.5 cm	S	regular	yellow	natural	SG5705
				over	79 cm	S	long	yellow	blue	SG5705
\	20-22 cm	32-37 cm	53-71 cm	to	73.5 cm	M	short	natural	yellow	SG5705
					73.5-78.5 cm	M	regular	natural	natural	SG5705
				over	79 cm	M	long	natural	blue	SG5705
	23-25 cm	37-42 cm	61-80 cm	to	73.5 cm	L	short	blue	yellow	SG5705
G					73.5-78.5 cm	L	regular	blue	natural	SG5705
SG White				over	79 cm	L	long	blue	blue	SG5705
	25-28 cm	42-47 cm	70-84 cm	to	78.5 cm	XL	regular	green	natural	SG5705
				over	79 cm	XL	long	green	blue	SG5706

CAROSKIN™



Made-to-Measure Flat Knit Compression Stockings for the treatment of Lymphoedema and Lipoedema

CAROSKIN™ compression stockings are a high quality made-to-measure flat knit compression stocking. Using flat knit technology with seam the CAROSKIN™ compression stockings are designed to sit smoothly on the skin and optimise patient comfort.

The unique CAROSKIN™ flat knit design exerts a flat uniform graduated compression on the limb area to be treated, which is positively acting for the preventive treatment of oedema.

FEATURES AND BENEFITS

- optimal fitting stocking
- excels in high sweat absorption
- increases wearing comfort
- uniform, graduated compression



CAROSKIN™ flat knit compression stockings offer a choice of two finishes, both available in beige or black colors:

CAROSKIN™ Comfort

The CAROSKIN™ Comfort flat knit compression stockings are available in Class I, Class II or Class III compression. This stocking is designed to significantly increase the patient's comfort.

CAROSKIN™ Strong

The CAROSKIN™ Strong flat knit compression stockings are available in Class I, Class II or Class III compression. If you're after a stocking with a shorter stretch feature then the CAROSKIN™ Strong is your choice.

Class I	Mild Compression	20 mmHG
Class II	Moderate Compression	27 mmHG
Class III	Firm Compression	35 mmHG



CAROSKIN™ Flat Knit Compression Stockings

- Compression Class I, II, III
- UV protection

CAROSKIN™ Comfort Material 79% Crabyon, 21% Lycra

CAROSKIN™ Strong Material 79% Polyamide, 21% Lycra

CAROSKIN™ Comfort / Strong Lower Limb


CODE	DESCRIPTION
4030001	CAROSKIN™ Panty AT
4031001	CAROSKIN™ Panty AT1
4034001	CAROSKIN™ Thigh High AG
4034001	CAROSKIN™ Thigh High and Waist AGG
4035001	CAROSKIN™ Below Knee AD
	-



	CODE	DESCRIPTION
	5018003	Grip Top 3cm AD - AG
ST	5018005	Grip Top 5cm AD - AG
IS LI	5018015	AGG Band
OPTIONS LIST	5019001	Slant end / Slipform
	5019002	Comfort Zone Knee
	5019003	Closed Toe
	5019004	Slant foot / Oblique

CAROSKIN™ Comfort / Strong Armsleeve

CODE	DESCRIPTION
4010001	CAROSKIN™ Armsleeve with shoulder cap
4012001	CAROSKIN™ Armsleeve CG



All-in-one sleeves available. Please advise at time of order.

	CODE	DESCRIPTION
LIST	5018003	Grip top band 3cm
SNS	5018005	Grip top band 5cm
OPTIONS	5019001	Slant end / Slipform
0	5019002	Comfort zone elbow

CAROSKIN™ Comfort / Strong Glove

CODE	DESCRIPTION
4016001	CAROSKIN™ Glove Open Finger
4017001	CAROSKIN™ Half Glove Open Thumb



All-in-one sleeves available. Please advise at time of order.

	CODE	DESCRIPTION
LIST	N/A	Wrist without pressure
SN	N/A	Glove and sleeve in one piece
PTIO		

LOWER LIMB





Patient name Contact Measured by				Customer	code	_ Date
AM	- ₄ - M RIGHT		▶	LEFT	(Circumference) Length	E D
AL AG AK	K G			G	AG	B1 B A2 AA2 AA2 AA2 Seen from the side
AE AD AC	E D C			E	AE AD AC	A1 H AA2 H AA2 Seen from the side (interior)
AB1 AB AA2	B1 B A2 A2	A1	H A1	B1 A2 A2	AB1 AB AA2	Right Left X X XX XX XXX XXX
Panty AT Special re	Panty AT1	Thigh High AG	Thigh and Waist AGG	Below Knee AD	Strong ST Class Cl. Toes Options Sil AC Sil Sil	001 Beige

ARMSLEEVE





Patient name Contact Measured by	DOB Customer code Signature	
GH RIGHT H CG F CF E CE Circumfe Leng	erence	G CG F CF CE CE
Armsleeve CG Armsleeve shoulder of Special requirements	Strong ST0 Class Class Quantity Righ Options Glov Silco Silco Slan	1 Beige C003 Black 1 Beige ST03 Black 1 Class II Class III 1 Left Pe and sleeve in one piece 1 band 3cm 1 one band 5cm 1 tend / Slipform 1 fort zone elbow

GLOVE





Patient name		Sex Male Female
Contact		Date
W2 V3 V4 V5 AC X1 Circumferer Length	x4 x3 x2	W4 W5 W5 A B
	olf Glove Strong Strong Class Coptions	C001 Beige C003 Black ST01 Beige ST03 Black Class I Class II Class III Wrist without pressure Glove and sleeve in one piece Left glove Right glove



Mediven Butler Off Compression Stocking Remover

This long handled device assists with the removal of compression stockings, socks and normal pantyhose. It is lightweight and is gentle on the skin as it has no sharp edges. It features a T-shaped thick grip handle, metal shaft and a curved tongue-like base with a hook. It can also double as a shoehorn.

Product Code: 900-A



Easy-Slide Compression Arm Stocking Applicator

An application aid for open fingered compression stockings. The easy slide is cone shaped and made of a slippery, paper like material with a pull loop at the tapered end. The Easy-Slide is pulled onto the arm and the stocking is then pulled over it. The Easy-Slide is removed by pulling the loop.

Product Codes: ES10 - S ES40 - XL



Eureka-On Compression Stocking Applicator

A compression stocking applicator for the leg or arm. It consists of a triangular shaped durable paper like slippery material with a central extended pull tab and eyelet holes on the sides for threading a thin tubular plastic lacing pin. Can also be used with pantyhose. Requires the use of two hands and for user to be able to reach their feet allows the compression garment to be rolled onto the 'doughnut' device ready for use.

Product Code: EUREKA - ON



Ezy-As Compression Stocking Applicator

A rigid plastic structure designed to assist with fitting closed and open pressure garments to upper and lower limbs. It is available in three sizes. Available in three sizes, with an optional handle attachment.

Product Codes: YELLOW SUT - SM 260 MED

BLUE SUT - SM 261 LG RED SUT - SM 259 SM

VENOSAN® Socks-Jet Donning and Doffing Aid

The patented donning and doffing aid for medical compression stockings. Features easy to use doffing and donning bars for stocking application. Simple for use by nursing staff, also for bedridden patients.

Product Code: AV002



Doff N' Donner

Designed to assist with the application of open and closed toe compression stockings. Suitable for fitting and removing garments on the upper and lower limbs below and above knee. Features two components: a 'doughnut' shaped soft silicon device, filled with a leak proof fluid and a cone with suction cup base that allows the compression garment to be rolled onto the 'doughnut' device ready for use.

Product Codes: DDCONE

DDDEVICE





VENOSAN® Silicone Dotted Donning Gloves

VENOSAN® silicone dotted donning gloves assist with applying stockings. Helps protect stockings from tears by finger nails. Available in sizes S/M and L/XL.

Product Codes: VC003 S/M

VC004 L/XL



VENOSAN® Coltex Roll-On Body Adhesive

The perfect way to hold up your hosiery - 24 hours a day. Use Coltex Adhesive, and you can be absolutely certain your stockings stay up at all times. Whatever the shape of your legs and whatever compression stockings you want to keep in place, this gentle but effective glue is a real confidence booster. Easy to apply and machine washable, Coltex Adhesive also contains hamamelis extracts to gently soothe your skin.



COMPRESSION THERAPY - SAMA®









Sama® Frottee Tubular Padding

- 95% cotton, 3% polyamide, 2% elastane
- Soft, time-saving tubular padding
- Multi-elastic applies with no wrinkling, creases or folds
- Snug and accurate fit contours to all shapes, sizes and lengths
- Re-usable, will endure up to 20 washes at 40 degrees celcius depending on exudation, sterilisable
- Cut to desired length, no fraying
- Available in 10 metre rolls in four sizes (S, M, L, XL) and in four colours (white, pink, blue, green)

Functional benefits:

- Easy to apply, time-saving regime
- Increased comfort from padding layer, will not crease or wrinkle
- Allows the skin to breathe
- All products are washable and re-usable
- Cut and only use what you need from padding roll

Code	Size	Colour	Quantity
DV88103	Small	White	10m roll
DV88104	Medium	White	10m roll
DV88105	Large	White	10m roll
DV88106	X-Large	White	10m roll
DV88113	Small	Pink	10m roll
DV88114	Medium	Pink	10m roll
DV88115	Large	Pink	10m roll
DV88123	Small	Blue	10m roll
DV88124	Medium	Blue	10m roll
DV88125	Large	Blue	10m roll
DV88133	Small	Green	10m roll
DV88134	Medium	Green	10m roll
DV88135	Large	Green	10m roll

Size Chart

Small	15 - 35 cm circumference
Medium	30 - 55 cm circumference
Large	40 - 75 cm circumference
X-Large	50 - 95 cm circumference (White only)

